


# DAQmx for HTBasic

*Professional National Instruments DAQ Driver*

## User's Guide






# **DAQmx for HTBasic**

---

## **User's Guide**

---



**Tech Soft**

[www.TechSoft.de](http://www.TechSoft.de)

Copyright 2005-2011 by Tech Soft GmbH, Berlin, Germany.

8th revision, January 2011

Printed in Germany. All rights reserved.

Information in this manual is subject to change without notice and does not represent a commitment on the part of the vendor. Contents of this publication may not be reproduced in any form without express written permission of Tech Soft GmbH.

The software described in this manual is furnished under a license agreement and may be used or copied only in accordance with the terms of the agreement.

High Tech Basic, HTBasic and TransEra are trademarks of TransEra corporation.  
HP, HP BASIC and Hewlett Packard are trademarks of Hewlett Packard Company.  
Windows, Windows 95, 98, Me, NT, 2000, XP, Windows 2003, Windows 2008, Vista and Windows 7 are trademarks of Microsoft Corporation. All other brand and product names are trademarks or registered trademarks of their respective companies.

Manual part number: H-TN0010

## **Table of contents**

<b>CONCEPT .....</b>	<b>2</b>
<b>INSTALLATION .....</b>	<b>2</b>
Getting started .....	3
<b>FUNCTION REFERENCE: OVERVIEW .....</b>	<b>4</b>
NI-DAQmx Function Overview .....	4
Determining HTBasic function names for NI-DAQmx functions.....	11
Special HTBasic Support Functions.....	11
Global variables for HTBasic.....	11
<b>FUNCTION DECLARATIONS FOR ALL IMPLEMENTED FUNCTIONS .....</b>	<b>13</b>
Task Configuration and Control.....	13
Channel Configuration / Creation .....	16
Timing Functions .....	36
Read Functions .....	40
Write Functions.....	45
Triggering Functions.....	51
Calibration.....	55
NI-DAQmx Specific Attribute Get/Set/Reset Function Declarations .....	77
Switch Functions .....	79
Signal Routing .....	82
TEDS .....	83
Events & Signals .....	84
Buffer .....	85
Other Get/Set specific functions (Part I).....	87
Other Get/Set specific functions (Part II).....	98
Callback services.....	165
Watchdog functions.....	167
Storage .....	168
HTBasic specific functions .....	169
<b>INDEX.....</b>	<b>173</b>

## **DAQmx Library for HTBasic**

© Tech Soft GmbH 2005-2011. All rights reserved.

### System requirements:

TransEra HTBasic 9.X

Windows 2000, Windows XP, Windows Vista

Hardware: Any National Instruments DAQ board which is supported by the NI-DAQmx driver (and simulated devices)

National Instruments system software "NI-DAQmx" version 8.0 or higher (tested up to NI-DAQmx 9.1)

Strongly recommended: NI-DAQmx 8.6.1

## **Concept**

HTBasic includes a National Instruments driver for older DAQ boards which uses the so called "Traditional NI-DAQ (Legacy)" driver model. This HTBasic driver is called "DAQNI".

However, many newer DAQ boards are no longer supported by this older driver model. Some time ago National Instruments introduced a new driver for all newer DAQ hardware called "NI-DAQmx. This driver simplifies the programming of the DAQ hardware significantly and provides a lot of new functionality.

The DAQmx Library for HTBasic allows you to use all the features of the new NI-DAQmx driver model under HTBasic including powerful functions like callback services which uses the built-in intelligence of the NI DAQ hardware allowing you to run multiple simultaneous measurements in the background (if the hardware supports it). The HTBasic library supports almost the entire function set of the NI-DAQmx driver (currently more than 200 functions) and provides a programming interface very similar to the ANSI C interface of the NI-DAQmx driver. In addition to the normal analog and digital input/output and counter functions the new driver also supports switches, special measurements like microphone measures, accelerator measures, special signal triggering, watchdog services, task management and much more (the hardware must support the function, e.g. not all National Instruments hardware includes a watchdog).

It is important that you understand the key concepts of the National Instruments NI-DAQmx driver before you start to use the DAQmx Library for HTBasic. Please read the National Instruments documentation "NI-DAQmx Help" and "NI-DAQmx C Reference Help" carefully.

## **Installation**

Install hardware and software from National Instruments as described in your National Instruments product manual. Note that you do not have to use physical hardware in order to test this library, it works with simulated devices too. Please refer to your MAX (Measurement & Automation Explorer) documentation to find out how to set up simulated devices in MAX.

After you have installed the National Instruments software and (if desired) hardware first test your devices in MAX (Measurement & Automation Explorer) using the test panel.

If the test finished successfully you can install the "DAQmx Library for HTBasic" by running the file `DAQmx_htb_setup.exe`. Follow the instructions on the screen. After successful installation start HTBasic and load the desired example program from the "File" menu (see below). Make sure that you have enabled the option "Change MSI on Open".

## **Example programs**

Many examples programs can be found in the subdirectory "HTB Examples" in the installation directory of your DAQmx for HTBasic Library. One of those examples is called "`combined_daqmx_example.prg`" which demonstrates the following basic functions:

- Analog output at port `ao0`
- Fast analog voltage input at port `ai0` and current input at Port `ai1`
- Digital output on bits 0..3 of digital port 0 (`port0`)
- Digital input on bits 0..3 of digital port 0 (`port1`)
- Counter I/O with counter 0 (`ctr0`)
- Optional: Usage of two stored tasks (`VoltageTask`, `DigitalInTask`)

You might connect the analog output port 0 (`ao0`) with analog input port 0 (`ai0`) in order to read back the analog output values in to your analog input port. Also, you might want to connect bits 0..3 of digital port 0 with bits 4..7 of digital port 1 for testing purposes.

Please note that the demo version only supports 2 simultaneous physical tasks (but it supports a virtually unlimited number of tasks using simulated devices).

If you have questions or problems please report them to:

### **Tech Soft GmbH**

Test & Measurement Division

Karmeliterweg 114

13465 Berlin, Germany

Email: [HTB@TechSoft.de](mailto:HTB@TechSoft.de)

<http://www.TechSoft.de>

## Introduction

The DAQmx Library for HTBasic was developed to use new DAQ interfaces like the PCI-6259 from National Instruments which are no longer supported by the current HTBasic DAQNI driver.

This library uses the Windows driver from National Instruments "NI-DAQmx" (`nicaiu.dll`). Due to the fact that the new driver model uses data types which are not available under HTBasic we developed a "wrapper DLL" which sits between the National Instruments driver and HTBasic.

The HTBasic driver consists of two files:

<code>daqmx_htb.dll</code>	Driver DLL ("Wrapper DLL")
<code>daqmx_htb.lib</code>	HTBasic library
<code>shfolder.dll</code>	Helper functions to determine special Windows directories

## Files after installation

<code>HTB Examples\*.*</code>	Many HTBasic example programs
<code>Driver\*.*</code>	Driver DLLs and HTBasic library
<code>Doc\DAQmx_htb_manual.pdf</code>	This document
<code>Source\*.*</code>	Source files for wrapper DLL for Visual C++ 6.0 *)

\*) Complete VC++ source code is only available in the special "Developer Version".

## Library usage in customer programs

In order to simplify the usage of the the function in customer programs we have developed several support functions. At first you should load the function `FNHx_manage_subs` from the main HTBasic library (make sure that your MSI points to the directory where the file "`daqmx_htb.lib`" is located):

```
IF NOT INMEM("FNHx_manage_subs") THEN LOADSUB FNHx_manage_subs FROM "daqmx_htb.lib"
```

After that call the loaded function and check the return value:

```
IF NOT FNHx_manage_subs("LOAD",Verbose) THEN STOP ! Load & initialize library
```

If the return value is 0 then the function failed for some reason. In this case the program will be stopped. If the function returned 1 all went fine and you can start to use the function library. In the example program there is a label called `User_routines:` where several user routines are called (analog I/O, digital I/O etc.).

After all functions are finished you should use the following lines at the very end of your program:

```
Dummy=FNHx_manage_subs("UNLOAD",Verbose) ! Unload and de-initialize library  
DELSUB FNHx_manage_subs
```

## Getting started

We recommend you to go through the source code of the example program in order to see how the library works and how to use it. It is important that you understand how the original functions of the National Instruments DAQmx library work, so always look into the documentation for this library since the HTBasic library is based on the original National Instruments documentation.

If you want to write your own applications a good starting point is the framework "`Program_framework.prg`". Make a copy of this file and use this as a starting point for your own programs. Go to label `User_routines:` and include your own code there.

## Device info

Use the function `FNHx_get_dev_info` in order to print some device specific information on the screen including device names, serial number, driver versions etc. This function is used in the example program and the application framework program too. When you set the variable `write_2_file` to 1 the function writes your current configuration into the file `daqmx_devinfo.txt` stored in the directory `[%APPDATA%]\HTB-DAQmx`. On English Windows systems this directory usually translates into `"c:\Documents and Settings\Username\Application Data"`.

## Function Reference: Overview

Function names under HTBasic are usually shorter than the original NI-DAQmx functions due to the fact that HTBasic only supports function names with up to 15 characters. All functions use `_cdecl` calling convention.

Functions without a specific NI-DAQmx version number in the following table can be used under NI-DAQmx 8.0. Some functions are only available under certain NI-DAQmx versions. Those functions require at least the NI-DAQmx version number installed on your system.

### NI-DAQmx Function Overview

NI-DAQmx function name *)	HTBasic function name	HTBasic Ret.type	NI-DAQmx Version
<b>1. Task Configuration and Control</b>			
DAQmxLoadTask	FNHx_loadtask	LONG	
DAQmxCreateTask	FNHx_createtask	LONG	
DAQmxStartTask	FNHx_starttask	LONG	
DAQmxStopTask	FNHx_stoptask	LONG	
DAQmxClearTask	FNHx_cleartask	LONG	
DAQmxWaitUntilTaskDone	FNHx_waituntaskdn	LONG	
DAQmxIsTaskDone	FNHx_istaskdone	LONG	
DAQmxTaskControl	FNHx_taskcontrol	LONG	
DAQmxGetNthTaskChannel	FNHx_getnthtaskch	LONG	
<b>2. Channel Configuration/Creation</b>			
DAQmxCreateAIVoltageChan	FNHx_craivoltchan	LONG	
DAQmxCreateAIVoltageRMSChan	FNHx_craivlrmschn	LONG	v8.6
DAQmxCreateAICurrentChan	FNHx_craicurrchan	LONG	
DAQmxCreateAICurrentRMSChan	FNHx_craicurmschn	LONG	v8.6
DAQmxCreateAIThrmcplChan	FNHx_craithcpchan	LONG	
DAQmxCreateAIRTDChan	FNHx_crairtchan	LONG	
DAQmxCreateAIThrmstrChanIex	FNHx_craithmchiex	LONG	
DAQmxCreateAIThrmstrChanVex	FNHx_craithmchvex	LONG	
DAQmxCreateAIFreqVoltageChan	FNHx_craifrvchan	LONG	
DAQmxCreateAIResistanceChan	FNHx_crairesichan	LONG	
DAQmxCreateAISTrainGageChan	FNHx_craistrgachn	LONG	
DAQmxCreateAIVoltageChanWithExcit	FNHx_craivochwexc	LONG	
DAQmxCreateAITempBuiltInSensorChan	FNHx_craitmpbschn	LONG	
DAQmxCreateAIACcelChan	FNHx_craiaaccelchn	LONG	
DAQmxCreateAIMicrophoneChan	FNHx_craimicrchan	LONG	
DAQmxCreateAIPosLVDTChan	FNHx_craiplvdtch	LONG	
DAQmxCreateAIPosRVDTChan	FNHx_craipsrvdtch	LONG	
DAQmxCreateTEDSAIVoltageChan	FNHx_crtaivoltchn	LONG	
DAQmxCreateTEDSAICurrentChan	FNHx_crtaicurrchn	LONG	
DAQmxCreateTEDSAIThrmcplChan	FNHx_crtaithcpchn	LONG	
DAQmxCreateTEDSAIRTDChan	FNHx_crtairtdchn	LONG	
DAQmxCreateTEDSAIThrmstrChanIex	FNHx_crtaitchchiex	LONG	
DAQmxCreateTEDSAIThrmstrChanVex	FNHx_crtaitchchvex	LONG	
DAQmxCreateTEDSAIResistanceChan	FNHx_crtairesichn	LONG	
DAQmxCreateTEDSAISTrainGageChan	FNHx_crtaistrgchn	LONG	
DAQmxCreateTEDSAIVoltageChanWithExcit	FNHx_crtaivochwexc	LONG	
DAQmxCreateTEDSAIACcelChan	FNHx_craitacclchn	LONG	
DAQmxCreateAIMicrophoneChan	FNHx_craimicrchan	LONG	
DAQmxCreateAIPosLVDTChan	FNHx_craiplvdtch	LONG	
DAQmxCreateAIPosRVDTChan	FNHx_craipsrvdtch	LONG	
DAQmxCreateAOVoltageChan	FNHx_craovoltchan	LONG	
DAQmxCreateAOCCurrentChan	FNHx_craocurrchan	LONG	
DAQmxCreateAOFuncGenChan	FNHx_craofungnchn	LONG	v8.6
DAQmxCreateDIChan	FNHx_createdichan	LONG	
DAQmxCreateDOChan	FNHx_createdochan	LONG	
DAQmxCreateCIFreqChan	FNHx_crcifreqchan	LONG	
DAQmxCreateCIPeriodChan	FNHx_crciperichan	LONG	
DAQmxCreateCICountEdgesChan	FNHx_crcicntedgch	LONG	
DAQmxCreateCIPulseWidthChan	FNHx_crciplswchan	LONG	
DAQmxCreateCISemiPeriodChan	FNHx_crcisemperch	LONG	
DAQmxCreateCITwoEdgeSepChan	FNHx_crcitwedsech	LONG	
DAQmxCreateCILinEncoderChan	FNHx_crcilinenccch	LONG	
DAQmxCreateCIAngEncoderChan	FNHx_crcianencchn	LONG	
DAQmxCreateCIGPSTimestampChan	FNHx_crcigpstschn	LONG	
DAQmxCreateCOPulseChanFreq	FNHx_crcopulchfrq	LONG	
DAQmxCreateCOPulseChanTime	FNHx_crcopulchtim	LONG	
DAQmxCreateCOPulseChanTicks	FNHx_crcopulchtck	LONG	

\*) As described in the "NI-DAQmx C Reference Help"



NI-DAQmx function name *)	HTBasic function name	HTBasic Ret.type	NI-DAQmx Version
<b>3. Timing Functions</b>			
DAQmxCfgSampClkTiming	FNHx_cfgsamclktim	LONG	
DAQmxCfgPipelinedSampClkTiming	FNHx_cfgplsclktim	LONG	v8.5
DAQmxCfgHandshakingTiming	FNHx_cfghshtiming	LONG	
DAQmxCfgBurstHandshakingTimingImportClock	FNHx_cfgbhtimiclk	LONG	
DAQmxCfgBurstHandshakingTimingExportClock	FNHx_cfgbhtimeclk	LONG	
DAQmxCfgChangeDetectionTiming	FNHx_cfgchgdettim	LONG	
DAQmxCfgImplicitTiming	FNHx_cfgimptiming	LONG	
<b>4. Read Functions</b>			
DAQmxReadAnalogF64	FNHx_readanalgf64	LONG	
DAQmxReadAnalogScalarF64	FNHx_readanalsf64	LONG	
DAQmxReadBinaryI16	FNHx_rdbinaryi16	LONG	
DAQmxReadBinaryI32	FNHx_rdbinaryi32	LONG	
DAQmxReadDigitalU16	FNHx_readdigitu16	LONG	
DAQmxReadDigitalU32	FNHx_readdigitu32	LONG	
DAQmxReadDigitalScalarU32	FNHx_readdigits32	LONG	
DAQmxReadDigitalLines	FNHx_readdiglines	LONG	
DAQmxReadCounterF64	FNHx_readcntrf64	LONG	
DAQmxReadCounterU32	FNHx_readcntru32	LONG	
DAQmxReadCounterScalarU32	FNHx_rdcntscalu32	LONG	
DAQmxReadCounterScalarF64	FNHx_rdcntscaLf64	LONG	
<b>5. Write Functions</b>			
DAQmxWriteAnalogF64	FNHx_writeanalf64	LONG	
DAQmxWriteAnalogScalarF64	FNHx_writeanasf64	LONG	
DAQmxWriteBinaryI16	FNHx_writebini16	LONG	
DAQmxWriteBinaryI32	FNHx_writebini32	LONG	
DAQmxWriteDigitalU16	FNHx_writdigitu16	LONG	
DAQmxWriteDigitalU32	FNHx_writdigitu32	LONG	
DAQmxWriteDigitalScalarU32	FNHx_writdigscu32	LONG	
DAQmxWriteDigitalLines	FNHx_writdiglines	LONG	
DAQmxWriteCtrFreq	FNHx_wrctrfreq	LONG	
DAQmxWriteCtrFreqScalar	FNHx_wrctrfrqscal	LONG	
DAQmxWriteCtrTime	FNHx_wrcrttime	LONG	
DAQmxWriteCtrTimeScalar	FNHx_wrcrtimscal	LONG	
DAQmxWriteCtrTicks	FNHx_wrcrtticks	LONG	
DAQmxWriteCtrTicksScalar	FNHx_wrcrttickssc	LONG	
<b>6. Triggering</b>			
DAQmxDisableStartTrig	FNHx_disstarttrig	LONG	
DAQmxCfgDigEdgeStartTrig	FNHx_cfdgedsttrig	LONG	
DAQmxCfgAnlgEdgeStartTrig	FNHx_cfganedstrig	LONG	
DAQmxCfgAnlgWindowStartTrig	FNHx_cfganwistrig	LONG	
DAQmxCfgDigPatternStartTrig	FNHx_cfgdiptstrig	LONG	
DAQmxDisableRefTrig	FNHx_disreftrig	LONG	
DAQmxCfgDigEdgeRefTrig	FNHx_cfgdiedrtrig	LONG	
DAQmxCfgAnlgEdgeRefTrig	FNHx_cfganedrtrig	LONG	
DAQmxCfgAnlgWindowRefTrig	FNHx_cfganwirtrig	LONG	
DAQmxCfgDigPatternRefTrig	FNHx_cfgdipartrig	LONG	
DAQmxDisableAdvTrig	FNHx_disadvtrig	LONG	
DAQmxCfgDigEdgeAdvTrig	FNHx_cfgdiedatrig	LONG	
DAQmxSendSoftwareTrigger	FNHx_sendsofttrig	LONG	
<b>7. Calibration</b>			
DAQmxSelfCal	FNHx_selfcal	LONG	
DAQmxPerformBridgeOffsetNullingCal	FNHx_perfbrofnlca	LONG	
DAQmxGetCalDevTemp	FNHx_getcadevtemp	LONG	
DAQmxGetSelfCalLastTemp	FNHx_getsfcaltemp	LONG	
DAQmxGetExtCalLastTemp	FNHx_getexcaltemp	LONG	
DAQmxESeriesCalAdjust	FNHx_esercaladjst	LONG	
DAQmxMSeriesCalAdjust	FNHx_msercaladjst	LONG	
DAQmxSSeriesCalAdjust	FNHx_ssercaladjst	LONG	
DAQmxSCBaseboardCalAdjust	FNHx_scbbcaldadjst	LONG	
DAQmxAOSeriesCalAdjust	FNHx_aosecaladjst	LONG	
DAQmxDeviceSupportsCal	FNHx_devsuppcal	LONG	
DAQmxGetSelfCalSupported	FNHx_selfcalsupp	LONG	
DAQmxInitExtCal	FNHx_initextcal	LONG	
DAQmxCloseExtCal	FNHx_closeextcal	LONG	
DAQmxGetSelfCalLastDateAndTime	FNHx_getscaldtati	LONG	
DAQmxGetExtCalLastDateAndTime	FNHx_getecaldtati	LONG	
DAQmxRestoreLastExtCalConst	FNHx_reslexcalcon	LONG	
DAQmxChangeExtCalPassword	FNHx_chnextcalpwd	LONG	
DAQmxGetExtCalRecommendedInterval	FNHx_getecalrecin	LONG	

NI-DAQmx function name *)	HTBasic function name	HTBasic Ret.type	NI-DAQmx Version
<b>7. Calibration (Continued)</b>			
DAQmxAdjustDSAAICal	FNHx_adjdsaaical	LONG	
DAQmxAdjustDSAAOCal	FNHx_adjdsaaocal	LONG	
DAQmxAdjustDSATimebaseCal	FNHx_adjdsatbcal	LONG	
DAQmxAdjust4204Cal	FNHx_adj4204cal	LONG	
DAQmxAdjust4220Cal	FNHx_adj4220cal	LONG	
DAQmxAdjust4224Cal	FNHx_adj4224cal	LONG	
DAQmxAdjust4225Cal	FNHx_adj4225cal	LONG	
DAQmxSetup1126Cal	FNHx_set1126cal	LONG	v8.3
DAQmxAdjust1126Cal	FNHx_adj1126cal	LONG	v8.3
DAQmxSetup1141Cal	FNHx_set1141cal	LONG	v8.3
DAQmxAdjust1141Cal	FNHx_adj1141cal	LONG	v8.3
DAQmxSetup1142Cal	FNHx_set1142cal	LONG	v8.3
DAQmxAdjust1142Cal	FNHx_adj1142cal	LONG	v8.3
DAQmxSetup1143Cal	FNHx_set1143cal	LONG	v8.3
DAQmxAdjust1143Cal	FNHx_adj1143cal	LONG	v8.3
DAQmxSetup153xCal	FNHx_set153xcal	LONG	v8.3
DAQmxAdjust153xCal	FNHx_adj153xcal	LONG	v8.3
DAQmxSetup1540Cal	FNHx_set1540cal	LONG	v8.3
DAQmxAdjust1540Cal	FNHx_adj1540cal	LONG	v8.3
DAQmxSetup1102Cal	FNHx_set1102cal	LONG	v8.5
DAQmxAdjust1102Cal	FNHx_adj1102cal	LONG	v8.5
DAQmxSetup1104Cal	FNHx_set1104cal	LONG	v8.5
DAQmxAdjust1104Cal	FNHx_adj1104cal	LONG	v8.5
DAQmxSetup1112Cal	FNHx_set1112cal	LONG	v8.5
DAQmxAdjust1112Cal	FNHx_adj1112cal	LONG	v8.5
DAQmxSetup1122Cal	FNHx_set1122cal	LONG	v8.7
DAQmxAdjust1122Cal	FNHx_adj1122cal	LONG	v8.7
DAQmxSetup1124Cal	FNHx_set1124cal	LONG	v8.5
DAQmxAdjust1124Cal	FNHx_adj1124cal	LONG	v8.5
DAQmxSetup1125Cal	FNHx_set1125cal	LONG	v8.5
DAQmxAdjust1125Cal	FNHx_adj1125cal	LONG	v8.5
DAQmxSetup1502Cal	FNHx_set1502cal	LONG	v8.5
DAQmxAdjust1502Cal	FNHx_adj1502cal	LONG	v8.5
DAQmxSetup1503Cal	FNHx_set1503cal	LONG	v8.5
DAQmxAdjust1503Cal	FNHx_adj1503cal	LONG	v8.5
DAQmxAdjust1503CurrentCal	FNHx_adj1503ccl	LONG	v8.5
DAQmxSetup1520Cal	FNHx_set1520cal	LONG	v8.5
DAQmxAdjust1520Cal	FNHx_adj1520cal	LONG	v8.5
DAQmxSetup1521Cal	FNHx_set1521cal	LONG	v8.7
DAQmxAdjust1521Cal	FNHx_adj1521cal	LONG	v8.7
<b>8. NI-DAQmx Specific Attribute Get/Set/Reset Function Declarations</b>			
DAQmxGetPhysicalChanName	FNHx_getphyschnam	LONG	
DAQmxSetCIPrescaler	FNHx_setciprescal	LONG	
<b>9. Device Control</b>			
DAQmxResetDevice	FNHx_resetdevice	LONG	
<b>10. Switch Functions</b>			
DAQmxSwitchCreateScanList	FNHx_swcrscanlist	LONG	
DAQmxSwitchConnect	FNHx_swconnect	LONG	
DAQmxSwitchConnectMulti	FNHx_swconnectmul	LONG	
DAQmxSwitchDisconnect	FNHx_swdisconnect	LONG	
DAQmxSwitchDisconnectMulti	FNHx_swdisconnmul	LONG	
DAQmxSwitchDisconnectAll	FNHx_swdisconnall	LONG	
DAQmxSwitchSetTopologyAndReset	FNHx_swsettopares	LONG	
DAQmxSwitchFindPath	FNHx_swfindpath	LONG	
DAQmxSwitchOpenRelays	FNHx_swopenrelays	LONG	
DAQmxSwitchCloseRelays	FNHx_swclosrelays	LONG	
DAQmxSwitchWaitForSettling	FNHx_swwaitforset	LONG	
<b>11. Signal Routing</b>			
DAQmxConnectTerms	FNHx_connterms	LONG	
DAQmxDisconnectTerms	FNHx_disconnterms	LONG	
DAQmxTristateOutputTerm	FNHx_tristoutterm	LONG	
<b>12. TEDS</b>			
DAQmxConfigureTEDS	FNHx_configteds	LONG	
DAQmxClearTEDS	FNHx_clearteds	LONG	

\*) As described in the "NI-DAQmx C Reference Help"

NI-DAQmx function name *)	HTBasic function name	HTBasic Ret.type	NI-DAQmx Version
<b>13. Events &amp; Signals</b>			
DAQmxExportSignal	FNHx_exportsignal	LONG	
<b>14. Buffer</b>			
DAQmxGetBufInputBufSize	FNHx_getbufinpsiz	LONG	
DAQmxSetBufInputBufSize	FNHx_setbufinpsiz	LONG	
DAQmxGetBufOutputBufSize	FNHx_getbufoutsiz	LONG	
DAQmxSetBufOutputBufSize	FNHx_setbufoutsiz	LONG	
DAQmxGetBufOutputOnbrdBufSize	FNHx_getobufoutsz	LONG	
DAQmxSetBufOutputOnbrdBufSize	FNHx_setobufoutsz	LONG	
<b>15. Other Get/Set specific functions</b>			
DAQmxGetAIMax	FNHx_getaimax	LONG	
DAQmxGetAIMin	FNHx_getaimin	LONG	
DAQmxGetAILowpassCutoffFreq	FNHx_getailpcofrq	LONG	
DAQmxSetAILowpassCutoffFreq	FNHx_setailpcofrq	LONG	
DAQmxGetAILowpassEnable	FNHx_getailwpenab	LONG	
DAQmxSetAILowpassEnable	FNHx_setailwpenab	LONG	
DAQmxGetAnlgEdgeStartTrigHyst	FNHx_getanedsthys	LONG	
DAQmxSetAnlgEdgeStartTrigHyst	FNHx_setanedsthys	LONG	
DAQmxGetAIACExcitSyncEnable	FNHx_getaiacexsen	LONG	
DAQmxSetAIACExcitSyncEnable	FNHx_setaiacexsen	LONG	
DAQmxGetReadOverloadedChansExist	FNHx_getrdovchexi	LONG	
DAQmxGetReadOverloadedChans	FNHx_getrdovchans	LONG	
DAQmxGetDigLvlPauseTrigSrc	FNHx_getdiglvptsc	LONG	
DAQmxSetDigLvlPauseTrigSrc	FNHx_setdiglvptsc	LONG	
DAQmxGetPauseTrigType	FNHx_getpaustrtyp	LONG	
DAQmxSetPauseTrigType	FNHx_setpaustrtyp	LONG	
DAQmxGetDigLvlPauseTrigWhen	FNHx_getdlpatrwhn	LONG	
DAQmxSetDigLvlPauseTrigWhen	FNHx_setdlpatrwhn	LONG	
DAQmxGetCIGPSSyncSrc	FNHx_getcigpssysr	LONG	
DAQmxSetCIGPSSyncSrc	FNHx_setcigpssysr	LONG	
DAQmxGetHshkDelayAfterXfer	FNHx_gethsdelaxfr	LONG	
DAQmxSetHshkDelayAfterXfer	FNHx_sethsdelaxfr	LONG	
DAQmxGetExportedHshkEventInterlockedDeassertDelay	FNHx_getexhseiddl	LONG	
DAQmxSetExportedHshkEventInterlockedDeassertDelay	FNHx_setexhseiddl	LONG	
DAQmxGetDIDigFltrEnable	FNHx_getdidgflten	LONG	
DAQmxSetDIDigFltrEnable	FNHx_setdidgflten	LONG	
DAQmxGetDIDigFltrMinPulseWidth	FNHx_getdidgflmpw	LONG	
DAQmxSetDIDigFltrMinPulseWidth	FNHx_setdidgflmpw	LONG	
DAQmxGetAdvTrigType	FNHx_getadvtrgtyp	LONG	
DAQmxSetAdvTrigType	FNHx_setadvtrgtyp	LONG	
DAQmxGetSwitchScanRepeatMode	FNHx_getswscrpmod	LONG	
DAQmxSetSwitchScanRepeatMode	FNHx_setswscrpmod	LONG	
DAQmxGetMasterTimebaseSrc	FNHx_getmastbasrc	LONG	
DAQmxSetMasterTimebaseSrc	FNHx_setmastbasrc	LONG	
DAQmxGetMasterTimebaseRate	FNHx_getmastbrate	LONG	
DAQmxSetMasterTimebaseRate	FNHx_setmastbrate	LONG	
DAQmxGetRefClkSrc	FNHx_getrefclksrc	LONG	
DAQmxSetRefClkSrc	FNHx_setrefclksrc	LONG	
DAQmxGetRefClkRate	FNHx_getrefclkrat	LONG	
DAQmxSetRefClkRate	FNHx_setrefclkrat	LONG	
DAQmxGetSampClkTimebaseSrc	FNHx_getsclktbsrc	LONG	
DAQmxSetSampClkTimebaseSrc	FNHx_setsclktbsrc	LONG	
DAQmxGetSyncPulseSrc	FNHx_getsynpulsrc	LONG	
DAQmxSetSyncPulseSrc	FNHx_setsynpulsrc	LONG	
DAQmxGetAIThrmcplScaleType	FNHx_getaitcsctyp	LONG	v8.3
DAQmxSetAIThrmcplScaleType	FNHx_setaitcsctyp	LONG	v8.3
DAQmxResetAIThrmcplScaleType	FNHx_rstaitcsctyp	LONG	v8.3
DAQmxGetAIIsTEDS	FNHx_getaiisteds	LONG	v8.3
DAQmxGetDILogicFamily	FNHx_getdilogifam	LONG	v8.3
DAQmxSetDILogicFamily	FNHx_setdilogifam	LONG	v8.3
DAQmxResetDILogicFamily	FNHx_rstdilogifam	LONG	v8.3
DAQmxGetDIMemMapEnable	FNHx_getdimmapena	LONG	v8.3
DAQmxSetDIMemMapEnable	FNHx_setdimmapena	LONG	v8.3
DAQmxResetDIMemMapEnable	FNHx_rstdimmapena	LONG	v8.3
DAQmxGetDIAcquireOn	FNHx_getdiacquion	LONG	v8.3
DAQmxSetDIAcquireOn	FNHx_setdiacquion	LONG	v8.3
DAQmxResetDIAcquireOn	FNHx_rstdiacquion	LONG	v8.3
DAQmxGetDOLineStatesStartState	FNHx_getdolsstrts	LONG	v8.3
DAQmxSetDOLineStatesStartState	FNHx_setdolsstrts	LONG	v8.3
DAQmxResetDOLineStatesStartState	FNHx_rstdolsstrts	LONG	v8.3

NI-DAQmx function name *)	HTBasic function name	HTBasic Ret.type	NI-DAQmx Version
<b>15. Other Get/Set specific functions (Continued 2)</b>			
DAQmxGetDOLineStatesPausedState	FNHx_getdolspauss	LONG	v8.3
DAQmxSetDOLineStatesPausedState	FNHx_setdolspauss	LONG	v8.3
DAQmxResetDOLineStatesPausedState	FNHx_rstdolspauss	LONG	v8.3
DAQmxGetDOLineStatesDoneState	FNHx_getdolsdone	LONG	v8.3
DAQmxSetDOLineStatesDoneState	FNHx_setdolsdone	LONG	v8.3
DAQmxResetDOLineStatesDoneState	FNHx_rstdolsdone	LONG	v8.3
DAQmxSetDigitalLogicFamilyPowerUpState	FNHx_setdlfampust	LONG	v8.3
DAQmxGetDOLogicFamily	FNHx_getdologifam	LONG	v8.3
DAQmxSetDOLogicFamily	FNHx_setdologifam	LONG	v8.3
DAQmxResetDOLogicFamily	FNHx_rstdologifam	LONG	v8.3
DAQmxGetDOMemMapEnable	FNHx_getdommapena	LONG	v8.3
DAQmxSetDOMemMapEnable	FNHx_setdommapena	LONG	v8.3
DAQmxResetDOMemMapEnable	FNHx_rstdommapena	LONG	v8.3
DAQmxGetDOGenerateOn	FNHx_getdogeneron	LONG	v8.3
DAQmxSetDOGenerateOn	FNHx_setdogeneron	LONG	v8.3
DAQmxResetDOGenerateOn	FNHx_rstdogeneron	LONG	v8.3
DAQmxSetPhysicalChanName	FNHx_setphyschnam	LONG	v8.3
DAQmxGetExportedPauseTrigOutputTerm	FNHx_getexptouttm	LONG	v8.3
DAQmxSetExportedPauseTrigOutputTerm	FNHx_setexptouttm	LONG	v8.3
DAQmxResetExportedPauseTrigOutputTerm	FNHx_rstexptouttm	LONG	v8.3
DAQmxGetExportedPauseTrigLvlActiveLvl	FNHx_getexptlaclv	LONG	v8.3
DAQmxSetExportedPauseTrigLvlActiveLvl	FNHx_setexptlaclv	LONG	v8.3
DAQmxResetExportedPauseTrigLvlActiveLvl	FNHx_rstexptlaclv	LONG	v8.3
DAQmxGetExportedRefTrigPulsePolarity	FNHx_getexrfrppo	LONG	v8.3
DAQmxSetExportedRefTrigPulsePolarity	FNHx_setexrfrppo	LONG	v8.3
DAQmxResetExportedRefTrigPulsePolarity	FNHx_rstexrfrppo	LONG	v8.3
DAQmxGetExportedStartTrigPulsePolarity	FNHx_getexsttrppo	LONG	v8.3
DAQmxSetExportedStartTrigPulsePolarity	FNHx_setexsttrppo	LONG	v8.3
DAQmxResetExportedStartTrigPulsePolarity	FNHx_rstexsttrppo	LONG	v8.3
DAQmxGetExportedRdyForXferEventDeassertCond	FNHx_getexrdfxedc	LONG	v8.3
DAQmxSetExportedRdyForXferEventDeassertCond	FNHx_setexrdfxedc	LONG	v8.3
DAQmxResetExportedRdyForXferEventDeassertCond	FNHx_rstexrdfxedc	LONG	v8.3
DAQmxGetExportedRdyForXferEventDeassertCondCustomThreshold	FNHx_getexrddyfxed	LONG	v8.3
DAQmxSetExportedRdyForXferEventDeassertCondCustomThreshold	FNHx_setexrddyfxed	LONG	v8.3
DAQmxResetExportedRdyForXferEventDeassertCondCustomThreshold	FNHx_rstexrddyfxed	LONG	v8.3
DAQmxGetExportedDataActiveEventOutputTerm	FNHx_getexdaceout	LONG	v8.3
DAQmxSetExportedDataActiveEventOutputTerm	FNHx_setexdaceout	LONG	v8.3
DAQmxResetExportedDataActiveEventOutputTerm	FNHx_rstexdaceout	LONG	v8.3
DAQmxGetExportedDataActiveEventLvlActiveLvl	FNHx_getexdaevlal	LONG	v8.3
DAQmxSetExportedDataActiveEventLvlActiveLvl	FNHx_setexdaevlal	LONG	v8.3
DAQmxResetExportedDataActiveEventLvlActiveLvl	FNHx_rstexdaevlal	LONG	v8.3
DAQmxGetExportedRdyForStartEventOutputTerm	FNHx_getexrfseotm	LONG	v8.3
DAQmxSetExportedRdyForStartEventOutputTerm	FNHx_setexrfseotm	LONG	v8.3
DAQmxResetExportedRdyForStartEventOutputTerm	FNHx_rstexrfseotm	LONG	v8.3
DAQmxGetExportedRdyForStartEventLvlActiveLvl	FNHx_getexrfselal	LONG	v8.3
DAQmxSetExportedRdyForStartEventLvlActiveLvl	FNHx_setexrfselal	LONG	v8.3
DAQmxResetExportedRdyForStartEventLvlActiveLvl	FNHx_rstexrfselal	LONG	v8.3
DAQmxGetDevProductCategory	FNHx_getdevprdcacat	LONG	v8.3
DAQmxGetDevChassisModuleDevNames	FNHx_getdevchmdnm	LONG	v8.3
DAQmxGetDevAnlgTrigSupported	FNHx_getdevatrsup	LONG	v8.3
DAQmxGetDevDigTrigSupported	FNHx_getdevdtrsups	LONG	v8.3
DAQmxGetDevAIMaxSingleChanRate	FNHx_getdevaimscr	LONG	v8.3
DAQmxGetDevAIMaxMultiChanRate	FNHx_getdevaimmcr	LONG	v8.3
DAQmxGetDevAIMinRate	FNHx_getdevaimrat	LONG	v8.3
DAQmxGetDevAISimultaneousSamplingSupported	FNHx_getdevasssup	LONG	v8.3
DAQmxGetDevAITrigUsage	FNHx_getdevaitrus	LONG	v8.3
DAQmxGetDevAIVoltageRngs	FNHx_getdevavrang	LONG	v8.3
DAQmxGetDevAIVoltageIntExcitDiscreteVals	FNHx_getdevaviedv	LONG	v8.3
DAQmxGetDevAIVoltageIntExcitRangeVals	FNHx_getdevavierv	LONG	v8.3
DAQmxGetDevAICurrentRngs	FNHx_getdevacurra	LONG	v8.3
DAQmxGetDevAICurrentIntExcitDiscreteVals	FNHx_getdevaciedv	LONG	v8.3
DAQmxGetDevAIFreqRngs	FNHx_getdevafrqrn	LONG	v8.3
DAQmxGetDevAIGains	FNHx_getdevaigain	LONG	v8.3
DAQmxGetDevAICouplings	FNHx_getdevaicoup	LONG	v8.3
DAQmxGetDevAILowpassCutoffFreqDiscreteVals	FNHx_getdevalcfdv	LONG	v8.3
DAQmxGetDevAILowpassCutoffFreqRangeVals	FNHx_getdevalcfrv	LONG	v8.3
DAQmxGetDevAOPhysicalChans	FNHx_getdevaophch	LONG	v8.3
DAQmxGetDevAOSampClkSupported	FNHx_getdevaoscsps	LONG	v8.3
DAQmxGetDevAOMaxRate	FNHx_getdevaomart	LONG	v8.3



NI-DAQmx function name *)	HTBasic function name	HTBasic Ret.type	NI-DAQmx Version
<b>15. Other Get/Set specific functions (Continued 3)</b>			
DAQmxGetDevAOMinRate	FNHx_getdevaomirt	LONG	v8.3
DAQmxGetDevAOTrigUsage	FNHx_getdevaotrus	LONG	v8.3
DAQmxGetDevAOVoltageRngs	FNHx_getdevaovorg	LONG	v8.3
DAQmxGetDevAOCurrentRngs	FNHx_getdevaocurg	LONG	v8.3
DAQmxGetDevAOGains	FNHx_getdevaogain	LONG	v8.3
DAQmxGetDevDIMaxRate	FNHx_getdevdimart	LONG	v8.3
DAQmxGetDevDITrigUsage	FNHx_getdevditrus	LONG	v8.3
DAQmxGetDevDOMaxRate	FNHx_getdevdomart	LONG	v8.3
DAQmxGetDevDOTrigUsage	FNHx_getdevdotrus	LONG	v8.3
DAQmxGetDevCITrigUsage	FNHx_getdevcitrus	LONG	v8.3
DAQmxGetDevCISampClkSupported	FNHx_getdevciscsp	LONG	v8.3
DAQmxGetDevCIMaxTimebase	FNHx_getdevcimatb	LONG	v8.3
DAQmxGetDevCOTrigUsage	FNHx_getdevcotrus	LONG	v8.3
DAQmxGetDevCOMaxTimebase	FNHx_getdevcomatb	LONG	v8.3
DAQmxGetDevNumDMACHans	FNHx_getdevnmdchn	LONG	v8.3
DAQmxGetDevCompactDAQChassisDevName	FNHx_getdevcdcdnm	LONG	v8.3
DAQmxGetDevCompactDAQSlotNum	FNHx_getdevcdslnm	LONG	v8.3
DAQmxGetTaskNumDevices	FNHx_gettsknumdev	LONG	v8.3
DAQmxGetSampClkUnderflowBehavior	FNHx_getsmclufbeh	LONG	v8.3
DAQmxSetSampClkUnderflowBehavior	FNHx_setsmclufbeh	LONG	v8.3
DAQmxResetSampClkUnderflowBehavior	FNHx_rstsmclufbeh	LONG	v8.3
DAQmxGetAIConvRateEx	FNHx_getaicvratex	LONG	v8.3
DAQmxSetAIConvRateEx	FNHx_setaicvratex	LONG	v8.3
DAQmxResetAIConvRateEx	FNHx_rstaicvratex	LONG	v8.3
DAQmxGetAIConvMaxRateEx	FNHx_getaicvmrtex	LONG	v8.3
DAQmxGetAIConvSrcEx	FNHx_getaicvsrcecx	LONG	v8.3
DAQmxSetAIConvSrcEx	FNHx_setaicvsrcecx	LONG	v8.3
DAQmxResetAIConvSrcEx	FNHx_rstaicvsrcecx	LONG	v8.3
DAQmxGetAIConvActiveEdgeEx	FNHx_getaicvaceex	LONG	v8.3
DAQmxSetAIConvActiveEdgeEx	FNHx_setaicvaceex	LONG	v8.3
DAQmxResetAIConvActiveEdgeEx	FNHx_rstaicvaceex	LONG	v8.3
DAQmxGetAIConvTimebaseDivEx	FNHx_getaicvtbdex	LONG	v8.3
DAQmxSetAIConvTimebaseDivEx	FNHx_setaicvtbdex	LONG	v8.3
DAQmxResetAIConvTimebaseDivEx	FNHx_rstaicvtbdex	LONG	v8.3
DAQmxGetAIConvTimebaseSrcEx	FNHx_getaicvtbsex	LONG	v8.3
DAQmxSetAIConvTimebaseSrcEx	FNHx_setaicvtbsex	LONG	v8.3
DAQmxResetAIConvTimebaseSrcEx	FNHx_rstaicvtbsex	LONG	v8.3
DAQmxGetDelayFromSampClkDelayUnits	FNHx_getdlfscdunt	LONG	v8.3
DAQmxSetDelayFromSampClkDelayUnits	FNHx_setdlfscdunt	LONG	v8.3
DAQmxResetDelayFromSampClkDelayUnits	FNHx_rstdlfscdunt	LONG	v8.3
DAQmxGetDelayFromSampClkDelayUnitsEx	FNHx_getdlfscduex	LONG	v8.3
DAQmxSetDelayFromSampClkDelayUnitsEx	FNHx_setdlfscduex	LONG	v8.3
DAQmxResetDelayFromSampClkDelayUnitsEx	FNHx_rstdlfscduex	LONG	v8.3
DAQmxGetDelayFromSampClkDelay	FNHx_getdlfscdlay	LONG	v8.3
DAQmxSetDelayFromSampClkDelay	FNHx_setdlfscdlay	LONG	v8.3
DAQmxResetDelayFromSampClkDelay	FNHx_rstdlfscdlay	LONG	v8.3
DAQmxGetDelayFromSampClkDelayEx	FNHx_getdlfscdlex	LONG	v8.3
DAQmxSetDelayFromSampClkDelayEx	FNHx_setdlfscdlex	LONG	v8.3
DAQmxResetDelayFromSampClkDelayEx	FNHx_rstdlfscdlex	LONG	v8.3
DAQmxGetDigPatternPauseTrigSrc	FNHx_getdigpptsrsrc	LONG	v8.3
DAQmxSetDigPatternPauseTrigSrc	FNHx_setdigpptsrsrc	LONG	v8.3
DAQmxResetDigPatternPauseTrigSrc	FNHx_rstdigpptsrsrc	LONG	v8.3
DAQmxGetDigPatternPauseTrigPattern	FNHx_getdigpptpat	LONG	v8.3
DAQmxSetDigPatternPauseTrigPattern	FNHx_setdigpptpat	LONG	v8.3
DAQmxResetDigPatternPauseTrigPattern	FNHx_rstdigpptpat	LONG	v8.3
DAQmxGetDigPatternPauseTrigWhen	FNHx_getdigpptwhn	LONG	v8.3
DAQmxSetDigPatternPauseTrigWhen	FNHx_setdigpptwhn	LONG	v8.3
DAQmxResetDigPatternPauseTrigWhen	FNHx_rstdigpptwhn	LONG	v8.3
DAQmxGetWriteNextWriteIsLast	FNHx_getwrnxwrla	LONG	v8.3
DAQmxSetWriteNextWriteIsLast	FNHx_setwrnxwrla	LONG	v8.3
DAQmxResetWriteNextWriteIsLast	FNHx_rstwrnxwrla	LONG	v8.3
DAQmxGetPhysicalChanAITermCfgs	FNHx_getphchaitcf	LONG	v8.3
DAQmxGetPhysicalChanAOTermCfgs	FNHx_getphchaotcf	LONG	v8.3
DAQmxGetPhysicalChanDIPortWidth	FNHx_getphchdipwd	LONG	v8.3
DAQmxGetPhysicalChanDISampClkSupported	FNHx_getphchdiscs	LONG	v8.3
DAQmxGetPhysicalChanDIChangeDetectSupported	FNHx_getphchdidsp	LONG	v8.3
DAQmxGetPhysicalChanDOPortWidth	FNHx_getphchdopwd	LONG	v8.3
DAQmxGetDevTerminals	FNHx_getdevtermin	LONG	v8.7

NI-DAQmx function name *)	HTBasic function name	HTBasic Ret.type	NI-DAQmx Version
DAQmxGetPhysicalChanDOSampClkSupported	FNHx_getphchdoscs	LONG	v8.3
DAQmxGetSampClkTimingResponseMode	FNHx_getsctresmod	LONG	v8.3
DAQmxSetSampClkTimingResponseMode	FNHx_setsctresmod	LONG	v8.3
DAQmxResetSampClkTimingResponseMode	FNHx_rstsctresmod	LONG	v8.3
DAQmxGetAIVoltagedBRef	FNHx_getaivolbref	LONG	v8.5
DAQmxSetAIVoltagedBRef	FNHx_setaivolbref	LONG	v8.5
DAQmxResetAIVoltagedBRef	FNHx_rstaivolbref	LONG	v8.5
DAQmxGetAISoundPressuredBRef	FNHx_getaispbyref	LONG	v8.5
DAQmxSetAISoundPressuredBRef	FNHx_setaispbyref	LONG	v8.5
DAQmxResetAISoundPressuredBRef	FNHx_rstaispbyref	LONG	v8.5
DAQmxGetAIAcceldBRef	FNHx_getaiacbyref	LONG	v8.5
DAQmxSetAIAcceldBRef	FNHx_setaiacbyref	LONG	v8.5
DAQmxResetAIAcceldBRef	FNHx_rstaiacbyref	LONG	v8.5
DAQmxGetAIADCTimingMode	FNHx_getaiadctimd	LONG	v8.5
DAQmxSetAIADCTimingMode	FNHx_setaiadctimd	LONG	v8.5
DAQmxResetAIADCTimingMode	FNHx_rstaiadctimd	LONG	v8.5
DAQmxGetCOConstrainedGenMode	FNHx_getcocongmod	LONG	v8.5
DAQmxSetCOConstrainedGenMode	FNHx_setcocongmod	LONG	v8.5
DAQmxResetCOConstrainedGenMode	FNHx_rstcocongmod	LONG	v8.5
DAQmxGetExportedSampClkDelayOffset	FNHx_getexscdloff	LONG	v8.5
DAQmxSetExportedSampClkDelayOffset	FNHx_setexscdloff	LONG	v8.5
DAQmxResetExportedSampClkDelayOffset	FNHx_rstexscdloff	LONG	v8.5
DAQmxGetDevCIMaxSize	FNHx_getdevcimasz	LONG	v8.5
DAQmxGetDevCOMaxSize	FNHx_getdevcomasz	LONG	v8.5
DAQmxGetReadOvercurrentChansExist	FNHx_getrdocuchex	LONG	v8.5
DAQmxGetReadOvercurrentChans	FNHx_getrdovcuchn	LONG	v8.5
DAQmxGetReadOpenCurrentLoopChansExist	FNHx_getrdoculcex	LONG	v8.5
DAQmxGetReadOpenCurrentLoopChans	FNHx_getrdopculch	LONG	v8.5
DAQmxGetWriteOvercurrentChansExist	FNHx_getwrovcchex	LONG	v8.5
DAQmxGetWriteOvercurrentChans	FNHx_getwrovccchns	LONG	v8.5
DAQmxGetWriteOpenCurrentLoopChansExist	FNHx_getwropcchex	LONG	v8.5
DAQmxGetWriteOpenCurrentLoopChans	FNHx_getwropccchns	LONG	v8.5
DAQmxGetWritePowerSupplyFaultChansExist	FNHx_getwrpsfchex	LONG	v8.5
DAQmxGetWritePowerSupplyFaultChans	FNHx_getwrpsfchns	LONG	v8.5
DAQmxGetAIVoltageACRMSUnits	FNHx_getaivoacrmu	LONG	v8.6
DAQmxSetAIVoltageACRMSUnits	FNHx_setaivoacrmu	LONG	v8.6
DAQmxResetAIVoltageACRMSUnits	FNHx_rstaivoacrmu	LONG	v8.6
DAQmxGetAICurrentACRMSUnits	FNHx_getaicuacrmu	LONG	v8.6
DAQmxSetAICurrentACRMSUnits	FNHx_setaicuacrmu	LONG	v8.6
DAQmxResetAICurrentACRMSUnits	FNHx_rstaicuacrmu	LONG	v8.6
DAQmxGetAOVoltageCurrentLimit	FNHx_getaoculimit	LONG	v8.6
DAQmxSetAOVoltageCurrentLimit	FNHx_setaoculimit	LONG	v8.6
DAQmxResetAOVoltageCurrentLimit	FNHx_rstaoculimit	LONG	v8.6
DAQmxGetAOFuncGenType	FNHx_getaofungtyp	LONG	v8.6
DAQmxSetAOFuncGenType	FNHx_setaofungtyp	LONG	v8.6
DAQmxResetAOFuncGenType	FNHx_rstaofungtyp	LONG	v8.6
DAQmxGetAOFuncGenFreq	FNHx_getaofungfrq	LONG	v8.6
DAQmxSetAOFuncGenFreq	FNHx_setaofungfrq	LONG	v8.6
DAQmxResetAOFuncGenFreq	FNHx_rstaofungfrq	LONG	v8.6
DAQmxGetAOFuncGenAmplitude	FNHx_getaofungamp	LONG	v8.6
DAQmxSetAOFuncGenAmplitude	FNHx_setaofungamp	LONG	v8.6
DAQmxResetAOFuncGenAmplitude	FNHx_rstaofungamp	LONG	v8.6
DAQmxGetAOFuncGenOffset	FNHx_getaofungoff	LONG	v8.6
DAQmxSetAOFuncGenOffset	FNHx_setaofungoff	LONG	v8.6
DAQmxResetAOFuncGenOffset	FNHx_rstaofungoff	LONG	v8.6
DAQmxGetAOFuncGenSquareDutyCycle	FNHx_getaofungsdc	LONG	v8.6
DAQmxSetAOFuncGenSquareDutyCycle	FNHx_setaofungsdc	LONG	v8.6
DAQmxResetAOFuncGenSquareDutyCycle	FNHx_rstaofungsdc	LONG	v8.6
DAQmxGetAOFuncGenModulationType	FNHx_getaofungmtp	LONG	v8.6
DAQmxSetAOFuncGenModulationType	FNHx_setaofungmtp	LONG	v8.6
DAQmxResetAOFuncGenModulationType	FNHx_rstaofungmtp	LONG	v8.6
DAQmxGetAOFuncGenFMDeviation	FNHx_getaofungfmd	LONG	v8.6
DAQmxSetAOFuncGenFMDeviation	FNHx_setaofungfmd	LONG	v8.6
DAQmxResetAOFuncGenFMDeviation	FNHx_rstaofungfmd	LONG	v8.6
DAQmxGetSampTimingEngine	FNHx_getsamtimg	LONG	v8.6
DAQmxSetSampTimingEngine	FNHx_setsamtimg	LONG	v8.6
DAQmxResetSampTimingEngine	FNHx_rstsamtimg	LONG	v8.6
DAQmxGetPhysicalChanAOManualControlEnable	FNHx_getphchaomce	LONG	v8.6
DAQmxSetPhysicalChanAOManualControlEnable	FNHx_setphchaomce	LONG	v8.6

NI-DAQmx function name *)	HTBasic function name	HTBasic Ret.type	NI-DAQmx Version
DAQmxResetPhysicalChanAOManualControlEnable	FNHx_rstphchaomce	LONG	v8.6
DAQmxGetPhysicalChanAOManualControlAmplitude	FNHx_getpcaomcamp	LONG	v8.6
DAQmxGetPhysicalChanAOManualControlFreq	FNHx_getpcaomcfrq	LONG	v8.6
DAQmxPerformBridgeOffsetNullingCalEx	FNHx_perfbrofncex	LONG	v8.7
DAQmxPerformStrainShuntCal	FNHx_perfstrshcal	LONG	v8.7
DAQmxPerformBridgeShuntCal	FNHx_perfbrishcal	LONG	v8.7
DAQmxSetCOPulseFreq	FNHx_setcopulfreq	LONG	
DAQmxSetCOPulseDutyCyc	FNHx_setcoputdcyc	LONG	
DAQmxSetBufferAttribute	FNHx_setbufattrib	LONG	
DAQmxSetTrigAttribute	FNHx_settriattrib	LONG	
<b>16. Callback services</b>			
DAQmxRegisterEveryNSamplesEvent	FNHx_regevnsmevnt	LONG	
DAQmxRegisterSignalEvent	FNHx_regsignevent	LONG	
DAQmxRegisterDoneEvent	FNHx_regdoneevent	LONG	
<b>17. Watchdog</b>			
DAQmxCreateWatchdogTimerTask	FNHx_crwdogtimtsk	LONG	
DAQmxControlWatchdogTask	FNHx_ctrlwdogtsk	LONG	
<b>18. Storage</b>			
DAQmxSaveTask	FNHx_savetask	LONG	
DAQmxDeleteSavedTask	FNHx_delsavedtask	LONG	
<b>19. Special calibration functions</b>			
DAQmxGetAIChanCalPolyForwardCoeff	FNHx_getaichcpfco	LONG	
DAQmxSetAIChanCalPolyForwardCoeff	FNHx_setaichcpfco	LONG	
DAQmxResetAIChanCalPolyForwardCoeff	FNHx_rstaichcpfco	LONG	
DAQmxGetAIDevScalingCoeff	FNHx_getaidevscco	LONG	

## Determining HTBasic function names for NI-DAQmx functions

In case you need to know the corresponding HTBasic function name for a given NI-DAQmx function you might want to use the HTBasic support function `FNHx_load_libs`:

```
INTEGER Debug, Dummy
Dummy=FNHx_load_libs(Debug, "DAQmxCreateTEDSAIPosRVDTChan")
```

### HTBasic Output:

```
NI-DAQmx function name: DAQmxCreateTEDSAIPosRVDTChan
HTBasic function name: FNHx_crtaiprvdtch
```

```
! this lists all available NI-DAQmx functions and their HTBasic function names
Dummy=FNHx_load_libs(Debug, "**")
```

**Info:** `FNHx_load_libs` is part of the HTBasic library file "daqmx\_htb.lib".

## Special HTBasic Support Functions

HTBasic function name	Description	HTBasic Return type
SUB Hx_manage_subs	Load / Init or Unload/De-init HTBasic lib	VOID
FNHx_getdaqmxvers\$	Retrieve DAQmx version numbers	STRING
FNHx_getdllvers\$	Retrieve DAQmx wrapper DLL version number	STRING
FNHx_getdevname\$	Extract device name (strips leading "/")	STRING
FNHx_gettnamwpref	Special function	LONG
SUB Hx_show_errstat	Show error status information	VOID
FNHx_get_constant	Return constant from nidaqmx.h	INT / LONG / REAL
FNHx_get_constant\$	Return constant from nidaqmx.h	STRING
SUB Hx_enter_kbd	Wait for key stroke	VOID
FNHx_get_dev_info	Get device and system information	LONG
FNHx_lib_version\$	Get HTBasic library version information	STRING
FNHx_numofdevices	Get numbers of NI devices defined in MAX	LONG
FNHx_getsfpath\$	Get path to special windows directories	STRING

## Global variables for HTBasic

The following global variables are used for callback services.

HTBasic variable name	Description	HTBasic Return type
g_Evnt_taskhandle	Taskhandle for which the callback occurred	LONG
g_Evnt_type	Event type for callback	LONG

g_Evnt_status	Special status information for event	LONG
g_Callback_locked	Callback Locking indicator	LONG
g_Callback_pending	Indicates that another callback is pending	LONG
g_nSamples	Number of read samples for callback	LONG

**Please note:** You must declare and initialize all variables in HTBasic before you call DLL functions, otherwise the system may become instable.

### Variable mapping

NI-DAQmx variable	HTBasic variable
short	INTEGER
long	LONG
uInt8	INTEGER
uInt16	INTEGER
uInt32	LONG
uInt64	N.A.
bool32	LONG
char	STRING
float64	REAL

### Use the National Instruments NI-DAQmx Documentation

It is important that you understand how the original functions of the National Instruments DAQmx library work, so always look into the documentation "NI-DAQmx C Reference Help" for a more detailed description of the parameters used and their meaning. You can find many example programs which can be easily ported to HTBasic in the following directory:

```
%Program Files%\National Instruments\NI-DAQ\Examples\DAQmx ANSI C
```

The function reference always shows the HTBasic function name and the original NI-DAQmx function name, see the example below:

<b>FNHx_loadtask</b>	NI-DAQmx: <b>DAQmxLoadTask</b>
----------------------	--------------------------------

**Purpose:** Loads an existing task

### Return variables of DLL functions

Due to the nature of DLL handling it is very important that you always use a return variable for the function calls, otherwise HTBasic may crash. Example:

#### Correct:

```
LONG Result
LONG Taskid2
Result=FNHx_stoptask((Taskid2))
!
IF Result=0 THEN
    !... subsequent commands...
END IF
```

#### Incorrect:

```
LONG Taskid2
!
IF FNHx_stoptask((Taskid2))=0 THEN ! this may crash HTBasic !!!
    !... subsequent commands...
END IF
```

The (incorrect) construction shown above should never be used, you always must use a return value for DLL functions, otherwise HTBasic may crash or become unstable.



# Function declarations for all implemented functions

## Task Configuration and Control

**FNHx\_loadtask**

NI-DAQmx: **DAQmxLoadTask**

**Purpose:** Loads an existing task

**HTBasic Definition:**

```
LONG Result=FNHx_loadtask((Taskname$),LONG Taskid)
```

Input:

```
Taskname$           Name of task to load
```

Return value: 0 = Function was successful  
                  **Taskid** now contains Task ID of loaded task  
                  <>0 = Error occurred

**Example:**

```
LONG Result
LONG Taskid
Result=FNHx_loadtask("VoltageTask",Taskid)
!
IF Result=0 THEN
!... subsequent commands...
END IF
```

**FNHx\_createtask**

NI-DAQmx: **DAQmxLoadTask**

**Purpose:** Creates a new task

**HTBasic Definition:**

```
LONG Result=FNHx_createtask(LONG Taskid,(Taskname$))
```

Input:

```
Taskname$           Name of task to be created (or NULL string "")
```

Return value: 0 = Function was successful  
                  **Taskid** now contains Task ID of newly created task  
                  <>0 = Error occurred

**Example:**

```
INTEGER Verbose
LONG Result
LONG Analogin_task
Result=FNHx_createtask(Analogin_task,"AnalogInTask")
!
IF Result=0 THEN
!... subsequent commands...
END IF
```

**FNHx\_starttask**

NI-DAQmx: **DAQmxStartTask**

**Purpose:** Starts a task

**HTBasic Definition:**

```
LONG Result=FNHx_starttask((LONG Taskid))
```

Input:

```
Taskid              Task ID of task to start
```

Return value: 0 = Function was successful  
                  <>0 = Error occurred

**Example:**

```
LONG Result
LONG Taskid
Result=FNHx_starttask((Taskid))
!
IF Result=0 THEN
!... subsequent commands...
END IF
```

**Purpose:** Stops a task

**HTBasic Definition:**

LONG Result=FNHx\_stoptask((LONG Taskid))

Input:

Taskid                    Task ID of task to stop

Return value:    0 = Function was successful  
                  <>0 = Error occurred

**Example:**

```
LONG Result
LONG Taskhandle
Result=FNHx_stoptask((Taskhandle))
!
IF Result=0 THEN
!... subsequent commands...
END IF
```

**Purpose:** Clears a task

**HTBasic Definition:**

LONG Result=FNHx\_cleartask((LONG Taskid))

Input:

Taskid                    Task ID of task to be deleted

Return value:    0 = Function was successful  
                  <>0 = Error occurred

**Example:**

```
LONG Result
LONG Taskid
Result=FNHx_cleartask((Taskid))
!
IF Result=0 THEN
!... subsequent commands...
END IF
```

**Purpose:** Waits until a task is done

**HTBasic Definition:**

LONG Result=FNHx\_waituntaskdn((LONG Taskid), REAL Timetowait)

Input:

Taskid                    Task ID of task to be checked

Timetowait                Wait time, until the task should be finished until a timeout occurs

Return value:    0 = Function was successful  
                  <>0 = Error occurred

**Example:**

```
LONG Result
LONG Taskid
REAL Timetowait
Timetowait=5                ! 5 wait time
Result=FNHx_waituntaskdn((Taskid),Timetowait)
!
IF Result=0 THEN
!... subsequent commands...
END IF
```

**Purpose:** Checks if a task is done

**HTBasic Definition:**

LONG Result=FNHx\_istaskdone((LONG Taskid), LONG Istaskdone)

Input:

Taskid Task ID of task to be checked

Return value: 0 = Function was successful  
**Istaskdone** now contains 0 (Task done) or 1 (Task still running)  
 <>0 = Error occurred

**Example:**

```
LONG Istaskdone
LONG Taskid
Result=FNHx_istaskdone((Taskid),Istaskdone)
!
IF Result=0 THEN
  PRINT "Task Status of";Taskid;";";Istaskdone
  !... subsequent commands...
END IF
```

**Purpose:** Control a task

**HTBasic Definition:**

LONG Result=FNHx\_taskcontrol((LONG Taskid), (LONG Action))

Input:

Taskid Task ID of task to be controlled

Action Control code (e.g. "DAQmx\_Val\_Task\_Start" or "DAQmx\_Val\_Task\_Abort")

Return value: 0 = Function was successful  
 <>0 = Error occurred

**Example:**

```
LONG Taskid, Action
INTEGER Iresult, Debug
Action=FNHx_get_constant("DAQmx_Val_Task_Abort",Iresult,Debug) ! get LONG constant
Result=FNHx_taskcontrol((Taskid),(Action))
!
IF Result=0 THEN
  !... subsequent commands...
END IF
```

**Purpose:** Returns the *n*th channel of a task

**HTBasic Definition:**

LONG Result=FNHx\_getnthtaskch((LONG Taskid), (LONG Index), Buf\$, (LONG Bufsize))

Input:

Taskid The task used in this function

Index The Nth channel you want to return. The index starts at 1

Bufsize The size, in bytes, of buffer. If you pass 0, this function returns the number of bytes needed to allocate.

Output:

Buf\$ The Nth channel in the index. If you pass NULL, this function returns the number of bytes needed to allocate.

Return value: 0 = Function was successful  
 <>0 = Error occurred

**Example:**

```
LONG Taskid, Index, Bufsize
DIM Buf$[1024]
ASize=SIZE(Buf$) ! Size of array
Result=FNHx_getnthtaskch((Taskid),(Index),Buf$, (Bufsize))
```

## Channel Configuration / Creation

**FNHx\_craivoltchan**

NI-DAQmx: **DAQmxCreateAIVoltageChan**

**Purpose:** Defines an analog Input Voltage channel

### HTBasic Definition:

```
LONG Result=FNHx_craivoltchan((LONG Taskid), Phys_chan$, My_chn$, (LONG Termcfg),  
                              (REAL Minval), (REAL Maxval), (LONG Units), Custscalname$)
```

#### Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Termcfg	The input terminal configuration for the channel
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the voltage measurements
Custscalname\$	The name of a custom scale to apply to the channel

**Return value:** 0 = Function was successful  
<>0 = Error occurred

### Example:

```
INTEGER Iresult1, Iresult2  
LONG Result, Taskid, Termcfg, Units  
REAL V_min, V_max  
DIM Tmp$[255], Devchn$[255]  
!  
Termcfg=FNHx_get_constant("DAQmx_Val_Cfg_Default",Iresult1) ! constant 1  
Units=FNHx_get_constant("DAQmx_Val_Volts",Iresult2) ! constant 2  
V_min=-10.0  
V_max=10.0  
Tmp$=""  
Devchn$="Dev1/ai0"  
Result=FNHx_craivoltchan((Taskid),Devchn$,Tmp$, (Termcfg), (V_min), (V_max), (Units), Tmp$)
```

**FNHx\_craicurrchan**

NI-DAQmx: **DAQmxCreateAICurrentChan**

**Purpose:** Defines an analog Current Input channel

### HTBasic Definition:

```
LONG Result=FNHx_craicurrchan((LONG Taskid), Phys_chan$, My_chn$, (LONG Termcfg),  
                              (REAL Minval), (REAL Maxval), (LONG Units),  
                              (REAL Shuntresloc), (REAL Extshuntresval), Custscalname$)
```

#### Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Termcfg	The input terminal configuration for the channel
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the voltage measurements
Shuntresloc	The location of the shunt resistor
Extshuntresval	The value, in ohms, of an external shunt resistor
Custscalname\$	The name of a custom scale to apply to the channel

**Return value:** 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Defines analog temperature input channel

**HTBasic Definition:**

```
LONG Result=FNHx_craithcpchan((LONG Taskid), Phys_chan$, My_chn$, (LONG Termcfg),
                              (REAL Minval), (REAL Maxval), (LONG Units),
                              (REAL Tcouptype), (REAL Cjcsrc), (REAL Cjcval), Cjcchan$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Termcfg	The input terminal configuration for the channel
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the voltage measurements
Tcouptype	The type of thermocouple connected to the channel
Cjcsrc	The source of cold junction compensation
Cjcval	The temperature of the cold junction of the thermocouple if you set <b>Cjcsrc</b> to DAQmx_Val_ConstVal
Cjcchan\$	The channel that acquires the temperature of the thermocouple cold-junction if you set <b>Cjcsrc</b> to DAQmx_Val_Chan

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates a channel that uses a RTD to measure temperature

**HTBasic Definition:**

```
LONG Result=FNHx_crairtchan((LONG Taskid), Phys_chan$, My_chn$, (LONG Termcfg),
                             (REAL Minval), (REAL Maxval), (LONG Units),
                             (REAL Rtdtype), (LONG Resiscfg), (LONG Curexcsrc),
                             (REAL Curexcval), (REAL R0))
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Termcfg	The input terminal configuration for the channel
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the measurements
Rtdtype	The type of RTD connected to the channel
Resiscfg	The configuration for resistance measurements
Curexcsrc	The source of excitation
Curexcval	The amount of excitation, in amperes, that the sensor requires
R0	The sensor resistance in ohms at 0 deg C for the Callendar-Van Dusen equation

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) that use a thermistor to measure temperature (thermistor requires voltage excitation)

**HTBasic Definition:**

```
LONG Result=FNHx_craithmchiex((LONG Taskid), Phys_chan$, My_chn$, (REAL Minval),
                              (REAL Maxval), (LONG Units), (LONG Resconfig),
                              (LONG Curexcsrc), (REAL Curexcval), (REAL A), (REAL B),
                              (REAL C))
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the measurements
Resconfig	The configuration for resistance measurements (2-, 3- or 4-wire)
Curexcsrc	The source of excitation
Curexcval	The amount of excitation, in amperes, that the sensor requires
A,B,C	The A, B and C constants from the Steinhart-Hart thermistor equation

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) that use a thermistor to measure temperature (thermistor requires voltage excitation)

**HTBasic Definition:**

```
LONG Result=FNHx_craithmchvex((LONG Taskid), Phys_chan$, My_chn$, (REAL Minval),
                              (REAL Maxval), (LONG Units), (LONG Resconfig),
                              (LONG Volexcsrc), (REAL Volexcval), (REAL A), (REAL B),
                              (REAL C), (REAL R1))
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the measurements
Resconfig	The configuration for resistance measurements (2-, 3- or 4-wire)
Volexcsrc	The source of excitation
Volexcval	The amount of excitation, in volts, that the sensor requires
A,B,C	The A, B and C constants from the Steinhart-Hart thermistor equation
R1	The value, in ohms, of the reference resistor

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) that use a frequency-to-voltage converter to measure frequency

**HTBasic Definition:**

```
LONG Result=FNHx_craifrvochan((LONG Taskid), Phys_chan$, My_chn$, (LONG Termcfg),
                              (REAL Minval), (REAL Maxval), (LONG Units),
                              (REAL Tholdlev), (REAL Hysteresis), Custscalename$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Termcfg	The input terminal configuration for the channel
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the measurements
Tholdlev	The voltage level at which to recognize waveform repetitions
Hysteresis	Specifies in volts a window below thresholdLevel
Custscalename\$	The name of a custom scale to apply to the channel

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Create channel(s) to measure resistance

**HTBasic Definition:**

```
LONG Result=FNHx_crairesichan((LONG Taskid), Phys_chan$, My_chn$, (LONG Termcfg),
                              (REAL Minval), (REAL Maxval), (LONG Units),
                              (REAL Resconfig), (REAL Curexcsrc), (REAL Curexcval),
                              Custscalename$)
```

**Input:**

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Termcfg	The input terminal configuration for the channel
Minval	The minimum value, in units, that you expect to measure
Maxval	The maximum value, in units, that you expect to measure
Units	The units to use to return the measurements
Resconfig	The configuration for resistance measurements
Curexcsrc	The source of excitation
Curexcval	The amount of excitation, in amperes, that the sensor requires
Custscalename\$	The name of a custom scale to apply to the channel

**Return value:** 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) to measure strain

**HTBasic Definition:**

```
LONG Result=FNHx_craistrgachn((LONG Taskid), Phys_chan$, My_chn$, (REAL Minval),
                              (REAL Maxval), (LONG Units), (LONG Strainconfig),
                              (LONG Volexcsrc), (REAL Volexcval), (REAL Gagefactor),
                              (REAL Inbrvol), (REAL Nomgageres), (REAL Poisrat),
                              (REAL Leadwireres), Custscalename$)
```

**Input:**

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Minval	The minimum value, in units, that you expect to measure
Maxval	The maximum value, in units, that you expect to measure
Units	The units to use to return the measurements
Strainconfig	The strain gage bridge configuration
Volexcsrc	The source of excitation
Volexcval	The amount of excitation, in volts, that the sensor requires
Gagefactor	The sensitivity of the strain gages (see DAQmx documentation)
Inbrvol	The bridge output voltage in the unloaded condition
Nomgageres	The resistance, in ohms, of the gages in an unstrained position
Poisrat	The ratio of lateral strain to axial strain in the material in which you measure strain
Leadwireres	The amount, in ohms, of resistance in the lead wires
Custscalename\$	The name of a custom scale to apply to the channel

**Return value:** 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) to measure voltage for custom sensors with excitation

**HTBasic Definition:**

```
LONG Result=FNHx_craivochwexc((LONG Taskid), Phys_chan$, My_chn$, (LONG Termconfig),
                              (REAL Minval), (REAL Maxval), (LONG Units), (LONG Bricfg),
                              (LONG Volexcsrc), (REAL Volexcval), (LONG Useexcforscal),
                              Custscalename$)
```

**Input:**

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Termconfig	The input terminal configuration for the channel
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the measurements

Bricfg                    The type of Wheatstone bridge the sensor is  
 Volexcsrc                The source of excitation  
 Volexcval                The amount of excitation, in volts, that the sensor requires  
 Useexcforscal            Specifies whether NI-DAQmx divides the measurement by the excitation  
 Custscalename\$         The name of a custom scale to apply to the channel

Return value:    0 = Function was successful  
                   <>0 = Error occurred

<b>FNHx_craitmpbschn</b>	NI-DAQmx: <b>DAQmxCreateAITempBuiltInSensorChan</b>
--------------------------	---

**Purpose:** Creates channel(s) to measure temperature with a built-in sensor of a terminal block or device

**HTBasic Definition:**

**LONG** Result=FNHx\_craitmpbschn(**LONG** Taskid, Phys\_chan\$, My\_chn\$, (**LONG** Units))

Input:

Taskid                    The task to which to add the channels that this function creates  
 Phys\_chan\$                The names of the physical channels to use to create virtual channels  
 My\_chn\$                    The name(s) to assign to the created virtual channel(s)  
 Units                      The units to use to return the measurements

Return value:    0 = Function was successful  
                   <>0 = Error occurred

<b>FNHx_craiaaccelchn</b>	NI-DAQmx: <b>DAQmxCreateAIAccelChan</b>
---------------------------	---

**Purpose:** Creates channel(s) that use an accelerometer to measure acceleration

**HTBasic Definition:**

**LONG** Result=FNHx\_craiaaccelchn(**LONG** Taskid), Phys\_chan\$, My\_chn\$, (**LONG** Termconfig),  
 (**REAL** Minval), (**REAL** Maxval), (**LONG** Units), (**LONG** Sensit),  
 (**LONG** Curexcsrc), (**REAL** Curexcval), Custscalename\$

Input:

Taskid                    The task to which to add the channels that this function creates  
 Phys\_chan\$                The names of the physical channels to use to create virtual channels  
 My\_chn\$                    The name(s) to assign to the created virtual channel(s)  
 Termconfig                The input terminal configuration for the channel  
 Minval                    The minimum value, in **Units**, that you expect to measure  
 Maxval                    The maximum value, in **Units**, that you expect to measure  
 Units                      The units to use to return the measurements  
 Sensit                    The sensitivity of the sensor  
 Curexcsrc                The source of excitation  
 Curexcval                The amount of excitation, in ampere, that the sensor requires  
 Custscalename\$         The name of a custom scale to apply to the channel

Return value:    0 = Function was successful  
                   <>0 = Error occurred

<b>FNHx_craimicrchan</b>	NI-DAQmx: <b>DAQmxCreateAIMicrophoneChan</b>
--------------------------	--

**Purpose:** Creates channel(s) that use a microphone to measure sound pressure

**HTBasic Definition:**

**LONG** Result=FNHx\_craimicrchan(**LONG** Taskid), Phys\_chan\$, My\_chn\$, (**LONG** Termconfig),  
 (**LONG** Units), (**REAL** Micsens), (**REAL** Maxsoundplevel),  
 (**LONG** Curexcsrc), (**REAL** Curexcval), Custscalename\$

Input:

Taskid                    The task to which to add the channels that this function creates  
 Phys\_chan\$                The names of the physical channels to use to create virtual channels  
 My\_chn\$                    The name(s) to assign to the created virtual channel(s)  
 Termconfig                The input terminal configuration for the channel  
 Units                      The units to use to return the measurements  
 Micsens                    The sensitivity of the microphone (in millivolts per pascal)  
 Maxsoundplevel         The maximum instantaneous sound pressure level you expect to measure.  
                               This value is in decibels, referenced to 20 micropascals  
 Curexcsrc                The source of excitation  
 Curexcval                The amount of excitation, in ampere, that the sensor requires  
 Custscalename\$         The name of a custom scale to apply to the channel

Return value:    0 = Function was successful  
                   <>0 = Error occurred



**Purpose:** Creates channel(s) that use an LVDT to measure linear position

**HTBasic Definition:**

```
LONG Result=FNHx_craipslvdтч((LONG Taskid), Phys_chan$, My_chn$, (REAL Minval),
                             (REAL Maxval),(LONG Units), (REAL Sens), (LONG Sensunits),
                             (LONG Volexcsrc), (REAL Volexcval), (REAL Volexcfrq),
                             (LONG Acexcwmode), Custscalename$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the measurements
Sens	The sensitivity of the sensor
Sensunits	The units of <b>Sens</b>
Volexcsrc	The source of excitation
Volexcval	The amount of excitation, in volts, that the sensor requires
Volexcfrq	The excitation frequency, in hertz, that the sensor requires
Acexcwmode	The number of leads on the sensor
Custscalename\$	The name of a custom scale to apply to the channel

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) that use an RVDT to measure angular position

**HTBasic Definition:**

```
LONG Result=FNHx_craipsrvdтч((LONG Taskid), Phys_chan$, My_chn$, (REAL Minval),
                             (REAL Maxval),(LONG Units), (REAL Sens), (LONG Sensunits),
                             (LONG Volexcsrc), (REAL Volexcval), (REAL Volexcfrq),
                             (LONG Acexcwmode), Custscalename$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the measurements
Sens	The sensitivity of the sensor
Sensunits	The units of <b>Sens</b>
Volexcsrc	The source of excitation
Volexcval	The amount of excitation, in volts, that the sensor requires
Volexcfrq	The excitation frequency, in hertz, that the sensor requires
Acexcwmode	The number of leads on the sensor
Custscalename\$	The name of a custom scale to apply to the channel

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) to measure voltage

**HTBasic Definition:**

```
LONG Result=FNHx_crtaivoltchn((LONG Taskid), Phys_chan$, My_chn$, (LONG Termcfg),
                              (REAL Minval), (REAL Maxval), (LONG Units),
                              Custscalename$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Termcfg	The input terminal configuration for the channel
Minval	The minimum value, in units, that you expect to measure
Maxval	The maximum value, in units, that you expect to measure
Units	The units to use to return the measurements
Custscalename\$	The name of a custom scale to apply to the channel

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) to measure current

**HTBasic Definition:**

```
LONG Result=FNHx_crtaicurrchn((LONG Taskid), Phys_chan$, My_chn$, (LONG Termcfg),
                              (REAL Minval), (REAL Maxval), (LONG Units),
                              (LONG Shuntresloc), (LONG Extshuntresval), Custscalename$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Termcfg	The input terminal configuration for the channel
Minval	The minimum value, in units, that you expect to measure
Maxval	The maximum value, in units, that you expect to measure
Units	The units to use to return the measurements
Shuntresloc	The location of the shunt resistor
Extshuntresval	The value, in ohms, of an external shunt resistor
Custscalename\$	The name of a custom scale to apply to the channel

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) that use a thermocouple to measure temperature

**HTBasic Definition:**

```
LONG Result=FNHx_crtaithcpchn((LONG Taskid), Phys_chan$, My_chn$, (LONG Termcfg),
                              (REAL Minval), (REAL Maxval), (LONG Units),
                              (LONG Cjcsource), (REAL Cjcval), Cjcchannel$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Termcfg	The input terminal configuration for the channel
Minval	The minimum value, in units, that you expect to measure
Maxval	The maximum value, in units, that you expect to measure
Units	The units to use to return the measurements
Cjcsource	The source of cold junction compensation
Cjcval	The temperature of the cold junction of the thermocouple if you set <b>Cjcsource</b> to "DAQmx_Val_ConstVal"
Cjcchannel\$	The channel that acquires the temperature of the thermocouple cold-junction if you set <b>Cjcsource</b> to "DAQmx_Val_Chan"

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) that use an RTD to measure temperature

**HTBasic Definition:**

```
LONG Result=FNHx_crtairtdchn((LONG Taskid), Phys_chan$, My_chn$, (REAL Minval),
                              (REAL Maxval), (LONG Units), (LONG Resconfig),
                              (LONG Curexcsrc), (REAL Curexcval))
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the measurements
Resconfig	The configuration for resistance measurements
Curexcsrc	The source of excitation
Curexcval	The amount of excitation, in amperes, that the sensor requires

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) that use a thermistor to measure temperature

**HTBasic Definition:**

```
LONG Result=FNHx_crtaiithchiex((LONG Taskid), Phys_chan$, My_chn$, (REAL Minval),
                                (REAL Maxval), (LONG Units), (LONG Resconfig),
                                (LONG Curexcsrc), (REAL Curexcval))
```

**Input:**

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the measurements
Resconfig	The configuration for resistance measurements (2-, 3- or 4-wire)
Curexcsrc	The source of excitation
Curexcval	The amount of excitation, in amperes, that the sensor requires

**Return value:** 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) that use a thermistor to measure temperature

**HTBasic Definition:**

```
LONG Result=FNHx_crtaiithchvex((LONG Taskid), Phys_chan$, My_chn$, (REAL Minval),
                                (REAL Maxval), (LONG Units), (LONG Resconfig),
                                (LONG Volexcsrc), (REAL Volexcval), (REAL R1))
```

**Input:**

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the measurements
Resconfig	The configuration for resistance measurements (2-, 3- or 4-wire)
Volexcsrc	The source of excitation
Volexcval	The amount of excitation, in volts, that the sensor requires
R1	The value, in ohms, of the reference resistor

**Return value:** 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) to measure resistance

**HTBasic Definition:**

```
LONG Result=FNHx_crtairesichn((LONG Taskid), Phys_chan$, My_chn$, (REAL Minval),
                                (REAL Maxval), (LONG Units), (REAL Resconfig),
                                (REAL Curexcsrc), (REAL Curexcval), Custscalename$)
```

**Input:**

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Minval	The minimum value, in units, that you expect to measure
Maxval	The maximum value, in units, that you expect to measure
Units	The units to use to return the measurements
Resconfig	The configuration for resistance measurements
Curexcsrc	The source of excitation
Curexcval	The amount of excitation, in amperes, that the sensor requires
Custscalename\$	The name of a custom scale to apply to the channel

**Return value:** 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) to measure strain

**HTBasic Definition:**

```
LONG Result=FNHx_crtaistrgchn((LONG Taskid), Physchn$, Mychn$, (REAL Minv), (REAL Maxv),
                              (LONG Units), (LONG Volexcsrc), (REAL Volexcval),
                              (REAL Inibrivol), (REAL Leadwireres), Custscalename$)
```

**Input:**

Taskid	The task to which to add the channels that this function creates
Physchn\$	The names of the physical channels to use to create virtual channels
Mychn\$	The name(s) to assign to the created virtual channel(s)
Minv	The minimum value, in units, that you expect to measure
Maxv	The maximum value, in units, that you expect to measure
Units	The units to use to return the measurements
Volexcsrc	The source of excitation
Volexcval	The amount of excitation, in volts, that the sensor requires
Inibrivol	The bridge output voltage in the unloaded condition
Leadwireres	The amount, in ohms, of resistance in the lead wires
Custscalename\$	The name of a custom scale to apply to the channel

**Return value:** 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) to measure voltage (for custom sensors that require excitation)

**HTBasic Definition:**

```
LONG Result=FNHx_crtaivochwex((LONG Taskid), Phys_chan$, My_chn$, (LONG Termconfig),
                              (REAL Minval), (REAL Maxval), (LONG Units),
                              (LONG Volexcsrc), (REAL Volexcval), Custscalename$)
```

**Input:**

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Termconfig	The input terminal configuration for the channel
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the measurements
Volexcsrc	The source of excitation
Volexcval	The amount of excitation, in volts, that the sensor requires
Custscalename\$	The name of a custom scale to apply to the channel

**Return value:** 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) that use an accelerometer to measure acceleration

**HTBasic Definition:**

```
LONG Result=FNHx_craitacclchn((LONG Taskid), Phys_chan$, My_chn$, (LONG Termconfig),
                              (REAL Minval), (REAL Maxval), (LONG Units),
                              (LONG Curexcsrc), (REAL Curexcval), Custscalename$)
```

**Input:**

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Termconfig	The input terminal configuration for the channel
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the measurements
Curexcsrc	The source of excitation
Curexcval	The amount of excitation, in ampere, that the sensor requires
Custscalename\$	The name of a custom scale to apply to the channel

**Return value:** 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) that use a microphone to measure sound pressure

**HTBasic Definition:**

```
LONG Result=FNHx_crtaimicchan((LONG Taskid), Phys_chan$, My_chn$, (LONG Termconfig),
                              (LONG Units), (REAL Maxsoundplevel), (LONG Curexcsrc),
                              (REAL Curexcval), Custscalename$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Termconfig	The input terminal configuration for the channel
Units	The units to use to return the measurements
Maxsoundplevel	The maximum instantaneous sound pressure level you expect to measure. This value is in decibels, referenced to 20 micropascals
Curexcsrc	The source of excitation
Curexcval	The amount of excitation, in ampere, that the sensor requires
Custscalename\$	The name of a custom scale to apply to the channel

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) that use an LVDT to measure linear position

**HTBasic Definition:**

```
LONG Result=FNHx_crtaiplvdtch((LONG Taskid), Physchn$, Mychn$, (REAL Minv), (REAL Maxv),
                              (LONG Units), (LONG Volexcsrc), (REAL Volexcval),
                              (REAL Volexcfrq), (LONG Acexcwmode), Custscalename$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Physchn\$	The names of the physical channels to use to create virtual channels
Mychn\$	The name(s) to assign to the created virtual channel(s)
Minv	The minimum value, in <b>Units</b> , that you expect to measure
Maxv	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the measurements
Volexcsrc	The source of excitation
Volexcval	The amount of excitation, in volts, that the sensor requires
Volexcfrq	The excitation frequency, in hertz, that the sensor requires
Acexcwmode	The number of leads on the sensor
Custscalename\$	The name of a custom scale to apply to the channel

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) that use an RVDT to measure angular position

**HTBasic Definition:**

```
LONG Result=FNHx_crtairvdtch((LONG Taskid), Physchn$, Mychn$, (REAL Minv), (REAL Maxv),
                              (LONG Units), (LONG Volexcsrc), (REAL Volexcval),
                              (REAL Volexcfrq), (LONG Acexcwmode), Custscalename$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Physchn\$	The names of the physical channels to use to create virtual channels
Mychn\$	The name(s) to assign to the created virtual channel(s)
Minv	The minimum value, in <b>Units</b> , that you expect to measure
Maxv	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the measurements
Volexcsrc	The source of excitation
Volexcval	The amount of excitation, in volts, that the sensor requires
Volexcfrq	The excitation frequency, in hertz, that the sensor requires
Acexcwmode	The number of leads on the sensor
Custscalename\$	The name of a custom scale to apply to the channel

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) to generate voltage

**HTBasic-Definition:**

```
LONG Result=FNIHx_craovoltchan((LONG Taskid), Phys_chan$, My_chn$, (REAL Minval),
                               (REAL Maxval), (LONG Units), Custscalename$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Minval	The minimum value, in units, that you expect to measure
Maxval	The maximum value, in units, that you expect to measure
Units	The units to use to return the measurements
Custscalename\$	The name of a custom scale to apply to the channel

Return value: 0 = Function was successful  
<>0 = Error occurred

**Example:**

```
INTEGER Iresult1
LONG Result, Taskid, Units
REAL V_min, V_max
DIM Tmp$[255], Devchn$[255]
!
Units=FNIHx_get_constant("DAQmx_Val_Volts",Iresult1)           ! constant 1
V_min=-10.0
V_max=10.0
Tmp$=""
Devchn$="Dev1/ao0"
Result=FNIHx_craovoltchan((Taskid),Devchn$,Tmp$, (V_min), (V_max), (Units), Tmp$)
```

**Purpose:** Creates channel(s) to generate current

**HTBasic-Definition:**

```
LONG Result=FNIHx_craocurrchan((LONG Taskid), Phys_chan$, My_chn$, (REAL Minval),
                                (REAL Maxval), (LONG Units), Custscalename$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Minval	The minimum value, in units, that you expect to measure
Maxval	The maximum value, in units, that you expect to measure
Units	The units to use to return the measurements
Custscalename\$	The name of a custom scale to apply to the channel

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) to measure digital signals

**HTBasic Definition:**

`LONG Result=FNHx_createdichan((LONG Taskid), Lines$, My_chn$, (LONG Linegrouping))`

Input:

Taskid	The task to which to add the channels that this function creates
Lines\$	The names of the digital lines used to create a virtual channel
My_chn\$	The name(s) to assign to the created virtual channel(s)
Linegrouping	Specifies whether to group digital lines into one or more virtual channels

Return value: 0 = Function was successful  
<>0 = Error occurred

**Example:**

```
INTEGER Iresult1
LONG Result, Taskid, Linegrouping
DIM Lines$(255), Tmp$(255)
!
Linegrouping=FNHx_get_constant("DAQmx_Val_ChanForAllLines",Iresult1)
Lines$="Dev1/port0/line0:3"
Tmp$=""
Result=FNHx_createdichan((Taskid),Lines$,Tmp$, (Linegrouping))
```

**Purpose:** Creates channel(s) to generate digital signals

**HTBasic Definition:**

`LONG Result=FNHx_createdochan((LONG Taskid), Lines$, My_chn$, (LONG Linegrouping))`

Input:

Taskid	The task to which to add the channels that this function creates
Lines\$	The names of the digital lines used to create a virtual channel
My_chn\$	The name(s) to assign to the created virtual channel(s)
Linegrouping	Specifies whether to group digital lines into one or more virtual channels

Return value: 0 = Function was successful  
<>0 = Error occurred

**Example:**

```
INTEGER Iresult1
LONG Result, Taskid, Linegrouping
DIM Lines$(255), Tmp$(255)
!
Linegrouping=FNHx_get_constant("DAQmx_Val_ChanForAllLines",Iresult1)
Lines$="Dev1/port0/line4:7"
Tmp$=""
Result=FNHx_createdochan((Taskid),Lines$,Tmp$, (Linegrouping))
```

**Purpose:** Creates a channel to measure the frequency of a digital signal and adds the channel to the task

**HTBasic Definition:**

```
LONG Result=FNHx_crcifreqchan((LONG Taskid), Counter$, My_chn$, (REAL Minval),
                              (REAL Maxval), (LONG Units), (LONG Edge),
                              (LONG Measmeth), (REAL Meastime), (LONG Divisor), Csn$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Counter\$	The name of the counter to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Minval	The minimum value, in units, that you expect to measure
Maxval	The maximum value, in units, that you expect to measure
Units	The units to use to return the measurements
Edge	Specifies between which edges to measure the frequency or period of the signal
Measmeth	The method used to calculate the period or frequency of the signal
Meastime	The length of time to measure the frequency or period of a digital signal, when measMethod is "DAQmx_Val_HighFreq2Ctr"
Divisor	The value by which to divide the input signal, when measMethod is "DAQmx_Val_LargeRng2Ctr"
Csn\$	The name of a custom scale to apply to the channel

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates a channel to measure the period of a digital signal and adds the channel to the task

**HTBasic Definition:**

```
LONG Result=FNHx_crciperichan((LONG Taskid), Counter$, My_chn$, (REAL Minval),
                              (REAL Maxval), (LONG Units), (LONG Edge),
                              (LONG Measmeth), (REAL Meastime), (LONG Divisor), Csn$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Counter\$	The name of the counter to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Minval	The minimum value, in units, that you expect to measure
Maxval	The maximum value, in units, that you expect to measure
Units	The units to use to return the measurements
Edge	Specifies between which edges to measure the frequency or period of the signal
Measmeth	The method used to calculate the period or frequency of the signal
Meastime	The length of time to measure the frequency or period of a digital signal, when measMethod is "DAQmx_Val_HighFreq2Ctr"
Divisor	The value by which to divide the input signal, when measMethod is "DAQmx_Val_LargeRng2Ctr"
Csn\$	The name of a custom scale to apply to the channel

Return value: 0 = Function was successful  
<>0 = Error occurred



**Purpose:** Creates a channel to count the number of rising or falling edges of a digital signal

**HTBasic Definition:**

```
LONG Result=FNHx_crcicntedgch((LONG Taskid), Counter$, My_chn$, (LONG Edge),
                               (LONG Initcount), (LONG Countdirection))
```

**Input:**

Taskid	The task to which to add the channels that this function creates
Lines\$	The names of the digital lines used to create a virtual channel
My_chn\$	The name(s) to assign to the created virtual channel(s)
Edge	Specifies on which edges of the input signal to increment or decrement the count
Initcount	The value from which to start counting
Countdirection	Specifies whether to increment or decrement the counter on each edge

**Return value:** 0 = Function was successful  
<>0 = Error occurred

**Example:**

```
INTEGER Iresult1, Iresult2
LONG Result, Taskid, Edge, Initcount, Countdir
DIM Counter$[255], Tmp$[255]
!
Counter$="Dev1/ctr0"           ! Counter 0
Tmp$=""
Edge=FNHx_get_constant("DAQmx_Val_Rising", Iresult1)
Initcount=0
Countdir=FNHx_get_constant("DAQmx_Val_CountUp", Iresult2)
Result=FNHx_crcicntedgch((Taskid), Counter$, Tmp$, (Edge), (Initcount), (Countdir))
```

**Purpose:** Creates a channel to measure the width of a digital pulse and adds the channel to the task

**HTBasic Definition:**

```
LONG Result=FNHx_crciplswchan((LONG Taskid), Counter$, My_chn$, (REAL Minval),
                               (REAL Maxval), (LONG Units), (LONG Startingedge), Csn$)
```

**Input:**

Taskid	The task to which to add the channels that this function creates
Counter\$	The name of the counter to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Minval	The minimum value, in units, that you expect to measure
Maxval	The maximum value, in units, that you expect to measure
Units	The units to use to return the measurements
Startingedge	Specifies on which edge to begin measuring pulse width
Csn\$	The name of a custom scale to apply to the channel

**Return value:** 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates a channel to measure the width of a digital pulse and adds the channel to the task

**HTBasic Definition:**

```
LONG Result=FNHx_crcisemperch((LONG Taskid), Counter$, My_chn$, (REAL Minval),
                               (REAL Maxval), (LONG Units), Csn$)
```

**Input:**

Taskid	The task to which to add the channels that this function creates
Counter\$	The name of the counter to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Minval	The minimum value, in units, that you expect to measure
Maxval	The maximum value, in units, that you expect to measure
Units	The units to use to return the measurements
Csn\$	The name of a custom scale to apply to the channel

**Return value:** 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates a channel that measures the amount of time between the rising or falling edge of one digital signal and the rising or falling edge of another digital signal

**HTBasic Definition:**

```
LONG Result=FNHx_crcitwedsech((LONG Taskid), Counter$, My_chn$, (REAL Minval),
                              (REAL Maxval), (LONG Units), (LONG Firstedge),
                              (LONG Secondedge), Csn$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Counter\$	The name of the counter to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Minval	The minimum value, in units, that you expect to measure
Maxval	The maximum value, in units, that you expect to measure
Units	The units to use to return the measurements
Firstedge	Specifies on which edge of the first signal to start each measurement
Secondedge	Specifies on which edge of the first signal to stop each measurement
Csn\$	The name of a custom scale to apply to the channel

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates a channel that uses a linear encoder to measure linear position

**HTBasic Definition:**

```
LONG Result=FNHx_crcilinencch((LONG Taskid), Counter$, My_chn$, (LONG Decodingtype),
                              (LONG Zidxenable), (REAL Zidxval), (LONG Zidxphase),
                              (LONG Units), (REAL Distppulse), (REAL Initpos), Csn$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Counter\$	The name of the counter to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Decodingtype	Specifies how to count and interpret the pulses that the encoder generates on signal A and signal B
Zidxenable	Specifies whether to enable z indexing for the measurement
Zidxval	The value, in units, to which to reset the measurement when signal Z is high and signal A and signal B are at the states you specify with <b>Zidxphase</b>
Zidxphase	The states at which signal A and signal B must be while signal Z is high for NI-DAQmx to reset the measurement
Units	The units to use to return the measurements
Distppulse	The distance measured for each pulse the encoder generates (in <b>Units</b> )
Initpos	The position of the encoder when the measurement begins (in <b>Units</b> )
Csn\$	The name of a custom scale to apply to the channel

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates a channel that uses an angular encoder to measure angular position

**HTBasic Definition:**

```
LONG Result=FNHx_crcianencchn((LONG Taskid), Counter$, My_chn$, (LONG Decodingtype),
                              (LONG Zidxenable), (REAL Zidxval), (LONG Zidxphase),
                              (LONG Units), (REAL Pulsesperref), (REAL Initang), Csn$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Counter\$	The name of the counter to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Decodingtype	Specifies how to count and interpret the pulses that the encoder generates on signal A and signal B
Zidxenable	Specifies whether to enable z indexing for the measurement
Zidxval	The value, in units, to which to reset the measurement when signal Z is high and signal A and signal B are at the states you specify with <b>Zidxphase</b>
Zidxphase	The states at which signal A and signal B must be while signal Z is high for NI-DAQmx to reset the measurement
Units	The units to use to return the measurements

Pulsesperref      The number of pulses the encoder generates per revolution  
 Initang            The starting angle of the encoder when measurement begins (in **Units**)  
 Csn\$                The name of a custom scale to apply to the channel

Return value:    0 = Function was successful  
                      <>0 = Error occurred

<b>FNHx_crcigpstschn</b>	<b>NI-DAQmx: DAQmxCreateCIGPSTimestampChan</b>
--------------------------	--

**Purpose:** Creates a channel that uses a special-purpose counter to take a timestamp and synchronizes that counter to a GPS receiver

**HTBasic Definition:**

**LONG** Result=FNHx\_crcigpstschn(**LONG** Taskid), Counter\$, My\_chn\$, (**LONG** Units),  
 (**LONG** Syncmethod), Csn\$)

Input:

Taskid              The task to which to add the channels that this function creates  
 Counter\$            The name of the counter to use to create virtual channels  
 My\_chn\$            The name(s) to assign to the created virtual channel(s)  
 Units                The units to use to return the measurements  
 Syncmethod         The method to use to synchronize the counter to a GPS receiver  
 Csn\$                The name of a custom scale to apply to the channel

Return value:    0 = Function was successful  
                      <>0 = Error occurred

<b>FNHx_crcopulchfrq</b>	<b>NI-DAQmx: DAQmxCreateCOPulseChanFreq</b>
--------------------------	---

**Purpose:** Creates channel(s) to generate digital pulses that **Freq** and **Dutycycle** define

**HTBasic Definition:**

**LONG** Result=FNHx\_crcopulchfrq(**LONG** Taskid), Counter\$, My\_chn\$, (**LONG** Units),  
 (**LONG** Idlestate), (**REAL** Initdelay), (**REAL** Freq),  
 (**REAL** Dutycycle))

Input:

Taskid              The task to which to add the channels that this function creates  
 Counter\$            The name of the counter to use to create virtual channels  
 My\_chn\$            The name(s) to assign to the created virtual channel(s)  
 Units                The units to use to return the measurements  
 Idlestate            The resting state of the output terminal  
 Initdelay            The amount of time in seconds to wait before generating the 1st pulse  
 Freq                 The frequency at which to generate pulses  
 Dutycycle            The width of the pulse divided by the pulse period.  
                        NI-DAQmx uses this ratio, combined with frequency, to determine pulse  
                        width and the interval between pulses

Return value:    0 = Function was successful  
                      <>0 = Error occurred

<b>FNHx_crcopulchtim</b>	<b>NI-DAQmx: DAQmxCreateCOPulseChanTime</b>
--------------------------	---

**Purpose:** Creates channel(s) to generate digital pulses defined by the amount of time the pulse is at a high state and the amount of time the pulse is at a low state

**HTBasic Definition:**

**LONG** Result=FNHx\_crcopulchtim(**LONG** Taskid), Counter\$, My\_chn\$, (**LONG** Units),  
 (**LONG** Idlestate), (**REAL** Initdelay), (**REAL** Lowtime),  
 (**REAL** Hightime))

Input:

Taskid              The task to which to add the channels that this function creates  
 Counter\$            The name of the counter to use to create virtual channels  
 My\_chn\$            The name(s) to assign to the created virtual channel(s)  
 Units                The units to use to return the measurements  
 Idlestate            The resting state of the output terminal  
 Initdelay            The amount of time in seconds to wait before generating the 1st pulse  
 Lowtime             The amount of time the pulse is low, in seconds  
 Hightime            The amount of time the pulse is high, in seconds

Return value:    0 = Function was successful  
                      <>0 = Error occurred

**Purpose:** Creates channel(s) to generate digital pulses defined by the number of timebase ticks that the pulse is at a high state and the number of timebase ticks that the pulse is at a low state

**HTBasic Definition:**

```
LONG Result=FNHx_crcopulchtck((LONG Taskid), Counter$, My_chn$, Srcterminal$,  
                               (LONG Idlestate), (REAL Initdelay), (REAL Lowticks),  
                               (REAL Highticks))
```

Input:

Taskid	The task to which to add the channels that this function creates
Counter\$	The name of the counter to use to create virtual channels
Srcterminal\$	The terminal to which you connect an external timebase
Idlestate	The resting state of the output terminal
Initdelay	The amount of time in seconds to wait before generating the 1st pulse
Lowticks	The number of timebase ticks that the pulse is low
Highticks	The number of timebase ticks that the pulse is high

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Creates channel(s) to measure RMS voltage and adds the channel(s) to the task you specify with **Taskid**

**HTBasic Definition:**

```
LONG Result=FNHx_craivlrmschn((LONG Taskid), Phys_chan$, My_chn$, (LONG Termcfg),
                               (REAL Minval), (REAL Maxval), (LONG Units), Custscalname$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Termcfg	The input terminal configuration for the channel
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the voltage measurements
Custscalname\$	The name of a custom scale to apply to the channel

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Creates channel(s) for RMS current measurement and adds the channel(s) to the task you specify with **Taskid**

**HTBasic Definition:**

```
LONG Result=FNHx_craicurmschn((LONG Taskid), Phys_chan$, My_chn$, (LONG Termcfg),
                               (REAL Minval), (REAL Maxval), (LONG Units),
                               (REAL Shuntresloc), (REAL Extshuntresval), Custscalname$)
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Termcfg	The input terminal configuration for the channel
Minval	The minimum value, in <b>Units</b> , that you expect to measure
Maxval	The maximum value, in <b>Units</b> , that you expect to measure
Units	The units to use to return the voltage measurements
Shuntresloc	The location of the shunt resistor
Extshuntresval	The value, in ohms, of an external shunt resistor
Custscalname\$	The name of a custom scale to apply to the channel

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Creates a channel for continually generating a waveform on the selected physical channel

**HTBasic Definition:**

```
LONG Result=FNHx_craofungnchn((LONG Taskid), Phys_chan$, My_chn$, (LONG Type),  
                               (REAL Freq), (REAL Amplitude), (LONG Offset))
```

Input:

Taskid	The task to which to add the channels that this function creates
Phys_chan\$	The names of the physical channels to use to create virtual channels
My_chn\$	The name(s) to assign to the created virtual channel(s)
Type	Specifies the kind of waveform to generate
Freq	The frequency of the waveform to generate in Hertz
Amplitude	The zero-to-peak amplitude of the waveform to generate in volts
Offset	The voltage offset of the waveform to generate

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

## Timing Functions

**FNHx\_cfgsamclktim**

NI-DAQmx: **DAQmxCfgSampClkTiming**

**Purpose:** Sets the source of the Sample Clock, the rate of the Sample Clock, and the number of samples to acquire or generate

### HTBasic Definition:

```
LONG Result=FNHx_cfgsamclktim((LONG Taskid), Source$, (REAL Rate), (LONG Activeedge),  
                               (LONG Samplemode), (LONG Sampsperchan))
```

#### Input:

Taskid	The task to which to add the channels that this function creates
Source\$	The source terminal of the Sample Clock
Rate	The sampling rate in samples per second
Activeedge	Specifies on which edge of the clock to acquire or generate samples
Samplemode	Specifies whether the task acquires or generates samples continuously or if it acquires or generates a finite number of samples
Sampsperchan	The number of samples to acquire or generate for each channel in the task if <b>Samplemode</b> is "DAQmx_Val_FiniteSamps". If <b>Samplemode</b> is "DAQmx_Val_ContSamps", NI-DAQmx uses this value to determine the buffer size.

**Return value:** 0 = Function was successful  
<>0 = Error occurred

### Example:

```
INTEGER Iresult1, Iresult2  
LONG Result, Taskid, Sampsperchn, Actedge, Smode  
REAL Rate  
DIM Tmp$[255], Devchn$[255]  
!  
Sampsperchn=100 ! Samples per channel (max. 32767)  
Actedge=FNHx_get_constant("DAQmx_Val_Rising", Iresult1) ! get constant 1  
Smode=FNHx_get_constant("DAQmx_Val_ContSamps", Iresult2) ! get constant 2  
Rate=20000.0 ! 20.000 Samples / s  
Result=FNHx_cfgsamclktim((Taskid), Tmp$, (Rate), (Actedge), (Smode), (Sampsperchn))
```

**FNHx\_cfgplscclktim**

NI-DAQmx 8.5: **DAQmxCfgPipelinedSampClkTiming**

**Purpose:** Sets the source of the Sample Clock, the rate of the Sample Clock, and the number of samples to acquire or generate

### HTBasic Definition:

```
LONG Result=FNHx_cfgplscclktim((LONG Taskid), Source$, (REAL Rate), (LONG Activeedge),  
                               (LONG Samplemode), (LONG Sampsperchan))
```

#### Input:

Taskid	The task to which to add the channels that this function creates
Source\$	The source terminal of the Sample Clock
Rate	The sampling rate in samples per second
Activeedge	Specifies on which edge of the clock to acquire or generate samples
Samplemode	Specifies whether the task acquires or generates samples continuously or if it acquires or generates a finite number of samples
Sampsperchan	The number of samples to acquire or generate for each channel in the task if <b>Samplemode</b> is "DAQmx_Val_FiniteSamps". If <b>Samplemode</b> is "DAQmx_Val_ContSamps", NI-DAQmx uses this value to determine the buffer size.

**Return value:** 0 = Function was successful  
<>0 = Error occurred

### Note:

This timing type allows handshaking using the Pause trigger, the Ready for Transfer event, or the Data Active event. Refer to the device documentation for more information.

This timing type is supported only by the NI 6536 and NI 6537 and is available with NI-DAQmx 8.5 and higher only.

**Purpose:** Determines the number of digital samples to acquire or generate using digital handshaking between the device and a peripheral device

**HTBasic Definition:**

```
LONG Result=FNHx_cfghshtiming((LONG Taskid), (LONG Samplemode), (LONG Sampsperchn_msb),  
                               (LONG Sampsperchn_lsb))
```

Input:

Taskid	The task to which to add the channels that this function creates
Samplemode	Specifies whether the task acquires or generates samples continuously or if it acquires or generates a finite number of samples
Sampsperchn_msb	Upper 32 Bit of Samples per channel variable (64 Bit)
Sampsperchn_lsb	Lower 32 Bit of Samples per channel variable (64 Bit)

The number of samples to acquire or generate for each channel in the task if **Samplemode** is "DAQmx\_Val\_FiniteSamps". If **Samplemode** is "DAQmx\_Val\_ContSamps", NI-DAQmx uses this value to determine the buffer size

Return value: 0 = Function was successful  
<>0 = Error occurred



**Purpose:** Configures when the DAQ device transfers data to a peripheral device, using an imported sample clock to control burst handshaking timing

**HTBasic Definition:**

```
LONG Result=FNHx_cfgbhtimiclk((LONG Taskid), (LONG Samplemode), (LONG Sampspcrhn_msb),
                               (LONG Sampspcrhn_lsb), (REAL Sampclockrate), Sampclocksrc$,
                               (LONG Sampclktactedge), (LONG Pausewhen), (LONG Rdyevntactlev))
```

Input:

Taskid	The task to which to add the channels that this function creates
Samplemode	Specifies whether the task acquires or generates samples continuously or if it acquires or generates a finite number of samples
Sampspcrhn_msb	Upper 32 Bit of Samples per channel variable (64 Bit)
Sampspcrhn_lsb	Lower 32 Bit of Samples per channel variable (64 Bit)
	The number of samples to acquire or generate for each channel in the task if <b>Samplemode</b> is "DAQmx_Val_FiniteSamps". If <b>Samplemode</b> is "DAQmx_Val_ContSamps", NI-DAQmx uses this value to determine the buffer size
Sampclockrate	Specifies the sampling rate in samples per channel per second
Sampclocksrc\$	Specifies the terminal of the signal to use as the Sample Clock
Sampclktactedge	Specifies on which edge of a clock pulse sampling takes place
Pausewhen	Specifies whether the task pauses while the signal is high or low
Rdyevntactlev	Specifies the polarity for the Ready for Transfer event

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Configures when the DAQ device transfers data to a peripheral device, using the DAQ device's onboard sample clock to control burst handshaking timing

**HTBasic Definition:**

```
LONG Result=FNHx_cfgbhtimeclk((LONG Taskid), (LONG Samplemode), (LONG Sampspcrhn_msb),
                               (LONG Sampspcrhn_lsb), (REAL Sampclockrate), Sampclkoutterm$,
                               (LONG Sampclkpulspol), (LONG Pausewhen), (LONG Rdyevntactlev))
```

Input:

Taskid	The task to which to add the channels that this function creates
Samplemode	Specifies whether the task acquires or generates samples continuously or if it acquires or generates a finite number of samples
Sampspcrhn_msb	Upper 32 Bit of Samples per channel variable (64 Bit)
Sampspcrhn_lsb	Lower 32 Bit of Samples per channel variable (64 Bit)
	The number of samples to acquire or generate for each channel in the task if <b>Samplemode</b> is "DAQmx_Val_FiniteSamps". If <b>Samplemode</b> is "DAQmx_Val_ContSamps", NI-DAQmx uses this value to determine the buffer size
Sampclockrate	Specifies the sampling rate in samples per channel per second
Sampclkoutterm\$	Specifies the terminal to which to route the Sample Clock
Sampclkpulspol	Specifies if the polarity for the exported sample clock is active high or active low
Pausewhen	Specifies whether the task pauses while the signal is high or low
Rdyevntactlev	Specifies the polarity for the Ready for Transfer event

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Configures the task to acquire samples on the rising and/or falling edges of the lines or ports you specify

**HTBasic Definition:**

```
LONG Result=FNHx_cfgchgdettim((LONG Taskid), Risedgechan$, Faledgechan$, (LONG Smode),
                              (LONG Sampsperchn_msb), (LONG Sampsperchn_lsb))
```

Input:

Taskid	The task to which to add the channels that this function creates
Risedgechan\$	The names of the digital lines or ports on which to detect rising edges
Faledgechan\$	The names of the digital lines or ports on which to detect falling edges
Smode	Specifies whether the task acquires or generates samples continuously or if it acquires or generates a finite number of samples
Sampsperchn_msb	Upper 32 Bit of Samples per channel variable (64 Bit)
Sampsperchn_lsb	Lower 32 Bit of Samples per channel variable (64 Bit)

The number of samples to acquire or generate for each channel in the task if **Samplemode** is "DAQmx\_Val\_FiniteSamps". If **Samplemode** is "DAQmx\_Val\_ContSamps", NI-DAQmx uses this value to determine the buffer size

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Sets only the number of samples to acquire or generate without specifying timing

**HTBasic Definition:**

```
LONG Result=FNHx_cfgimptiming((LONG Taskid), (LONG Samplemode), (LONG Sampchtoacquire))
```

Input:

Taskid	The task to which to add the channels that this function creates
Samplemode	Specifies whether the task acquires or generates samples continuously or if it acquires or generates a finite number of samples
Sampchtoacquire	The number of samples to acquire or generate for each channel in the task if <b>Samplemode</b> is "DAQmx_Val_FiniteSamps". If <b>Samplemode</b> is "DAQmx_Val_ContSamps", NI-DAQmx uses this value to determine the buffer size

Return value: 0 = Function was successful  
<>0 = Error occurred

**Example:**

```
INTEGER Iresult1
LONG Result, Taskid, Samplemode, Sampchtoacquire
!
Samplemode=FNHx_get_constant("DAQmx_Val_FiniteSamps", Iresult1)      ! Sample mode
Sampchtoacquire=10
Result=FNHx_cfgimptiming((Taskid), (Samplemode), (Sampchtoacquire))
```

## Read Functions

**FNHx\_readanalf64**

NI-DAQmx: **DAQmxReadAnalogF64**

**Purpose:** Reads multiple floating-point samples from a task that contains one or more analog input channels

### HTBasic Definition:

```
LONG Result=FNHx_readanalf64((LONG Taskid), (LONG Numpersam), (REAL T_out),  
                              (LONG Fillmode), REAL Data_in(*), (LONG Asize),  
                              LONG Sperchnread, LONG Reserved)
```

### Input:

Taskid	The task to read samples from
Numpersam	The number of samples, per channel, to read
T_out	The amount of time, in seconds, to wait for the function to read the sample(s)
Fillmode	Specifies whether or not the samples are interleaved
Asize	The size of the array, in samples, into which samples are read
Reserved	Reserved for future use. Pass 0 to this parameter

**Return value:** 0 = Function was successful  
                  **Sperchnread** = Number of read REAL values  
                  **Data\_in(\*)** = Measurements  
                  <>0 = Error occurred

### Example:

```
INTEGER Iresult1  
LONG Result, Taskid, Numpersam, Sperchnread, Reserved  
REAL T_out  
REAL Data_in(32767)  
DIM Tmp$(255)  
!  
Fillmode=FNHx_get_constant("DAQmx_Val_GroupByScanNumber",Iresult1)  
T_out=5.0  
Asize=SIZE(Data_in,1)           ! Size of array  
Numpersam=100                  ! read 100 samples  
Reserved=0  
Result=FNHx_readanalf64((Taskid), (Numpersam), (T_out), (Fillmode), Data_in(*), (Asize),  
                        Sperchnread, Reserved)  
!  
IF Result=0 THEN  
    !... subsequent commands...  
END IF
```

**FNHx\_readanalsf64**

NI-DAQmx: **DAQmxReadAnalogScalarF64**

**Purpose:** Reads a single floating-point sample from a task that contains a single analog input channel

### HTBasic Definition:

```
LONG Result=FNHx_readanalsf64((LONG Taskid), (REAL T_out), REAL Data_in, LONG Reserved)
```

### Input:

Taskid	The task to read samples from
T_out	The amount of time, in seconds, to wait for the function to read the sample(s)
Reserved	Reserved for future use. Pass 0 to this parameter

**Return value:** 0 = Function was successful  
                  REAL variable **Data\_in** contains the measurement  
                  <>0 = Error occurred

### Example:

```
LONG Result, Taskid, Reserved  
REAL T_out, Avalue  
!  
T_out=5.0  
Reserved=0  
Result=FNHx_readanalsf64((Taskid), (T_out), Avalue, Reserved)
```

**Purpose:** Reads multiple unscaled, signed 16-bit integer samples

**HTBasic Definition:**

```
LONG Result=FNHx_rdbinaryi16((LONG Taskid), (LONG Numpersam), (REAL T_out),
                             (LONG Fillmode), INTEGER Data_in(*), (LONG Asize),
                             LONG Sperchnread, LONG Reserved)
```

Input:

Taskid	The task to read samples from
Numpersam	The number of samples, per channel, to read
T_out	The amount of time, in seconds, to wait for the function to read the sample(s)
Fillmode	Specifies whether or not the samples are interleaved
Asize	The size of the array, in samples, into which samples are read
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
                   **Sperchnread** = Number of read INTEGER values from each channel  
                   **Data\_in(\*)** = 16 Bit unscaled measurements  
 <>0 = Error occurred

**Purpose:** Reads multiple unscaled, signed 32-bit integer samples

**HTBasic Definition:**

```
LONG Result=FNHx_rdbinaryi32((LONG Taskid), (LONG Numpersam), (REAL T_out),
                              (LONG Fillmode), LONG Data_in(*), (LONG Asize),
                              LONG Sperchnread, LONG Reserved)
```

Input:

Taskid	The task to read samples from
Numpersam	The number of samples, per channel, to read
T_out	The amount of time, in seconds, to wait for the function to read the sample(s)
Fillmode	Specifies whether or not the samples are interleaved
Asize	The size of the array, in samples, into which samples are read
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
                   **Sperchnread** = Number of read INTEGER values from each channel  
                   **Data\_in(\*)** = 32 Bit unscaled measurements  
 <>0 = Error occurred

**Purpose:** Reads multiple 16-bit integer samples

**HTBasic Definition:**

```
LONG Result=FNHx_readdigitu16((LONG Taskid), (LONG Sampsperchan), (REAL T_out),
                               (LONG Fillmode), INTEGER Dig_in(*), (LONG Asize),
                               LONG Samsreaddig, LONG Reserved)
```

Input:

Taskid	The task to read samples from
Numpersam	The number of samples, per channel, to read
T_out	The amount of time, in seconds, to wait for the function to read the sample(s)
Fillmode	Specifies whether or not the samples are interleaved
Asize	The size of the array, in samples, into which samples are read
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
                   **Dig\_in** = 16 Bit digital measurements  
                   **Samsreaddig** = Number of read INTEGER values from each channel  
 <>0 = Error occurred

**Purpose:** Reads multiple 32-bit integer samples

**HTBasic Definition:**

```
LONG Result=FNHx_readdigitu32((LONG Taskid), (LONG Sampsperchan), (REAL T_out),
                              (LONG Fillmode), INTEGER Dig_in(*), (LONG Asize),
                              LONG Samsreaddig, LONG Reserved)
```

Input:

Taskid	The task to read samples from
Numpersam	The number of samples, per channel, to read
T_out	The amount of time, in seconds, to wait for the function to read the sample(s)
Fillmode	Specifies whether or not the samples are interleaved
Asize	The size of the array, in samples, into which samples are read
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
                   **Dig\_in** = 32 Bit digital measurements  
                   **Samsreaddig** = Number of read INTEGER values from each channel  
 <>0 = Error occurred

**Example:**

```
INTEGER Iresult1, Dig_in(2000)
LONG Result, Taskid, Sampsperchan, Asize, Samsreaddig, Reserved
REAL T_out
!
Sampsperchan=4           ! number of samples per channel
T_out=10.0               ! 10 seconds timeout
Asize=SIZE(Dig_in,1)    ! size of array
MAT Dig_in=(0)
Reserved=0
!
Fillmode=FNHx_get_constant("DAQmx_Val_GroupByChannel",Iresult1) ! constant 1
Result=DEF FNHx_readdigitu32((Taskid), (Sampsperchan), (T_out), (Fillmode), Dig_in(*),
                              (Asize), Samsreaddig, Reserved)
```

**Purpose:** Reads a single 32-bit integer sample

**HTBasic Definition:**

```
LONG Result=FNHx_readdigits32((LONG Taskid), (REAL T_out), LONG Value, LONG Reserved)
```

Input:

Taskid	The task to read samples from
T_out	The amount of time, in seconds, to wait for the function to read the sample(s)
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
                   **Value** = Measurement  
 <>0 = Error occurred

**Example:**

```
LONG Result, Taskid, Value, Reserved
REAL T_out
!
T_out=10.0               ! 10s Timeout
Reserved=0
!
Result=DEF FNHx_readdigitu32((Taskid), (T_out), Value, Reserved)
```

**Purpose:** Reads multiple samples from each digital line in a task

**HTBasic Definition:**

```
LONG Result=FNHx_readdiglines((LONG Taskid), (LONG Numpersam), (REAL T_out),
                               (LONG Fillmode), INTEGER Dig_in(*), (LONG Asize),
                               LONG Samsreaddig, LONG Numbytespersamp, LONG Reserved)
```

Input:

Taskid	The task to read samples from
Numpersam	The number of samples, per channel, to read
T_out	The amount of time, in seconds, to wait for the function to read the sample(s)
Fillmode	Specifies whether or not the samples are interleaved
Asize	The size of the array, in samples, into which samples are read
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful

<b>Dig_in</b>	= Measurements
<b>Samsreaddig</b>	= Number of read samples from each channel
<b>Numbytespersamp</b>	= The number of elements in readArray that constitutes a sample per channel. For each sample per channel, numBytesPerSamp is the number of bytes that channel consists of.

<>0 = Error occurred

**Example:**

```
INTEGER Iresult1, Dig_in(2000)
LONG Result, Taskid, Sampsperchan, Asize, Samsreaddig, Numbytespersamp, Reserved
REAL T_out
!
T_out=10.0
Sampsperchan=4           ! number of samples (Bits) per channel
Asize=SIZE(Dig_in,1)     ! array size
MAT Dig_in=(0)
!
Fillmode=FNHx_get_constant("DAQmx_Val_GroupByChannel",Iresult1) ! constant 1
Result=FNHx_readdiglines((Taskid),(Sampsperchan),(T_out),(Fillmode),Dig_in(*),(Asize),
                          Samsreaddig,Numbytespersamp,Reserved)
```

**Purpose:** Reads multiple floating-point samples from a counter task

**HTBasic Definition:**

```
LONG Result=FNHx_readcntrf64((LONG Taskid), (LONG Numpersam), (REAL T_out),
                              REAL Data(*), (LONG Asize), LONG Samsread, LONG Reserved)
```

Input:

Taskid	The task to read samples from
Numpersam	The number of samples, per channel, to read
T_out	The amount of time, in seconds, to wait for the function to read the sample(s)
Asize	The size of the array, in samples, into which samples are read
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful

<b>Data</b>	= Array filled with samples
<b>Samsread</b>	= Number of read INTEGER values from each channel

<>0 = Error occurred

**Purpose:** Reads multiple 32-bit integer samples from a counter task

**HTBasic Definition:**

```
LONG Result=FNHx_readcntru32((LONG Taskid), (LONG Numpersam), (REAL T_out),
                             LONG Data(*), (LONG Asize), LONG Samsread, LONG Reserved)
```

Input:

Taskid	The task to read samples from
Numpersam	The number of samples, per channel, to read
T_out	The amount of time, in seconds, to wait for the function to read the sample(s)
Asize	The size of the array, in samples, into which samples are read
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
                   **Data** = Array filled with samples  
                   **Samsread** = Number of read INTEGER values from each channel  
 <>0 = Error occurred

**Purpose:** Reads a single floating-point sample from a counter task

**HTBasic Definition:**

```
LONG Result=FNHx_rdcntscal64((LONG Taskid), (REAL Timeout), REAL Value, LONG Reserved)
```

Input:

Taskid	The task to read samples from
T_out	The amount of time, in seconds, to wait for the function to read the sample(s)
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
                   **Value** = Measurement  
 <>0 = Error occurred

**Example:**

```
LONG Result, Taskid, Reserved
REAL T_out, Value
T_out=10.0           ! 10s Timeout
Reserved=0
Result=FNHx_rdcntscal64((Taskid), T_out, Value, Reserved)
```

**Purpose:** Reads a 32-bit integer sample from a counter task

**HTBasic Definition:**

```
LONG Result=FNHx_rdcntscal32((LONG Taskid), (REAL Timeout), LONG Value, LONG Reserved)
```

Input:

Taskid	The task to read samples from
T_out	The amount of time, in seconds, to wait for the function to read the sample(s)
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
                   **Value** = Measurement  
 <>0 = Error occurred

**Example:**

```
LONG Result, Taskid, Value, Reserved
REAL T_out
T_out=10.0           ! 10s Timeout
Reserved=0
Result=FNHx_rdcntscal32((Taskid), (T_out), Value, Reserved)
```

## Write Functions

**FNHx\_writeanalf64**

NI-DAQmx: **DAQmxWriteAnalogF64**

**Purpose:** Writes multiple floating-point samples to a task

### HTBasic Definition:

```
LONG Result=FNHx_writeanalf64((LONG Taskid), (LONG Numpersam), (LONG Autostart),  
                              (REAL T_out), (LONG Dlayout), REAL Data_out(*),  
                              LONG Smpchwrit, LONG Reserved)
```

### Input:

Taskid	The task to write the samples to
Numpersam	The number of samples, per channel, to read
Autostart	Specifies whether or not this function automatically starts the task if you do not start it
T_out	The amount of time, in seconds, to wait for the function to read the sample(s)
Dlayout	Specifies how the samples are arranged (interleaved / noninterleaved)
Data_out(*)	The array of 64-bit samples to write to the task
Reserved	Reserved for future use. Pass 0 to this parameter

**Return value:** 0 = Function was successful  
Data from **Data\_out(\*)** written to the buffer  
**Smpchwrit** The actual number of samples per channel successfully written to the buffer  
<>0 = Error occurred

### Example:

```
INTEGER Iresult1  
LONG Result, Taskid, Numpersam, Autostart, Sampsperchan, Dlayout, Reserved, Smpchwrit  
REAL T_out, Data_out(32767)  
DIM Tmp$(255)  
!  
Sampsperchan=1  
Autostart=1  
T_out=10.0  
Data_out(1)=PI  
Dlayout=FNHx_get_constant("DAQmx_Val_GroupByChannel",Iresult1) ! get constant  
Result=FNHx_writeanalf64((Taskid), (Sampsperchan), (Autostart), (T_out), (Dlayout),  
                          Data_out(1), Smpchwrit, Reserved)
```

**FNHx\_writeanasf64**

NI-DAQmx: **DAQmxWriteAnalogScalarF64**

**Purpose:** Writes a floating-point sample to a task

### HTBasic Definition:

```
LONG Result=FNHx_writeanasf64((LONG Taskid), (LONG Autostart), (REAL T_out),  
                              (REAL Value), LONG Reserved)
```

### Input:

Taskid	The task to write the samples to
Autostart	Specifies whether or not this function automatically starts the task if you do not start it
T_out	The amount of time, in seconds, to wait for the function to read the sample(s)
Value	The 64-bit REAL variable to write to the task
Reserved	Reserved for future use. Pass 0 to this parameter

**Return value:** 0 = Function was successful  
Data from **Data\_out(\*)** written to the buffer  
<>0 = Error occurred

### Example:

```
LONG Result, Taskid, Autostart, Reserved  
REAL T_out, Value  
!  
Autostart=1  
T_out=10.0  
Value=1.654  
Reserved=0  
Result=FNHx_writeanasf64((Taskid), (Autostart), (T_out), (Value), Reserved)
```



**Purpose:** Writes unscaled signed integer samples to the task

**HTBasic Definition:**

```
LONG Result=FNHx_writebini16((LONG Taskid), (LONG Numpersam), (LONG Autostart),
                              (REAL T_out), (LONG Dlayout), INTEGER Data(*),
                              LONG Sperchnwrit, LONG Reserved)
```

Input:

Taskid	The task to write the samples to
Numpersam	The number of samples, per channel, to read
Autostart	Specifies whether or not this function automatically starts the task
T_out	The amount of time, in seconds, to wait for the function to write the sample(s)
Dlayout	Specifies whether or not the samples are interleaved
Data(*)	Array of 16-Bit samples to write to the task
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
**Sperchnwrit** = Number of INTEGER values written to buffer  
 <>0 = Error occurred

**Purpose:** Writes unscaled 32-bit integer samples to the task

**HTBasic Definition:**

```
LONG Result=FNHx_writebini32((LONG Taskid), (LONG Numpersam), (LONG Autostart),
                              (REAL T_out), (LONG Dlayout), LONG Data(*),
                              LONG Sperchnwrit, LONG Reserved)
```

Input:

Taskid	The task to write the samples to
Numpersam	The number of samples, per channel, to write
Autostart	Specifies whether or not this function automatically starts the task
T_out	The amount of time, in seconds, to wait for the function to write the sample(s)
Dlayout	Specifies whether or not the samples are interleaved
Data(*)	Array of 32-Bit samples to write to the task
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
**Sperchnwrit** = Number of LONG values written to buffer  
 <>0 = Error occurred

**Purpose:** Writes multiple 16-bit unsigned integer samples to a task that contains one or more digital output channels. Use this format for devices with up to 16 lines per port

**HTBasic Definition:**

```
LONG Result=FNHx_writdigitu16((LONG Taskid), (LONG Sampsperschan), (LONG Autostart),
                              (REAL T_out), (LONG Datalayout), INTEGER Dig_out(*),
                              LONG Sampswritten, LONG Reserved)
```

Input:

Taskid	The task to write the samples to
Sampsperschan	The number of samples, per channel, to write
Autostart	Specifies whether or not this function automatically starts the task
T_out	The amount of time, in seconds, to wait for the function to write the sample(s)
Datalayout	Specifies whether or not the samples are interleaved
Dig_out(*)	INTEGER array with data to be written to the buffer
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
**Sampswritten** = Number of values written to buffer for each channel  
 <>0 = Error occurred

**Purpose:** Writes multiple 32-bit unsigned integer samples to a task that contains one or more digital output channels. Use this format for devices with up to 32 lines per port

**HTBasic Definition:**

```
LONG Result=FNHx_writdigitu32((LONG Taskid), (LONG Sampsperchan), (LONG Autostart),
                              (REAL T_out), (LONG Datalayout), (LONG Asize),
                              LONG Dig_out(*), LONG Sampswritten, LONG Reserved)
```

Input:

Taskid	The task to write the samples to
Sampsperchan	The number of samples, per channel, to write
Autostart	Specifies whether or not this function automatically starts the task
T_out	The amount of time, in seconds, to wait for the function to write the sample(s)
Datalayout	Specifies whether or not the samples are interleaved
Dig_out(*)	LONG array with data to be written to the buffer
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
                   **Sampswritten** = Number of values written to buffer for each channel  
                   <>0 = Error occurred

**Purpose:** Writes a single 32-bit unsigned integer sample to a task

**HTBasic Definition:**

```
LONG Result=FNHx_writdigscu32((LONG Taskid), (LONG Autostart), (REAL T_out), LONG Value,
                              LONG Reserved)
```

Input:

Taskid	The task to write the samples to
Autostart	Specifies whether or not this function automatically starts the task
T_out	The amount of time, in seconds, to wait for the function to write
Value	32-Bit sample to write to the task
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Example:**

```
LONG Result, Taskid, Value, Reserved
REAL T_out
!
Reserved=0
Autostart=1
T_out=5
Value=255
Result=FNHx_writdigscu32((Digitalout_task), (Autostart), (T_out), Value, Reserved)
```

**Purpose:** Writes multiple samples to each digital line in a task

**HTBasic Definition:**

```
LONG Result=FNHx_writdiglines((LONG Taskid), (LONG Sampsperchan), (LONG Autostart),
                              (REAL T_out), (LONG Datalayout), (LONG Asize),
                              INTEGER Dig_out(*), LONG Sampswritten, LONG Reserved)
```

Input:

Taskid	The task to write the samples to
Sampsperchan	The number of samples, per channel, to write
Autostart	Specifies whether or not this function automatically starts the task
T_out	The amount of time, in seconds, to wait for the function to write the sample(s)
Datalayout	Specifies whether or not the samples are interleaved
Asize	Size of data array
Dig_out(*)	INTEGER array with data to be written to the buffer
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
                   **Sampswritten** = Number of values written to buffer for each channel  
                   <>0 = Error occurred

**Example:**

```

INTEGER Iresult1
LONG Result, Taskid, Sampsperchan, Autostart, Asize, Sampswritten, Datalayout
LONG Dig_out(32767), Reserved
REAL T_out
!
Sampsperchan=2           ! number of sample (Bits) per channel
Autostart=1
T_out=10.0
Asize=SIZE(Dig_out,1)   ! array size
Dig_out(1)=255          ! 1. output value
Dig_out(2)=0            ! 2. output value
Reserved=0
!
Datalayout=FNHx_get_constant("DAQmx_Val_GroupByChannel",Iresult1) ! constant 1
Result=FNHx_writdiglines((Taskid),(Sampsperchan),(Autostart),(T_out),
                        (Datalayout),(Asize),Dig_out(*),Sampswritten,Reserved)

```

**FNHx\_wrctrfreq**NI-DAQmx: **DAQmxWriteCtrFreq****Purpose:** Writes a new pulse frequency and duty cycle to each channel**HTBasic Definition:**

```

LONG Result=FNHx_wrctrfreq((LONG Taskid), (LONG Sampsperchan), (LONG Autostart),
                          (REAL T_out), (LONG Dlayout), (REAL Freq), (REAL Dutycycle),
                          LONG Sampswritten, LONG Reserved)

```

Input:

Taskid	The task to write the values to
Sampsperchan	The number of samples, per channel, to write
Autostart	Specifies whether or not this function automatically starts the task
T_out	The amount of time, in seconds, to wait for the function to write the sample(s)
Dlayout	Specifies whether or not the samples are interleaved
Freq	The frequency at which to generate pulses
Dutycycle	The width of the pulse divided by the pulse period. NI-DAQmx uses this ratio, combined with frequency, to determine pulse width and the interval between pulses
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
**Sampswritten** = Number of values written to buffer for each channel  
 <>0 = Error occurred

**Purpose:** Writes a new pulse frequency and duty cycle to each channel

**HTBasic Definition:**

```
LONG Result=FNHx_wrctrfrqscal((LONG Taskid), (LONG Autostart), (REAL T_out),
                              (REAL Freq), (REAL DutyCycle), LONG Reserved)
```

Input:

Taskid	The task to write the values to
Autostart	Specifies whether or not this function automatically starts the task
T_out	The amount of time, in seconds, to wait for the function to write the sample(s)
Freq	The frequency at which to generate pulses
DutyCycle	The width of the pulse divided by the pulse period. NI-DAQmx uses this ratio, combined with frequency, to determine pulse width and the interval between pulses
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Writes a new pulse high time and low time to each channel

**HTBasic Definition:**

```
LONG Result=FNHx_wrctrtime((LONG Taskid), (LONG Sampsperschan), (LONG Autostart),
                           (REAL T_out), (LONG Dlayout), (REAL Hightime),
                           (REAL Lowtime), LONG Sampswritten, LONG Reserved)
```

Input:

Taskid	The task to write the values to
Sampsperschan	The number of samples, per channel, to write
Autostart	Specifies whether or not this function automatically starts the task
T_out	The amount of time, in seconds, to wait for the function to write the sample(s)
Dlayout	Specifies whether or not the samples are interleaved
Hightime	The amount of time the pulse is high, in seconds
Lowtime	The amount of time the pulse is low, in seconds
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
                  **Sampswritten** = Number of values written to buffer for each channel  
<>0 = Error occurred

**Purpose:** Writes a new pulse high time and low time to each channel

**HTBasic Definition:**

```
LONG Result=FNHx_wrctrtime scal((LONG Taskid), (LONG Autostart), (REAL T_out),
                                 (REAL Hightime), (REAL Lowtime), LONG Reserved)
```

Input:

Taskid	The task to write the values to
Autostart	Specifies whether or not this function automatically starts the task
T_out	The amount of time, in seconds, to wait for the function to write the sample(s)
Hightime	The amount of time the pulse is high, in seconds
Lowtime	The amount of time the pulse is low, in seconds
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Writes new pulse high tick counts and low tick counts to each channel

**HTBasic Definition:**

```
LONG Result=FNHx_wrcrticks((LONG Taskid), (LONG Sampsperschan), (LONG Autostart),
                           (REAL T_out), (LONG Dlayout), (REAL Highticks),
                           (REAL Lowticks), LONG Sampswritten, LONG Reserved)
```

Input:

Taskid	The task to write the values to
Sampsperschan	The number of samples, per channel, to write
Autostart	Specifies whether or not this function automatically starts the task
T_out	The amount of time, in seconds, to wait for the function to write the sample(s)
Dlayout	Specifies whether or not the samples are interleaved
Highticks	The number of timebase ticks the pulse is high
Lowticks	The number of timebase ticks the pulse is low
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
                   **Sampswritten** = Number of values written to buffer for each channel  
                   <>0 = Error occurred

**Purpose:** Writes a new pulse high time and low time to each channel

**HTBasic Definition:**

```
LONG Result=FNHx_wrcrtickssc((LONG Taskid), (LONG Autostart), (REAL T_out),
                             (REAL Highticks), (REAL Lowticks), LONG Reserved)
```

Input:

Taskid	The task to write the values to
Autostart	Specifies whether or not this function automatically starts the task
T_out	The amount of time, in seconds, to wait for the function to write the sample(s)
Highticks	The number of timebase ticks the pulse is high
Lowticks	The number of timebase ticks the pulse is low
Reserved	Reserved for future use. Pass 0 to this parameter

Return value: 0 = Function was successful  
                   <>0 = Error occurred

## Triggering Functions

**FNHx\_disstarttrig**

NI-DAQmx: **DAQmxDisableStartTrig**

**Purpose:** Configures the task to start acquiring or generating samples immediately upon starting the task

**HTBasic Definition:**

`LONG Result=FNHx_disstarttrig((LONG Taskid))`

Input:

Taskid                    The task used in this function

Return value:    0 = Function was successful  
                  <>0 = Error occurred

**FNHx\_cfdgedsttrig**

NI-DAQmx: **DAQmxCfgDigEdgeStartTrig**

**Purpose:** Configures the task to start acquiring or generating samples on a rising or falling edge of a digital signal

**HTBasic Definition:**

`LONG Result=FNHx_cfdgedsttrig((LONG Taskid), (Triggersrc$), (LONG Triggeredge))`

Input:

Taskid                    The task used in this function  
Triggersrc\$              The name of a terminal where there is a digital signal to use as the source of the trigger  
Triggeredge              Specifies on which edge of a digital signal to start acquiring or generating samples

Return value:    0 = Function was successful  
                  <>0 = Error occurred

**FNHx\_cfganedstrig**

NI-DAQmx: **DAQmxCfgAnlgEdgeStartTrig**

**Purpose:** Configures the task to start acquiring or generating samples when an analog signal crosses the level you specify

**HTBasic Definition:**

`LONG Result=FNHx_cfganedstrig((LONG Taskid), (Triggersrc$), (LONG Triggerslop), (REAL Triggerlev))`

Input:

Taskid                    The task used in this function  
Triggersrc\$              The name of a channel or terminal where there is an analog signal to use as the source of the trigger  
Triggerslop              Specifies on which slope of the signal to start acquiring or generating samples when the signal crosses **Triggerlev**  
Triggerlev                The threshold at which to start acquiring or generating samples

Return value:    0 = Function was successful  
                  <>0 = Error occurred

**FNHx\_cfganwistrig**

NI-DAQmx: **DAQmxCfgAnlgWindowStartTrig**

**Purpose:** Configures the task to start acquiring or generating samples when an analog signal enters or leaves a range you specify

**HTBasic Definition:**

`LONG Result=FNHx_cfganwistrig((LONG Taskid), (Triggersrc$), (LONG Triggerwhen), (REAL Windowtop), (REAL Windowbottom))`

Input:

Taskid                    The task used in this function  
Triggersrc\$              The name of a channel or terminal where there is an analog signal to use as the source of the trigger  
Triggerwhen              Specifies whether the task starts measuring or generating samples when the signal enters the window or when it leaves the window  
Windowtop                The upper limit of the window  
Windowbottom            The lower limit of the window

Return value:    0 = Function was successful  
                  <>0 = Error occurred

**Purpose:** Configures a task to start acquiring or generating samples when a digital pattern is matched

**HTBasic Definition:**

```
LONG Result=FNHx_cfgdiptrig((LONG Taskid), (Triggersrc$), (Triggerpat$),
                             (LONG Triggerwhen))
```

Input:

Taskid	The task used in this function
Triggersrc\$	Specifies the physical channels to use for pattern matching
Triggerpat\$	Specifies the digital pattern that must be met for the trigger to occur
Triggerwhen	Specifies the conditions under which the trigger occurs

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Disables reference triggering for the measurement or generation

**HTBasic Definition:**

```
LONG Result=FNHx_disreftrig((LONG Taskid))
```

Input:

Taskid	The task used in this function
--------	--------------------------------

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Configures the task to stop the acquisition when the device acquires all pretrigger samples, detects a rising or falling edge of a digital signal, and acquires all posttrigger samples

**HTBasic Definition:**

```
LONG Result=FNHx_cfgdiedrtrig((LONG Taskid), (Triggersrc$), (LONG Triggeredge),
                               (LONG Pretrigsam))
```

Input:

Taskid	The task used in this function
Triggersrc\$	Specifies the name of a terminal where there is a digital signal to use as the source of the trigger
Triggeredge	Specifies on which edge of the digital signal the Reference Trigger occurs
Pretrigsam	The minimum number of samples per channel to acquire before recognizing the Reference Trigger.

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Configures the task to start acquiring or generating samples when an analog signal crosses the level you specify

**HTBasic Definition:**

```
LONG Result=FNHx_cfganedrtrig((LONG Taskid), (Triggersrc$), (LONG Triggerslop),
                               (REAL Triggerlev), (LONG Pretrigsam))
```

Input:

Taskid	The task used in this function
Triggersrc\$	The name of a channel or terminal where there is an analog signal to use as the source of the trigger
Triggerslop	Specifies on which slope of the signal to start acquiring or generating samples when the signal crosses <b>Triggerlev</b>
Triggerlev	The threshold at which to start acquiring or generating samples
Pretrigsam	The minimum number of samples per channel to acquire before recognizing the Reference Trigger

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Configures the task to stop the acquisition when the device acquires all pretrigger samples, an analog signal enters or leaves a range you specify, and the device acquires all post-trigger samples

**HTBasic Definition:**

```
LONG Result=FNHx_cfganwirtrig((LONG Taskid), (Triggersrc$), (LONG Triggerwhen),
                               (REAL Windowtop), (REAL Windowbottom), (LONG Pretrigsam))
```

Input:

Taskid	The task used in this function
Triggersrc\$	The name of a virtual channel or terminal where there is an analog signal to use as the source of the trigger
Triggerwhen	Specifies whether the task starts measuring or generating samples when the signal enters the window or when it leaves the window
Windowtop	The upper limit of the window
Windowbottom	The lower limit of the window
Pretrigsam	The minimum number of samples per channel to acquire before recognizing the Reference Trigger

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Configures the task to stop the acquisition when the device acquires all pretrigger samples, matches or does not match a digital pattern, and acquires all posttrigger samples

**HTBasic Definition:**

```
LONG Result=FNHx_cfgdipartrig((LONG Taskid), (Triggersrc$), (Triggerpat$),
                               (LONG Triggerwhen), (LONG Pretrigsam))
```

Input:

Taskid	The task used in this function
Triggersrc\$	Specifies the physical channels to use for pattern matching
Triggerpat\$	Specifies the digital pattern that must be met for the trigger to occur
Triggerwhen	Specifies the conditions under which the trigger occurs
Pretrigsam	The minimum number of samples per channel to acquire before recognizing the Reference Trigger

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Disables the advance triggering for the task

**HTBasic Definition:**

```
LONG Result=FNHx_disadvtrig((LONG Taskid))
```

Input:

Taskid	The task used in this function
--------	--------------------------------

Return value: 0 = Function was successful  
<>0 = Error occurred



**Purpose:** Configures a switch task to advance to the next entry in a scan list on a rising or falling edge of a digital signal

**HTBasic Definition:**

`LONG Result=FNHx_cfgdiedatrig((LONG Taskid), (Triggersrc$), (LONG Triggeredge))`

Input:

Taskid	The task used in this function
Triggersrc\$	Specifies the name of a terminal where there is a digital signal to use as the source of the trigger
Triggeredge	Specifies on which edge of the digital signal the Reference

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Generates the specified software trigger

**HTBasic Definition:**

`LONG Result=FNHx_sendsofttrig((LONG Taskid), (LONG Triggerid))`

Input:

Taskid	The task used in this function
Triggerid	Specifies which software trigger to generate

Return value: 0 = Function was successful  
<>0 = Error occurred

## Calibration

**FNHx\_selfcal**

NI-DAQmx: **DAQmxSelfCal**

**Purpose:** Measures the onboard reference voltage of the device and adjusts the self-calibration constants to account for any errors caused by short-term fluctuations in the operating environment

### HTBasic Definition:

`LONG Result=FNHx_selfcal(Devicename$)`

#### Input:

`Devicename$`            The name of the device, as configured in MAX, to which this operation applies

Return value:    0 = Function was successful  
                  <>0 = Error occurred

### Example:

```
LONG Result
DIM Devicename$[255]
!
Devicename$="Dev1"
Result=FNHx_selfcal(Devicename$)
!
IF Result=0 THEN
    PRINT "Device '&Devicename$&' has been self-calibrated."
END IF
```

**FNHx\_perfbrofnlca**

NI-DAQmx: **DAQmxPerformBridgeOffsetNullingCal**

**Purpose:** Measures the onboard reference voltage of the device and adjusts the self-calibration constants to account for any errors caused by short-term fluctuations in the operating environment

### HTBasic Definition:

`LONG Result=FNHx_perfbrofnlca((LONG Taskid), Channel$)`

#### Input:

`Taskid`                The task used in this function  
`Channel$`              A subset of virtual channels in the task that you want to calibrate. Use this parameter if you do not want to calibrate all the channels in the task or if some channels in the task measure non-bridge-based sensors

Return value:    0 = Function was successful  
                  <>0 = Error occurred

**FNHx\_getcdevtemp**

NI-DAQmx: **DAQmxGetCalDevTemp**

**Purpose:** Returns in degrees Celsius the current temperature of the device

### HTBasic Definition:

`LONG Result=FNHx_getcdevtemp((Devicename$), REAL Data)`

#### Input:

`Devicename$`            The name of the device for this function

Return value:    0 = Function was successful  
                  **Data:** Current temperature of the device in degrees  
                  <>0 = Error occurred

**Purpose:** Returns in degrees Celsius the temperature of the device at the time of the last self calibration. Compare this temperature to the current onboard temperature to determine if you should perform another calibration.

**HTBasic Definition:**

LONG Result=FNHx\_getsfcaltemp((Devicename\$), REAL Data)

Input:

Devicename\$            The name of the device for this function

Return value:    0 = Function was successful

**Data:** Temperature of the device at the time of last calibration  
<>0 = Error occurred

**Purpose:** Returns in degrees Celsius the temperature of the device at the time of the last external calibration.

**HTBasic Definition:**

LONG Result=FNHx\_getexcaltemp((Devicename\$), REAL Data)

Input:

Devicename\$            The name of the device for this function

Return value:    0 = Function was successful

**Data:** Temperature of the device at the time of last external calibration  
<>0 = Error occurred

**Purpose:** Starts an external calibration session on a device

**HTBasic Definition:**

LONG Result=FNHx\_initextcal((Devicename\$), (Password\$), LONG Calhandle)

Input:

Devicename\$            The name of the device, as configured in MAX, to which this operation applies

Password\$              The current calibration password for the device (case sensitive)

Return value:    0 = Function was successful

**Calhandle** contains a reference to the calibration session  
<>0 = Error occurred

**Purpose:** Closes an open external calibration session

**HTBasic Definition:**

LONG Result=FNHx\_closeextcal((LONG Calhandle), (LONG Action))

Input:

Calhandle              A reference to the calibration session that you created using the DAQmxInitExtCal function

Action                 Specifies how to close the calibration session

Return value:    0 = Function was successful

<>0 = Error occurred

**Purpose:** Indicates the last date and time that the device underwent a self calibration

**HTBasic Definition:**

LONG Result=FNHx\_getscaldtati((Devname\$), LONG Year, LONG Month, LONG Day, LONG Hour, LONG Minute)

Input:

Devname\$	The name of the device for this function
Year	The last year that the device underwent a self-calibration
Month	The last month that the device underwent a self-calibration
Day	The last day that the device underwent a self-calibration
Hour	The last hour, on a 24-hour clock, that the device underwent a self-calibration
Minute	The last minute that the device underwent a self-calibration

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Indicates the last date and time that the device underwent an external calibration

**HTBasic Definition:**

LONG Result=FNHx\_getecaldtati((Devname\$), LONG Year, LONG Month, LONG Day, LONG Hour, LONG Minute)

Input:

Devname\$	The name of the device for this function
Year	The last year that the device underwent an external calibration
Month	The last month that the device underwent an external calibration
Day	The last day that the device underwent an external calibration
Hour	The last hour, on a 24-hour clock, that the device underwent an external calibration
Minute	The last minute that the device underwent an external calibration

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Sets the self calibration constants of the device to the the current external calibration constants

**HTBasic Definition:**

LONG Result=FNHx\_reslexcalcon((Devname\$))

Input:

Devname\$	The name of the device for this function
-----------	--

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Changes the external calibration password of the device

**HTBasic Definition:**

LONG Result=FNHx\_chnextcalpwd((Devname\$), (Oldpassword\$), (Newpassword\$))

Input:

Devname\$	The name of the device for this function
Oldpassword\$	The current calibration password for the device
Newpassword\$	The new calibration password for the device (four chars max.)

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Indicates in months the National Instruments recommended interval between each external calibration

**HTBasic Definition:**

LONG Result=FNHx\_getecalrecin((Devicename\$), LONG Data)

Input:

Devicename\$            The name of the device, as configured in MAX, to which this operation applies

Return value:    0 = Function was successful  
                               **Data** contains the recommended interval between each external calibration  
                               <>0 = Error occurred

**Purpose:** Adjusts the external calibration constants for the analog input section of a DSA device

**HTBasic Definition:**

LONG Result=FNHx\_adjdsaaical((LONG Calhandle), (REAL Refvoltage))

Input:

Calhandle            A reference to the calibration handle that you created using the DAQmxInitExtCal function  
 Refvoltage            The known voltage, in volts, to use as a reference for calibration

Return value:    0 = Function was successful  
                               <>0 = Error occurred

**Purpose:** Adjusts the external calibration constants for the analog output section of a DSA device

**HTBasic Definition:**

LONG Result=FNHx\_adjdsaaocal((LONG Calhandle), (LONG Channel), (REAL Reqlowvolt), (REAL Actlowvolt), (REAL Reqhigvolt), (REAL Acthigvolt), (REAL Gainset))

Input:

Calhandle            A reference to the calibration handle that you created using the DAQmxInitExtCal function  
 Channel            The number of the channel to calibrate. This number is the numeric portion of the physical channel name, not the full physical channel name  
 Reqlowvolt            The low voltage you attempted to generate at the gain setting you specified  
 Actlowvolt            The actual low voltage as measured by an external sensor  
 Reqhigvolt            The high voltage you attempted to generate at the gain setting you specified  
 Acthigvolt            The actual high voltage as measured by an external sensor  
 Gainset            The gain setting you used when you attempted to generate the requested high voltage and requested low voltage

Return value:    0 = Function was successful  
                               <>0 = Error occurred

**Purpose:** Adjusts the external calibration constant for the timebase of a DSA device with an adjustable oscillator

**HTBasic Definition:**

`LONG Result=FNHx_adjdsatbcal((LONG Calhandle), (REAL Reffreq))`

Input:

Calhandle	A reference to the calibration handle that you created using the DAQmxInitExtCal function
Reffreq	The frequency, in hertz, of the signal to use as a reference for calibration

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Adjusts the internal and external calibration constants for the SCMP pod on the PXI-4204 device

**HTBasic Definition:**

`LONG Result=FNHx_adj4204cal((LONG Calhandle), (Channames$), (REAL Lpfreq), (LONG Trackholden, (REAL Inputval))`

Input:

Calhandle	A reference to the calibration handle that you created using the DAQmxInitExtCal function
Channames\$	The physical channel(s) to calibrate
Lpfreq	The low pass cutoff frequency, in hertz, (6 or 10000) on the SCMP pod to calibrate
Trackholden	Specifies whether calibrating for trackHold is enabled or disabled
Inputval	The known voltage, in volts, to use as a reference for calibration

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Adjusts the internal and external calibration constants for the SCMP pod on the PXI-4220 device

**HTBasic Definition:**

`LONG Result=FNHx_adj4220cal((LONG Calhandle), (Channames$), (REAL Gain), (REAL Inputval))`

Input:

Calhandle	A reference to the calibration handle that you created using the DAQmxInitExtCal function
Channames\$	The physical channel(s) to calibrate
Gain	The gain value on the SCMP pod to calibrate
Inputval	The known voltage, in volts, to use as a reference for calibration

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Adjusts the internal and external calibration constants for the SCMP pod on the PXI-4224 device

**HTBasic Definition:**

```
LONG Result=FNHx_adj4224cal((LONG Calhandle), (Channames$), (REAL Gain),
                             (REAL Inputval))
```

Input:

Calhandle	A reference to the calibration handle that you created using the DAQmxInitExtCal function
Channames\$	The physical channel(s) to calibrate
Gain	The gain value on the SCMP pod to calibrate
Inputval	The known voltage, in volts, to use as a reference for calibration

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Adjusts the internal and external calibration constants for the SCMP pod on the PXI-4225 device

**HTBasic Definition:**

```
LONG Result=FNHx_adj4225cal((LONG Calhandle), (Channames$), (REAL Gain),
                             (REAL Inputval))
```

Input:

Calhandle	A reference to the calibration handle that you created using the DAQmxInitExtCal function
Channames\$	The physical channel(s) to calibrate
Gain	The gain value on the SCMP pod to calibrate
Inputval	The known voltage, in volts, to use as a reference for calibration

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Adjusts the external calibration constants on an E-Series device

**HTBasic Definition:**

```
LONG Result=FNHx_esercaladjst((LONG Calhandle), (REAL Refvoltage))
```

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration. This voltage should be between +6.000 V and +9.999 V

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Adjusts the external calibration constants for an M Series device

**HTBasic Definition:**

LONG Result=FNHx\_msercaladjst((LONG Calhandle), (REAL Refvoltage))

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration. This voltage should be between +6.000 V and +9.999 V

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Adjusts the external calibration constants for an S Series device

**HTBasic Definition:**

LONG Result=FNHx\_ssercaladjst((LONG Calhandle), (REAL Refvoltage))

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration. This voltage should be between +6.000 V and +9.999 V

Return value: 0 = Function was successful  
<>0 = Error occurred



**Purpose:** Adjusts the external calibration constants on for the baseboard of an SC Series device

**HTBasic Definition:**

LONG Result=FNHx\_scbcaladjst((LONG Calhandle), (REAL Refvoltage))

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration. This voltage should be between +6.000 V and +9.999 V

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Adjusts the external calibration constants on an AO Series device

**HTBasic Definition:**

LONG Result=FNHx\_aosecaladjst((LONG Calhandle), (REAL Refvoltage))

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration. This voltage should be between +6.000 V and +9.999 V

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Checks whether or not calibration is supported by device

**HTBasic Definition:**

LONG Result=FNHx\_devsuppcal(DeviceName\$, (LONG Calsupported))

Input:

DeviceName\$	The name of the device, as configured in MAX, to which this operation applies
--------------	---

Return value: 0 = Function was successful  
**Calsupported** returns 1 if calibration is supported, otherwise 0  
<>0 = Error occurred

**Purpose:** Checks whether or not self calibration is supported by device

**HTBasic Definition:**

LONG Result=FNHx\_selfcalsupp(DeviceName\$, (LONG Calsupported))

Input:

DeviceName\$	The name of the device, as configured in MAX, to which this operation applies
--------------	---

Return value: 0 = Function was successful  
**Calsupported** returns 1 if calibration is supported, otherwise 0  
<>0 = Error occurred

**Purpose:** Specifies the channel on the SCXI-1104 module for calibration

**HTBasic Definition:**

`LONG Result=FNHx_set1104cal((LONG Calhandle), (Channel$))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Adjusts the calibration constants for the SCXI-1104 module

**HTBasic Definition:**

`LONG Result=FNHx_adj1104cal((LONG Calhandle), (REAL Refvoltage), (REAL Measoutput))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration
Measoutput	The voltage measured at the output of the module

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Specifies the channel on the SCXI-1112 module for calibration

**HTBasic Definition:**

`LONG Result=FNHx_set1112cal((LONG Calhandle), (Channel$))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Adjusts the calibration constants for the SCXI-1112 module

**HTBasic Definition:**

`LONG Result=FNHx_adj1112cal((LONG Calhandle), (REAL Refvoltage), (REAL Measoutput))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration
Measoutput	The voltage measured at the output of the module

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Writes the specified binary value to the D/A circuitry on the specified channel at the specified range.

**HTBasic Definition:**

`LONG Result=FNHx_set1124cal((LONG Calhandle), (Channel$), (REAL Range), (REAL Daqvalue))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate
Range	The range to calibrate
Daqvalue	The binary number to write to the DAC circuitry

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Adjusts the calibration constants for the SCXI-1124 module.

**HTBasic Definition:**

`LONG Result=FNHx_adj1124cal((LONG Calhandle), (REAL Measoutput))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Measoutput	The voltage measured at the output of the module

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Sets the SCXI-1125 module to the specified gain value.

**HTBasic Definition:**

`LONG Result=FNHx_set1125cal((LONG Calhandle), (Channel$), (REAL Gain))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate
Gain	The gain value to calibrate

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Adjusts the external calibration constants for the SCXI-1125 module

**HTBasic Definition:**

`LONG Result=FNHx_adj1125cal((LONG Calhandle), (REAL Refvoltage), (REAL Measoutput))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration
Measoutput	The voltage measured at the output of the module

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Specifies the channel and upper frequency limit on the SCXI-1126 module for calibration

**HTBasic Definition:**

`LONG Result=FNHx_set1126cal((LONG Calhandle), (Channel$), (REAL Upperfreqlimit))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate
Upperfreqlimit	The high frequency limit in hertz, with 0 Hz as the low frequency limit, which most closely encapsulates the ranges to be calibrated

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Adjusts the calibration constants for the SCXI-1126 module

**HTBasic Definition:**

`LONG Result=FNHx_adj1126cal((LONG Calhandle), (REAL Reffreq), (REAL Measoutput))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Reffreq	The known frequency, in Hertz, to use as a reference for calibration
Measoutput	The voltage measured at the output of the module

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Sets the SCXI-1141 module to the specified gain value.

**HTBasic Definition:**

`LONG Result=FNHx_set1141cal((LONG Calhandle), (Channel$), (REAL Gain))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate
Gain	The gain value to calibrate

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Adjusts the external calibration constants for the SCXI-1141 module

**HTBasic Definition:**

`LONG Result=FNHx_adj1141cal((LONG Calhandle), (REAL Refvoltage), (REAL Measoutput))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration
Measoutput	The voltage measured at the output of the module

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Sets the SCXI-1142 module to the specified gain value.

**HTBasic Definition:**

`LONG Result=FNHx_set1142cal((LONG Calhandle), (Channel$), (REAL Gain))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate
Gain	The gain value to calibrate

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Adjusts the external calibration constants for the SCXI-1142 module

**HTBasic Definition:**

`LONG Result=FNHx_adj1142cal((LONG Calhandle), (REAL Refvoltage), (REAL Measoutput))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration
Measoutput	The voltage measured at the output of the module

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Sets the SCXI-1143 module to the specified gain value.

**HTBasic Definition:**

`LONG Result=FNHx_set1143cal((LONG Calhandle), (Channel$), (REAL Gain))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate
Gain	The gain value to calibrate

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Adjusts the external calibration constants for the SCXI-1143 module

**HTBasic Definition:**

`LONG Result=FNHx_adj1143cal((LONG Calhandle), (REAL Refvoltage), (REAL Measoutput))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration
Measoutput	The voltage measured at the output of the module

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Sets the SCXI-1502 module to the specified gain value.

**HTBasic Definition:**

`LONG Result=FNHx_set1502cal((LONG Calhandle), (Channel$), (REAL Gain))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate
Gain	The gain value to calibrate

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Adjusts the external calibration constants for the SCXI-1502 module

**HTBasic Definition:**

`LONG Result=FNHx_adj1502cal((LONG Calhandle), (REAL Refvoltage), (REAL Measoutput))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration
Measoutput	The voltage measured at the output of the module

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Sets the SCXI-1503 module to the specified gain value.

**HTBasic Definition:**

`LONG Result=FNHx_set1503cal((LONG Calhandle), (Channel$), (REAL Gain))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate
Gain	The gain value to calibrate

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Adjusts the external calibration constants for the SCXI-1503 module

**HTBasic Definition:**

`LONG Result=FNHx_adj1503cal((LONG Calhandle), (REAL Refvoltage), (REAL Measoutput))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration
Measoutput	The voltage measured at the output of the module

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Adjusts the external calibration constants for the SCXI-1503 module

**HTBasic Definition:**

`LONG Result=FNHx_adj1503cc1((LONG Calhandle), (Channel$), (REAL Meascurrent))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate
Meascurrent	The current measured at the output of the module

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Sets the SCXI-1520 module to the specified gain value.

**HTBasic Definition:**

`LONG Result=FNHx_set1503cal((LONG Calhandle), (Channel$), (REAL Gain))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate
Gain	The gain value to calibrate

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Adjusts the external calibration constants for the SCXI-1520 module

**HTBasic Definition:**

`LONG Result=FNHx_adj1520cal((LONG Calhandle), (REAL Refvoltage), (REAL Measoutput))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration
Measoutput	The voltage measured at the output of the module

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Sets the SCXI-1102 module to the specified gain value.

**HTBasic Definition:**

`LONG Result=FNHx_set1102cal((LONG Calhandle), (Channel$), (REAL Gain))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate
Gain	The gain value to calibrate

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Adjusts the external calibration constants for the SCXI-1102 module

**HTBasic Definition:**

`LONG Result=FNHx_adj1102cal((LONG Calhandle), (REAL Refvoltage), (REAL Measoutput))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration
Measoutput	The voltage measured at the output of the module

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Sets the SCXI-153x module to the specified gain value.

**HTBasic Definition:**

`LONG Result=FNHx_set153xcal((LONG Calhandle), (Channel$), (REAL Gain))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate
Gain	The gain value to calibrate

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Adjusts the external calibration constants for the SCXI-153x module

**HTBasic Definition:**

`LONG Result=FNHx_adj153xcal((LONG Calhandle), (REAL Refvoltage), (REAL Measoutput))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration
Measoutput	The voltage measured at the output of the module

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.



**Purpose:** Sets the SCXI-1540 module to the specified gain value.

**HTBasic Definition:**

```
LONG Result=FNHx_set1540cal((LONG Calhandle), (Channel$), (REAL Excvolt),
                             (REAL Excfreq))
```

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate
Excvolt	The RMS value of the internal AC excitation voltage
Excfreq	The frequency of the internal AC excitation voltage

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Adjusts the external calibration constants for the SCXI-1540 module

**HTBasic Definition:**

```
LONG Result=FNHx_adj1540cal((LONG Calhandle), (REAL Refvoltage), (REAL Measoutput),
                             (LONG Inputcalsrc))
```

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The excitation RMS voltage measured from the front of the module
Measoutput	The voltage measured at the output of the module
Inputcalsrc	The calibration input source selection

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Sets the SCXI-1122 module to the specified gain value.

**HTBasic Definition:**

```
LONG Result=FNHx_set1122cal((LONG Calhandle), (Channel$), (REAL Gain))
```

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate
Gain	The gain value to calibrate

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.7 and higher.

**Purpose:** Adjusts the external calibration constants for the SCXI-1122 module

**HTBasic Definition:**

`LONG Result=FNHx_adj1122cal((LONG Calhandle), (REAL Refvoltage), (REAL Measoutput))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration
Measoutput	The voltage measured at the output of the module

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.7 and higher.

**Purpose:** Sets the SCXI-1521 module to the specified gain value.

**HTBasic Definition:**

`LONG Result=FNHx_set1521cal((LONG Calhandle), (Channel$), (REAL Gain))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Channel\$	The physical channel to calibrate
Gain	The gain value to calibrate

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.7 and higher.

**Purpose:** Adjusts the external calibration constants for the SCXI-1521 module

**HTBasic Definition:**

`LONG Result=FNHx_adj1521cal((LONG Calhandle), (REAL Refvoltage), (REAL Measoutput))`

Input:

Calhandle	A reference to the calibration session that you created using the DAQmxInitExtCal function
Refvoltage	The known voltage, in volts, to use as a reference for calibration
Measoutput	The voltage measured at the output of the module

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.7 and higher.

**Purpose:** Performs a bridge offset nulling calibration on the channels in the task

**HTBasic Definition:**

`LONG Result=FNHx_perfbrofncex((LONG Taskid), Channel$, (LONG Skipsupchans))`

Input:

Taskid	The task used in this function
Channel\$	A subset of virtual channels in the task that you want to calibrate
Skipsupchans	Specifies whether or not to skip channels that do not support calibration

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.7 and higher.

**Purpose:** Performs shunt calibration for the specified channels using a strain gage sensor

**HTBasic Definition:**

```
LONG Result=FNHx_perfstrshcal((LONG Taskid), Chan$, (REAL Shuntresval),
                               (LONG Shuntresloc), (LONG Skipsupchans))
```

Input:

Taskid	The task used in this function
Chan\$	A subset of virtual channels in the task that you want to calibrate
Shuntresval	The resistance, in ohms, of the shunt resistor
Shuntresloc	The location of the shunt resistor
Skipsupchans	Specifies whether or not to skip channels that do not support calibration

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.7 and higher.

**Purpose:** Performs shunt calibration for the specified channels using a bridge sensor

**HTBasic Definition:**

```
LONG Result=FNHx_perfbrishcal((LONG Taskid), Chan$, (REAL Shuntresval),
                               (LONG Shuntresloc), (REAL Bridgeres), (LONG Skipsupchans))
```

Input:

Taskid	The task used in this function
Chan\$	A subset of virtual channels in the task that you want to calibrate
Shuntresval	The resistance, in ohms, of the shunt resistor
Shuntresloc	The location of the shunt resistor
Bridgeres	The resistance, in ohms, of the bridge sensor
Skipsupchans	Specifies whether or not to skip channels that do not support calibration

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.7 and higher.

**Purpose:** Retrieves the forward polynomial values used for calibrating the channel

**HTBasic Definition:**

```
LONG Result=FNHx_getaichcpfco((LONG Taskid), Channel$, REAL Data, (LONG Arraysizesam))
```

Input:

Taskid	The task ID for which the function should be executed
Channel\$	The physical channel to calibrate
Data	The value
Arraysizesam	Array size in samples

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Specifies the forward polynomial values used for calibrating the channel

**HTBasic Definition:**

`LONG Result=FNHx_setaihcpfco((LONG Taskid), Channel$, REAL Data, (LONG Arraysizesam))`

Input:

Taskid	The task ID for which the function should be executed
Channel\$	The physical channel to calibrate
Data	The value to calibrate
Arraysizesam	Array size in samples

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Resets the forward polynomial values used for calibrating the channel

**HTBasic Definition:**

`LONG Result=FNHx_rstaihcpfco((LONG Taskid), Channel$)`

Input:

Taskid	The task ID for which the function should be executed
Channel\$	The physical channel to calibrate

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Indicates the coefficients of a polynomial equation that NI-DAQmx uses to scale values from the native format of the device to volts. Each element of the array corresponds to a term of the equation. For example, if index two of the array is 4, the third term of the equation is  $4x^2$ . Scaling coefficients do not account for any custom scales or sensors contained by the channel.

**HTBasic Definition:**

`LONG Result=FNHx_getaidevscco((LONG Taskid), Channel$, REAL Data, (LONG Arraysizesam))`

Input:

Taskid	The task ID for which the function should be executed
Channel\$	The physical channel
Data	The value
Arraysizesam	Array size in samples

Return value: 0 = Function was successful  
<>0 = Error occurred

## NI-DAQmx Specific Attribute Get/Set/Reset Function Declarations

**FNHx\_getphyschnam**

NI-DAQmx: **DAQmxGetPhysicalChanName**

**Purpose:** Returns the physical channel name of the task

**HTBasic Definition:**

```
LONG Result=FNHx_getphyschnam((LONG Taskid), Channel$, Dat$, (LONG Bufsize))
```

**Input:**

Taskid	The task used in this function
Channel\$	Logical channel name
Bufsize	Buffer size of Dat\$

**Return value:** 0 = Function was successful  
                  **Dat\$** contains the physical channel name  
                  <>0 = Error occurred

**Example:**

```
LONG Result, Taskid, Bufsize
DIM Channel$[255], Dat$[255]
!
Bufsize=MAXLEN(Dat$)                ! buffersize
Channel$=""
Result=FNHx_getphyschnam((Taskid), Channel$, Dat$, (Bufsize))
!
IF Result=0 THEN
    PRINT "Physical channel name: ";Dat$
END IF
```

**FNHx\_setciprescal**

NI-DAQmx: **DAQmxSetCIPrescaler**

**Purpose:** Sets the Counter Input Prescaler property

**HTBasic Definition:**

```
LONG Result=FNHx_setciprescal((LONG Taskid), Channel$, (LONG Prescaler))
```

**Input:**

Taskid	The task used in this function
Channel\$	Logical channel name
Prescaler	Specifies the divisor to apply to the signal you connect to the counter source terminal

**Return value:** 0 = Function was successful  
                  <>0 = Error occurred

**Example:**

```
LONG Result, Taskid, Prescaler
DIM Channel$[255]
!
Channel$="Dev1/ctr0"
Prescaler=2                ! sets prescaler (=divisor) to 2
Result=FNHx_setciprescal((Taskid), Channel$, (Prescaler))
```

## Device Control

**FNHx\_resetdevice**

NI-DAQmx: **DAQmxResetDevice**

**Purpose:** Immediately aborts all tasks associated with a device and returns the device to an initialized state. Aborting a task stops and releases any resources the task reserved

### HTBasic Definition:

```
LONG Result=FNHx_resetdevice(Devicename$)
```

#### Input:

Devicename\$            The name of the device, as configured in MAX, to which this operation applies

Return value:    0 = Function was successful  
                  <>0 = Error occurred

### Example:

```
LONG Result
DIM Devicename$[255]
!
Devicename$="Dev1"
Result=FNHx_resetdevice(Devicename$)
!
IF Result=0 THEN
    PRINT "Device '&Devicename$&' has been reset."
END IF
```

## Switch Functions

**FNHx\_swcrscanlist**

NI-DAQmx: **DAQmxSwitchCreateScanList**

**Purpose:** Creates a new switch scanning task that uses the specified scan list and applies it to the specified task

### HTBasic Definition:

`LONG Result=FNHx_swcrscanlist((Scanlist$), LONG Taskid)`

#### Input:

Scanlist\$                Uses a special syntax to specify the sequence of connections and disconnections for the task

Return value:    0 = Function was successful  
                              **Taskid** now contains Task ID of newly created task  
                              <>0 = Error occurred

**FNHx\_swconnect**

NI-DAQmx: **DAQmxSwitchConnect**

**Purpose:** Makes a connection between two switch channels. When using this function, specify only the two connection endpoints using **Swchan1** and **Swchan2**

### HTBasic Definition:

`LONG Result=FNHx_swconnect((Swchan1$), (Swchan2$), (LONG Waitforsettling))`

#### Input:

Swchan1\$                The first channel to connect  
Swchan2\$                The second channel to connect  
Waitforsettling        If TRUE, this function waits for the switches to settle before returning. If FALSE, the function returns immediately after the operation

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**FNHx\_swconnectmul**

NI-DAQmx: **DAQmxSwitchConnectMulti**

**Purpose:** Makes one or more connections specified by a connection list. You can specify only the two endpoints, or you can specify the explicit path between two endpoints

### HTBasic Definition:

`LONG Result=FNHx_swconnectmul((Connlist$), (LONG Waitforsettling))`

#### Input:

Connlist\$                A list of connections, using a special syntax, to make between switch channels  
Waitforsettling        If TRUE, this function waits for the switches to settle before returning. If FALSE, the function returns immediately after the operation

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**FNHx\_swdisconnect**

NI-DAQmx: **DAQmxSwitchDisconnect**

**Purpose:** Disconnects two switch channels

### HTBasic Definition:

`LONG Result=FNHx_swdisconnect((Swchan1$), (Swchan2$), (LONG Waitforsettling))`

#### Input:

Swchan1\$                The first channel to disconnect  
Swchan2\$                The second channel to disconnect  
Waitforsettling        If TRUE, this function waits for the switches to settle before returning. If FALSE, the function returns immediately after the operation

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**Purpose:** Terminates all active connections on the device, which places the relays into the topology reset state

**HTBasic Definition:**

`LONG Result=FNHx_swdisconnectall((Devicename$), (LONG Waitforsettling))`

Input:

Devicename\$	The name of the device, as configured in MAX, to which this operation applies
Waitforsettling	If TRUE, this function waits for the switches to settle before returning. If FALSE, the function returns immediately after the operation

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Resets a switch device and sets its topology to the one specified with newTopology

**HTBasic Definition:**

`LONG Result=FNHx_swsettopares((Devicename$), (Newtopology$))`

Input:

Devicename\$	The name of the device, as configured in MAX, to which this operation applies
Newtopology\$	The switch topology to use on the device

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Returns information about the path between Swchan1\$ and Swchan2\$

**HTBasic Definition:**

`LONG Result=FNHx_swfindpath((Swchan1$), (Swchan2$), Path$, (LONG Pathbufsize), (LONG Pathstatus))`

Input:

Swchan1\$	The first channel to connect
Swchan2\$	The second channel to connect
Pathbufsize	The size of path. If you pass 0, this function returns the size of the buffer needed to allocate

Return value: 0 = Function was successful

<b>Pathstatus</b>	The status of the requested path
<b>Path\$</b>	The existing path or an available path between Swchan1\$ and Swchan2\$

<>0 = Error occurred

**Purpose:** Opens the specified relays

**HTBasic Definition:**

`LONG Result=FNHx_swopenrelays((Relays$), (LONG Waitforsettling))`

Input:

Relays\$	A set of relays to open
Waitforsettling	If TRUE, this function waits for the switches to settle before returning. If FALSE, the function returns immediately after the operation

Return value: 0 = Function was successful  
<>0 = Error occurred



**Purpose:** Closes the specified relays

**HTBasic Definition:**

`LONG Result=FNHx_swclosrelays((Relays$), (LONG Waitforsettling))`

Input:

Relays\$	A set of relays to close
Waitforsettling	If TRUE, this function waits for the switches to settle before returning. If FALSE, the function returns immediately after the operation

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Waits for the settling time on the device to expire

**HTBasic Definition:**

`LONG Result=FNHx_swaitforset((Devicename$))`

Input:

Devicename\$	The name of the device, as configured in MAX, to which this operation applies
--------------	---

Return value: 0 = Function was successful  
<>0 = Error occurred

## Signal Routing

**FNHx\_connterms**

NI-DAQmx: **DAQmxConnectTerms**

**Purpose:** Creates a route between a source and destination terminal

**HTBasic Definition:**

**LONG** Result=FNHx\_connterms((Srcterm\$), (Destterm\$), (**LONG** Signalmodifiers))

Input:

Srcterm\$	The originating terminal of the route. You can specify a terminal name
Destterm\$	The receiving terminal of the route. You can specify a terminal name
Signalmodifiers	Specifies whether or not to invert the signal routed from the sourceTerminal to the destinationTerminal

Return value: 0 = Function was successful  
<>0 = Error occurred

**FNHx\_disconnterms**

NI-DAQmx: **DAQmxDisconnectTerms**

**Purpose:** Removes signal routes previously created using DAQmxConnectTerms

**HTBasic Definition:**

**LONG** Result=FNHx\_disconnterms((Srcterm\$), (Destterm\$))

Input:

Srcterm\$	The originating terminal of the route. You can specify a terminal name
Destterm\$	The receiving terminal of the route. You can specify a terminal name

Return value: 0 = Function was successful  
<>0 = Error occurred

**FNHx\_tristoutterm**

NI-DAQmx: **DAQmxTristateOutputTerm**

**Purpose:** Removes signal routes previously created using DAQmxConnectTerms

**HTBasic Definition:**

**LONG** Result=FNHx\_tristoutterm((Outputterm\$))

Input:

Outputterm\$	The terminal on the I/O connector to set to high-impedance state. You can specify a terminal name
--------------	---

Return value: 0 = Function was successful  
<>0 = Error occurred

## TEDS

**FNHx\_configteds**

NI-DAQmx: **DAQmxConfigureTEDS**

**Purpose:** Associates TEDS information with the physical channel you specify. If you do not specify the filename of a data sheet in the filePath parameter, this function attempts to find a TEDS sensor connected to the physical channel. This function temporarily overrides any TEDS configuration for the physical channel that you performed in MAX

### HTBasic Definition:

**LONG** Result=FNHx\_configteds(Physchan\$, FilePath\$)

#### Input:

Physchan\$	The name of the physical channel you want to configure
FilePath\$	The path to a Virtual TEDS data sheet that you want to associate with the physical channel. If you do not specify the filename of a data sheet, this function attempts to find a TEDS sensor connected to the physical channel

Return value: 0 = Function was successful  
<>0 = Error occurred

**FNHx\_clearteds**

NI-DAQmx: **DAQmxClearTEDS**

**Purpose:** Removes TEDS information from the physical channel you specify. This function temporarily overrides any TEDS configuration for the physical channel that you performed in MAX

### HTBasic Definition:

**LONG** Result=FNHx\_clearteds(Physchan\$)

#### Input:

Physchan\$	The name of the physical channel you want to clear
------------	--

Return value: 0 = Function was successful  
<>0 = Error occurred

## Events & Signals

**FNHx\_exportsignal**

NI-DAQmx: **DAQmxExportSignal**

**Purpose:** Creates a new switch scanning task that uses the specified scan list and applies it to the specified task

### HTBasic Definition:

`LONG Result=FNHx_exportsignal((LONG Taskid), (LONG Signalid), (Outputterminal$))`

### Input:

Taskid	The task used in this function
Signalid	The name of the trigger, clock, or event to export
Outputterminal\$	The destination terminal of the exported signal

**Return value:** 0 = Function was successful  
<>0 = Error occurred

## Buffer

**FNHx\_getbufinsiz**

NI-DAQmx: **DAQmxGetBufInputBufSize**

**Purpose:** Determines the number of samples the input buffer can hold for each channel in the task. Zero indicates to allocate no buffer. Use a buffer size of 0 to perform a hardware-timed operation without using a buffer. Setting this property overrides the automatic input buffer allocation that NI-DAQmx performs

### HTBasic Definition:

`LONG Result=FNHx_getbufinsiz((LONG Taskid), LONG Bufsize)`

Input:  
Taskid                    The task used in this function

Return value:    0 = Function was successful  
                              **Bufsize**    Size of buffer in samples  
                              <>0 = Error occurred

**FNHx\_setbufinsiz**

NI-DAQmx: **DAQmxSetBufInputBufSize**

**Purpose:** Specifies the number of samples the input buffer can hold for each channel in the task. Zero indicates to allocate no buffer. Use a buffer size of 0 to perform a hardware-timed operation without using a buffer. Setting this property overrides the automatic input buffer allocation that NI-DAQmx performs

### HTBasic Definition:

`LONG Result=FNHx_setbufinsiz((LONG Taskid), LONG Bufsize)`

Input:  
Taskid                    The task used in this function  
Bufsize                   Size of buffer in samples

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**FNHx\_getbufoutsiz**

NI-DAQmx: **DAQmxGetBufOutputBufSize**

**Purpose:** Determines the number of samples the output buffer can hold for each channel in the task. Zero indicates to allocate no buffer. Use a buffer size of 0 to perform a hardware-timed operation without using a buffer. Setting this property overrides the automatic output buffer allocation that NI-DAQmx performs

### HTBasic Definition:

`LONG Result=FNHx_getbufoutsiz((LONG Taskid), LONG Bufsize)`

Input:  
Taskid                    The task used in this function

Return value:    0 = Function was successful  
                              **Bufsize**    Size of buffer in samples  
                              <>0 = Error occurred

**FNHx\_setbufoutsiz**

NI-DAQmx: **DAQmxSetBufOutputBufSize**

**Purpose:** Specifies the number of samples the output buffer can hold for each channel in the task. Zero indicates to allocate no buffer. Use a buffer size of 0 to perform a hardware-timed operation without using a buffer. Setting this property overrides the automatic output buffer allocation that NI-DAQmx performs

### HTBasic Definition:

`LONG Result=FNHx_setbufoutsiz((LONG Taskid), LONG Bufsize)`

Input:  
Taskid                    The task used in this function  
Bufsize                   Size of buffer in samples

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**Purpose:** Determines in samples per channel the size of the onboard output buffer of the device

**HTBasic Definition:**

`LONG Result=FNHx_getobufoutsz((LONG Taskid), LONG Bufsize)`

Input:

Taskid                    The task used in this function

Return value:    0 = Function was successful  
                              **Bufsize**    Size of buffer in samples  
                              <>0 = Error occurred

**Purpose:** Specifies in samples per channel the size of the onboard output buffer of the device

**HTBasic Definition:**

`LONG Result=FNHx_setobufoutsz((LONG Taskid), LONG Bufsize)`

Input:

Taskid                    The task used in this function  
Bufsize                    Size of buffer in samples

Return value:    0 = Function was successful  
                              <>0 = Error occurred

## Other Get/Set specific functions (Part I)

**FNHx\_getaimax**

NI-DAQmx: **DAQmxGetAIMax**

**Purpose:** Returns the coerced maximum value that the device can measure with the current settings

**HTBasic Definition:**

`LONG Result=FNHx_getaimax((LONG Taskid), (Channel$), REAL Aidata)`

Input:

Taskid                    The task used in this function  
Channel\$                 The name of the channel you want to query

Return value:    0 = Function was successful  
                              **Aidata**    Maximum value you expect to measure (in Units)  
                              <>0 = Error occurred

**FNHx\_getaimin**

NI-DAQmx: **DAQmxGetAIMin**

**Purpose:** Returns the coerced minimum value that the device can measure with the current settings

**HTBasic Definition:**

`LONG Result=FNHx_getaimin((LONG Taskid), (Channel$), REAL Aidata)`

Input:

Taskid                    The task used in this function  
Channel\$                 The name of the channel you want to query

Return value:    0 = Function was successful  
                              **Aidata**    Minimum value you expect to measure (in Units)  
                              <>0 = Error occurred

**FNHx\_getailpcofrq**

NI-DAQmx: **DAQmxGetAILowpassCutoffFreq**

**Purpose:** Returns the frequency in Hertz that corresponds to the -3dB cutoff of the filter

**HTBasic Definition:**

`LONG Result=FNHx_getailpcofrq((LONG Taskid), (Channel$), REAL Fdata)`

Input:

Taskid                    The task used in this function  
Channel\$                 The name of the channel you want to query

Return value:    0 = Function was successful  
                              **Fdata**    Frequency that corresponds to the -3dB cutoff of the filter  
                              <>0 = Error occurred

**FNHx\_setailpcofrq**

NI-DAQmx: **DAQmxSetAILowpassCutoffFreq**

**Purpose:** Specifies the frequency in Hertz that corresponds to the -3dB cutoff of the filter

**HTBasic Definition:**

`LONG Result=FNHx_setailpcofrq((LONG Taskid), (Channel$), (REAL Fdata))`

Input:

Taskid                    The task used in this function  
Channel\$                 The name of the channel you want to set  
Fdata                     Frequency that corresponds to the -3dB cutoff of the filter

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**Purpose:** Determines whether to enable the lowpass filter of the channel

**HTBasic Definition:**

`LONG Result=FNHx_getailwpenab((LONG Taskid), (Channel$), LONG Ldata)`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to set

Return value: 0 = Function was successful  
                   **Ldata** Specifies whether to enable the lowpass filter  
                                   of the channel  
 <>0 = Error occurred

**Purpose:** Specifies whether to enable the lowpass filter of the channel

**HTBasic Definition:**

`LONG Result=FNHx_setailwpenab((LONG Taskid), (Channel$), (LONG Ldata))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to set
Ldata	Specifies whether to enable the lowpass filter of the channel

Return value: 0 = Function was successful  
 <>0 = Error occurred

**Purpose:** Determines the hysteresis level in the units of the measurement or generation

**HTBasic Definition:**

`LONG Result=FNHx_getanedsthys((LONG Taskid), (Channel$), REAL Hdata)`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to set

Return value: 0 = Function was successful  
                   **Hdata** Hysteresis level in the units of the measurement  
                                   or generation  
 <>0 = Error occurred

**Purpose:** Specifies a hysteresis level in the units of the measurement or generation

**HTBasic Definition:**

`LONG Result=FNHx_setanedsthys((LONG Taskid), (Channel$), (REAL Hdata))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to set
Hdata	Specifies whether to enable the lowpass filter of the channel

Return value: 0 = Function was successful  
 <>0 = Error occurred



**Purpose:** Determines whether to synchronize the AC excitation source of the channel to that of another channel

**HTBasic Definition:**

`LONG Result=FNHx_getaiacexsen((LONG Taskid), (Channel$), LONG Edata)`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to get

Return value: 0 = Function was successful  
                   **Edata** Determines whether to synchronize the AC excitation source of the channel to that of another channel  
 <>0 = Error occurred

**Purpose:** Specifies whether to synchronize the AC excitation source of the channel to that of another channel

**HTBasic Definition:**

`LONG Result=FNHx_setaiacexsen((LONG Taskid), (Channel$), (LONG Edata))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to set
Edata	Specifies whether to synchronize the AC excitation source of the channel to that of another channel

Return value: 0 = Function was successful  
 <>0 = Error occurred

**Purpose:** Indicates if the device detected an overload in any channel in the task

**HTBasic Definition:**

`LONG Result=FNHx_getrdovchexi((LONG Taskid), LONG Odata)`

Input:

Taskid	The task used in this function
--------	--------------------------------

Return value: 0 = Function was successful  
                   **Odata** Indicates if the device detected an overload in any channel in the task  
 <>0 = Error occurred

**Purpose:** Indicates the names of any overloaded virtual channels in the task. You must read Overloaded Channels Exist before you read this property. Otherwise, you will receive an error

**HTBasic Definition:**

`LONG Result=FNHx_getrdovchans((LONG Taskid), Odata$, (LONG BufferSize))`

Input:

Taskid	The task used in this function
BufferSize	Size of the buffer <b>Odata\$</b>

Return value: 0 = Function was successful  
                   **Odata\$** Names of overloaded virtual channels in the task  
 <>0 = Error occurred

**Purpose:** Determines the name of a terminal where there is a digital signal to use as the source of the Pause Trigger

**HTBasic Definition:**

`LONG Result=FNHx_getdiglvptsc((LONG Taskid), Ldata$, (LONG Buffersize))`

Input:

Taskid                    The task used in this function  
Buffersize                Size of the buffer **Odata\$**

Return value:    0 = Function was successful  
                              **Ldata\$**    Determines the name of the terminal where there is a digital signal to use as the source of the Pause Trigger  
                              <>0 = Error occurred

**Purpose:** Specifies the name of a terminal where there is a digital signal to use as the source of the Pause Trigger

**HTBasic Definition:**

`LONG Result=FNHx_setdiglvptsc((LONG Taskid), Ldata$)`

Input:

Taskid                    The task used in this function  
Ldata\$                    Specifies the name of the terminal where there is a digital signal to use as the source of the Pause Trigger

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**Purpose:** Determines the type of trigger to use to pause a task

**HTBasic Definition:**

`LONG Result=FNHx_getpaustrtyp((LONG Taskid), LONG Tdata)`

Input:

Taskid                    The task used in this function

Return value:    0 = Function was successful  
                              **Tdata**    Returns the type of trigger to use to pause a task  
                              <>0 = Error occurred

**Purpose:** Specifies the type of trigger to use to pause a task

**HTBasic Definition:**

`LONG Result=FNHx_setpaustrtyp((LONG Taskid), (LONG Tdata))`

Input:

Taskid                    The task used in this function  
Tdata                      Specifies the type of trigger to use to pause a task

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**Purpose:** Determines whether the task pauses while the signal is high or low

**HTBasic Definition:**

`LONG Result=FNHx_getdlpatrwhn((LONG Taskid), LONG Wdata)`

Input:

Taskid                    The task used in this function

Return value:    0 = Function was successful

**Wdata**    Determines whether the task pauses while the signal is high or low  
<>0 = Error occurred

**Purpose:** Specifies whether the task pauses while the signal is high or low

**HTBasic Definition:**

`LONG Result=FNHx_setdlpatrwhn((LONG Taskid), (LONG Wdata))`

Input:

Taskid                    The task used in this function  
Wdata                     Specifies whether the task pauses while the signal is high or low

Return value:    0 = Function was successful

                              <>0 = Error occurred

**Purpose:** Determines the terminal to which the GPS synchronization signal is connected

**HTBasic Definition:**

`LONG Result=FNHx_getcigpssysr((LONG Taskid), (Channel$), Gdata$, (LONG BufferSize))`

Input:

Taskid                    The task used in this function  
Channel\$                  The name of the channel you want to work with  
BufferSize                Size of the buffer **Gdata\$**

Return value:    0 = Function was successful

**Gdata\$**    Determines the terminal to which the GPS synchronization signal is connected  
<>0 = Error occurred

**Purpose:** Specifies the terminal to which the GPS synchronization signal is connected

**HTBasic Definition:**

`LONG Result=FNHx_setcigpssysr((LONG Taskid), (Channel$), Gdata$)`

Input:

Taskid                    The task used in this function  
Channel\$                  The name of the channel you want to set  
Gdata\$                    Specifies the terminal to which the GPS synchronization signal is connected

Return value:    0 = Function was successful

                              <>0 = Error occurred

**Purpose:** Determines the number of seconds to wait after a handshake cycle before starting a new handshake cycle

**HTBasic Definition:**

LONG Result=FNHx\_gethsdelaxfr((LONG Taskid), REAL Ddata)

Input:

Taskid                    The task used in this function

Return value:

0 = Function was successful  
       **Ddata**    Number of seconds to wait after a handshake cycle before starting a new handshake cycle  
 <>0 = Error occurred

**Purpose:** Specifies the number of seconds to wait after a handshake cycle before starting a new handshake cycle

**HTBasic Definition:**

LONG Result=FNHx\_sethsdelaxfr((LONG Taskid), (REAL Ddata))

Input:

Taskid                    The task used in this function  
 Ddata                    Number of seconds to wait after a handshake cycle before starting a new handshake cycle

Return value:

0 = Function was successful  
 <>0 = Error occurred

**Purpose:** Determines in seconds the amount of time to wait after the Handshake Trigger asserts before deasserting the Handshake Event if Output Behavior is DAQmx\_Val\_Interlocked

**HTBasic Definition:**

LONG Result=FNHx\_getexhseiddl((LONG Taskid), REAL Idata)

Input:

Taskid                    The task used in this function

Return value:

0 = Function was successful  
       **Idata**    Returns in seconds the amount of time to wait after the Handshake Trigger asserts before deasserting the Handshake Event if Output Behavior is DAQmx\_Val\_Interlocked  
 <>0 = Error occurred

**Purpose:** Specifies in seconds the amount of time to wait after the Handshake Trigger asserts before deasserting the Handshake Event if Output Behavior is DAQmx\_Val\_Interlocked

**HTBasic Definition:**

LONG Result=FNHx\_setexhseiddl((LONG Taskid), (REAL Idata))

Input:

Taskid                    The task used in this function  
 Idata                    Specifies in seconds the amount of time to wait after the Handshake Trigger asserts before deasserting the Handshake Event if Output Behavior is DAQmx\_Val\_Interlocked

Return value:

0 = Function was successful  
 <>0 = Error occurred

**Purpose:** Determines whether to enable the digital filter for the line(s) or port(s). You can enable the filter on a line-by-line basis. You do not have to enable the filter for all lines in a channel

**HTBasic Definition:**

`LONG Result=FNHx_getdidgflten((LONG Taskid), (Channel$), LONG Fdata)`

Input:

Taskid                      The task used in this function  
Channel\$                    The name of the channel you want to work with

Return value:     0 = Function was successful  
                           **Fdata**     Determines whether to enable the digital filter for the line(s) or port(s)  
                           <>0 = Error occurred

**Purpose:** Specifies whether to enable the digital filter for the line(s) or port(s). You can enable the filter on a line-by-line basis. You do not have to enable the filter for all lines in a channel

**HTBasic Definition:**

`LONG Result=FNHx_setdidgflten((LONG Taskid), (Channel$), (LONG Fdata))`

Input:

Taskid                      The task used in this function  
Channel\$                    The name of the channel you want to work with  
Fdata                        Specifies whether to enable the digital filter for line(s) / port(s)

Return value:     0 = Function was successful  
                           <>0 = Error occurred

**Purpose:** Determines in seconds the minimum pulse width the filter recognizes as a valid high or low state transition

**HTBasic Definition:**

`LONG Result=FNHx_getdidgflmpw((LONG Taskid), (Channel$), REAL Mdata)`

Input:

Taskid                      The task used in this function  
Channel\$                    The name of the channel you want to work with

Return value:     0 = Function was successful  
                           **Mdata**     Minimum pulse width (in seconds) the filter recognizes as a valid high or low state transition  
                           <>0 = Error occurred

**Purpose:** Specifies in seconds the minimum pulse width the filter recognizes as a valid high or low state transition

**HTBasic Definition:**

`LONG Result=FNHx_setdidgflmpw((LONG Taskid), (Channel$), (REAL Mdata))`

Input:

Taskid                      The task used in this function  
Channel\$                    The name of the channel you want to work with  
Mdata                        Minimum pulse width (in seconds) the filter recognizes as a valid high or low state transition

Return value:     0 = Function was successful  
                           <>0 = Error occurred

**Purpose:** Determines the type of trigger to use to advance to the next entry in a switch scan list

**HTBasic Definition:**

`LONG Result=FNHx_getadvtrgtyp((LONG Taskid), LONG Adata)`

Input:

Taskid                    The task used in this function

Return value:

0 = Function was successful  
       **Adata**    Type of trigger to use to advance to the next entry in a  
                   switch scan list  
 <>0 = Error occurred

**Purpose:** Specifies the type of trigger to use to advance to the next entry in a switch scan list

**HTBasic Definition:**

`LONG Result=FNHx_setadvtrgtyp((LONG Taskid), (LONG Adata))`

Input:

Taskid                    The task used in this function  
 Adata                    Type of trigger to use to advance to the next entry in a switch  
                           scan list

Return value:

0 = Function was successful  
 <>0 = Error occurred

**Purpose:** Determines if the task advances through the scan list multiple times

**HTBasic Definition:**

`LONG Result=FNHx_getswscrpmo((LONG Taskid), LONG Rdata)`

Input:

Taskid                    The task used in this function

Return value:

0 = Function was successful  
       **Rdata**    Determines if the task advances through the scan list  
                   multiple times  
 <>0 = Error occurred

**Purpose:** Specifies if the task advances through the scan list multiple times

**HTBasic Definition:**

`LONG Result=FNHx_setswscrpmo((LONG Taskid), (LONG Rdata))`

Input:

Taskid                    The task used in this function  
 Rdata                    Specifies if the task advances through the scan list multiple times

Return value:

0 = Function was successful  
 <>0 = Error occurred

**Purpose:** Determines the terminal of the signal to use as the Master Timebase. On an E Series device, you can choose only between the onboard 20MHz Timebase or the RTSI7 terminal

**HTBasic Definition:**

`LONG Result=FNHx_getmastbasrc((LONG Taskid), Tdata$, (LONG Buffersize))`

Input:

Taskid	The task used in this function
Buffersize	Size of the buffer <b>Tdata\$</b>

Return value: 0 = Function was successful  
**Tdata\$** Determines the terminal of the signal to use as the Master Timebase  
 <>0 = Error occurred

**Purpose:** Specifies the terminal of the signal to use as the Master Timebase. On an E Series device, you can choose only between the onboard 20MHz Timebase or the RTSI7 terminal

**HTBasic Definition:**

`LONG Result=FNHx_setmastbasrc((LONG Taskid), Tdata$)`

Input:

Taskid	The task used in this function
Tdata\$	Specifies the terminal of the signal to use as the Master Timebase

Return value: 0 = Function was successful  
 <>0 = Error occurred

**Purpose:** Determines the rate of the Master Timebase

**HTBasic Definition:**

`LONG Result=FNHx_getmastbrate((LONG Taskid), (Channel$), REAL Mdata)`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to work with

Return value: 0 = Function was successful  
**Mdata** Determines the rate of the Master Timebase  
 <>0 = Error occurred

**Purpose:** Specifies the rate of the Master Timebase

**HTBasic Definition:**

`LONG Result=FNHx_setmastbrate((LONG Taskid), (REAL Mdata))`

Input:

Taskid	The task used in this function
Mdata	Specifies the rate of the Master Timebase

Return value: 0 = Function was successful  
 <>0 = Error occurred

**Purpose:** Determines the terminal of the signal to use as the Reference Clock

**HTBasic Definition:**

`LONG Result=FNHx_getrefclksrc((LONG Taskid), Tdata$, (LONG BufferSize))`

Input:

Taskid                    The task used in this function  
 BufferSize                Size of the buffer **Tdata\$**

Return value:    0 = Function was successful  
                               **Tdata\$**    The terminal of the signal to use as the Reference Clock  
                               <>0 = Error occurred

**Purpose:** Specifies the terminal of the signal to use as the Reference Clock

**HTBasic Definition:**

`LONG Result=FNHx_setrefclksrc((LONG Taskid), Tdata$)`

Input:

Taskid                    The task used in this function  
 Tdata\$                    The terminal of the signal to use as the Reference Clock

Return value:    0 = Function was successful  
                               <>0 = Error occurred

**Purpose:** Determines the frequency of the Reference Clock

**HTBasic Definition:**

`LONG Result=FNHx_getrefclkrat((LONG Taskid), (Channel$), REAL Fdata)`

Input:

Taskid                    The task used in this function  
 Channel\$                 The name of the channel you want to work with

Return value:    0 = Function was successful  
                               **Fdata**    Determines the frequency of the Reference Clock  
                               <>0 = Error occurred

**Purpose:** Specifies the frequency of the Reference Clock

**HTBasic Definition:**

`LONG Result=FNHx_setrefclkrat((LONG Taskid), (REAL Fdata))`

Input:

Taskid                    The task used in this function  
 Fdata                     Specifies the frequency of the Reference Clock

Return value:    0 = Function was successful  
                               <>0 = Error occurred



**Purpose:** Determines the terminal of the signal to use as the Sample Clock Timebase

**HTBasic Definition:**

LONG Result=FNIx\_getsclktbsrc((LONG Taskid), Tdata\$, (LONG BufferSize))

Input:

Taskid                   The task used in this function  
 BufferSize               Size of the buffer Tdata\$

Return value:   0 = Function was successful  
                   Tdata\$   The terminal of the signal to use as the Sample Clock  
                                   Timebase  
                   <>0 = Error occurred

**Purpose:** Specifies the terminal of the signal to use as the Sample Clock Timebase

**HTBasic Definition:**

LONG Result=FNIx\_setsclktbsrc((LONG Taskid), Tdata\$)

Input:

Taskid                   The task used in this function  
 Tdata\$                   The terminal of the signal to use as the Sample Clock Timebase

Return value:   0 = Function was successful  
                   <>0 = Error occurred

**Purpose:** Determines the terminal of the signal to use as the synchronization pulse. The synchronization pulse resets the clock dividers and the ADCs/DACs on the device

**HTBasic Definition:**

LONG Result=FNIx\_getsynpulsr((LONG Taskid), Tdata\$, (LONG BufferSize))

Input:

Taskid                   The task used in this function  
 BufferSize               Size of the buffer Tdata\$

Return value:   0 = Function was successful  
                   Tdata\$   The terminal of the signal to use as the synchronization  
                                   pulse  
                   <>0 = Error occurred

**Purpose:** Specifies the terminal of the signal to use as the synchronization pulse. The synchronization pulse resets the clock dividers and the ADCs/DACs on the device

**HTBasic Definition:**

LONG Result=FNIx\_setsynpulsr((LONG Taskid), Tdata\$)

Input:

Taskid                   The task used in this function  
 Tdata\$                   The terminal of the signal to use as the synchronization pulse

Return value:   0 = Function was successful  
                   <>0 = Error occurred

## Other Get/Set specific functions (Part II)

**FNHx\_getaivolbref**

**NI-DAQmx 8.5: DAQmxGetAIVoltagedBRef**

**Purpose:** Determines the decibel reference level in the units of the channel

**HTBasic Definition:**

`LONG Result=FNHx_getaivolbref((LONG Taskid), (Channel$), REAL Data)`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful

**Data:** Specifies the decibel reference level in the units of the channel. When you read samples as a waveform, the decibel reference level is included in the waveform attributes  
<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.5 and higher.

**FNHx\_setaivolbref**

**NI-DAQmx 8.5: DAQmxSetAIVoltagedBRef**

**Purpose:** Specifies the decibel reference level in the units of the channel

**HTBasic Definition:**

`LONG Result=FNHx_setaivolbref((LONG Taskid), (Channel$), (REAL Data))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with
Data	Specifies the decibel reference level in the units of the channel. When you read samples as a waveform, the decibel reference level is included in the waveform attributes

Return value: 0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.5 and higher.

**FNHx\_rstaivolbref**

**NI-DAQmx 8.5: DAQmxResetAIVoltagedBRef**

**Purpose:** Resets the decibel reference level in the units of the channel

**HTBasic Definition:**

`LONG Result=FNHx_rstaivolbref((LONG Taskid), (Channel$))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Determines the method or equation form that the thermocouple scale uses

**HTBasic Definition:**

`LONG Result=FNHx_getaitcsctyp((LONG Taskid), (Channel$), LONG Data)`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful  
**Data:** Specifies the method or equation form that the thermocouple scale uses

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the method or equation form that the thermocouple scale uses

**HTBasic Definition:**

`LONG Result=FNHx_setaitcsctyp((LONG Taskid), (Channel$), (LONG Data))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with
Data	Specifies the method or equation form that the thermocouple scale uses. Possible values: DAQmx_Val_Polynomial, DAQmx_Val_Table

Return value: 0 = Function was successful  
 <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the method or equation form that the thermocouple scale uses

**HTBasic Definition:**

`LONG Result=FNHx_rstaitcsctyp((LONG Taskid), (Channel$))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful  
 <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the decibel reference level in the units of the channel

**HTBasic Definition:**

`LONG Result=FNIHx_getaispbyref((LONG Taskid), (Channel$), REAL Data)`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful

**Data:** Specifies the decibel reference level in the units of the channel. When you read samples as a waveform, the decibel reference level is included in the waveform attributes

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Specifies the decibel reference level in the units of the channel

**HTBasic Definition:**

`LONG Result=FNIHx_setaispbyref((LONG Taskid), (Channel$), (REAL Data))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with
Data	Specifies the decibel reference level in the units of the channel. When you read samples as a waveform, the decibel reference level is included in the waveform attributes

Return value: 0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Resets the decibel reference level in the units of the channel

**HTBasic Definition:**

`LONG Result=FNIHx_rstaispbyref((LONG Taskid), (Channel$))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Determines the decibel reference level in the units of the channel

**HTBasic Definition:**

`LONG Result=FNHx_getaiacbyref((LONG Taskid), (Channel$), REAL Data)`

Input:

Taskid                    The task to which to add the channels that this function creates  
Channel\$                 A subset of virtual channels in the task that you want to work with

Return value:    0 = Function was successful

**Data:** Specifies the decibel reference level in the units of the channel. When you read samples as a waveform, the decibel reference level is included in the waveform attributes

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Specifies the decibel reference level in the units of the channel

**HTBasic Definition:**

`LONG Result=FNHx_setaiacbyref((LONG Taskid), (Channel$), (REAL Data))`

Input:

Taskid                    The task to which to add the channels that this function creates  
Channel\$                 A subset of virtual channels in the task that you want to work with  
Data                       Specifies the decibel reference level in the units of the channel. When you read samples as a waveform, the decibel reference level is included in the waveform attributes

Return value:    0 = Function was successful

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Resets the decibel reference level in the units of the channel

**HTBasic Definition:**

`LONG Result=FNHx_rstaiacbyref((LONG Taskid), (Channel$))`

Input:

Taskid                    The task to which to add the channels that this function creates  
Channel\$                 A subset of virtual channels in the task that you want to work with

Return value:    0 = Function was successful

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Indicates if the virtual channel was initialized using a TEDS bitstream from the corresponding physical channel

**HTBasic Definition:**

`LONG Result=FNHx_getaiisted((LONG Taskid), (Channel$), LONG Data)`

Input:

Taskid                    The task to which to add the channels that this function creates  
Channel\$                 The source terminal of the Sample Clock

Return value:    0 = Function was successful

**Data:** Indicates if the virtual channel was initialized using a TEDS bitstream from the corresponding physical channel

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the ADC Timing mode, controlling the tradeoff between speed and effective resolution

**HTBasic Definition:**

`LONG Result=FNIHx_getaiadctimd((LONG Taskid), (Channel$), LONG Data)`

Input:

Taskid                    The task to which to add the channels that this function creates  
Channel\$                 A subset of virtual channels in the task that you want to work with

Return value:    0 = Function was successful

**Data:** Specifies the ADC timing mode, controlling the tradeoff between speed and effective resolution. Some ADC timing modes provide increased powerline noise rejection.

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Specifies the ADC Timing mode, controlling the tradeoff between speed and effective resolution

**HTBasic Definition:**

`LONG Result=FNIHx_setaiadctimd((LONG Taskid), (Channel$), (LONG Data))`

Input:

Taskid                    The task to which to add the channels that this function creates  
Channel\$                 A subset of virtual channels in the task that you want to work with  
Data                       Specifies the ADC timing mode, controlling the tradeoff between speed and effective resolution. Some ADC timing modes provide increased powerline noise rejection

Return value:    0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Resets the ADC Timing mode, controlling the tradeoff between speed and effective resolution

**HTBasic Definition:**

`LONG Result=FNIHx_rstaiadctimd((LONG Taskid), (Channel$))`

Input:

Taskid                    The task to which to add the channels that this function creates  
Channel\$                 A subset of virtual channels in the task that you want to work with

Return value:    0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Specifies the logic family to use for acquisition

**HTBasic Definition:**

`LONG Result=FNIHx_getdilogifam((LONG Taskid), (Channel$), LONG Data)`

Input:

Taskid                    The task to which to add the channels that this function creates  
Channel\$                 A subset of virtual channels in the task that you want to work with

Return value:    0 = Function was successful

**Data:** Specifies the logic family to use for acquisition. A logic family corresponds to voltage thresholds that are compatible with a group of voltage standards.

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the logic family to use for acquisition

**HTBasic Definition:**

`LONG Result=FNIHx_setdilogifam((LONG Taskid), (Channel$), (LONG Data))`

Input:

Taskid                    The task to which to add the channels that this function creates  
Channel\$                 A subset of virtual channels in the task that you want to work with  
Data                       Specifies the logic family to use for acquisition. A logic family corresponds to voltage thresholds that are compatible with a group of voltage standards

Return value:    0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the logic family to use for acquisition

**HTBasic Definition:**

`LONG Result=FNIHx_rstdilogifam((LONG Taskid), (Channel$))`

Input:

Taskid                    The task to which to add the channels that this function creates  
Channel\$                 A subset of virtual channels in the task that you want to work with

Return value:    0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines if NI-DAQmx maps hardware registers to the memory space of the customer process, if possible

**HTBasic Definition:**

`LONG Result=FNHx_getdimmapena((LONG Taskid), (Channel$), LONG Data)`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful

**Data:** Specifies if NI-DAQmx maps hardware registers to the memory space of the customer process, if possible.

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies if NI-DAQmx maps hardware registers to the memory space of the customer process, if possible

**HTBasic Definition:**

`LONG Result=FNHx_setdimmapena((LONG Taskid), (Channel$), (LONG Data))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with
Data	Specifies if NI-DAQmx maps hardware registers to the memory space of the customer process, if possible.

Return value: 0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets NI-DAQmx hardware register mappings, if possible

**HTBasic Definition:**

`LONG Result=FNHx_rstdimmapena((LONG Taskid), (Channel$))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.



**Purpose:** Determines on which edge of the sample clock to acquire samples

**HTBasic Definition:**

`LONG Result=FNHx_getdiacquion((LONG Taskid), (Channel$), LONG Data)`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful  
                   Data: Specifies on which edge of the sample clock to acquire samples.

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies on which edge of the sample clock to acquire samples

**HTBasic Definition:**

`LONG Result=FNHx_setdiacquion((LONG Taskid), (Channel$), (LONG Data))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with
Data	Specifies on which edge of the sample clock to acquire samples.

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets behaviour on which edge of the sample clock to acquire samples

**HTBasic Definition:**

`LONG Result=FNHx_rstdiacquion((LONG Taskid), (Channel$))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the state of the lines in a digital output task when the task starts

**HTBasic Definition:**

`LONG Result=FNIx_getdolsstrts((LONG Taskid), (Channel$), LONG Data)`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful

**Data:** Specifies the state of the lines in a digital output task when the task starts.

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the state of the lines in a digital output task when the task starts

**HTBasic Definition:**

`LONG Result=FNIx_setdolsstrts((LONG Taskid), (Channel$), (LONG Data))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with
Data	Specifies the state of the lines in a digital output task when the task starts.

Return value: 0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the state of the lines in a digital output task when the task starts to default value

**HTBasic Definition:**

`LONG Result=FNIx_rstdolsstrts((LONG Taskid), (Channel$))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the state of the lines in a digital output task when the task pauses

**HTBasic Definition:**

`LONG Result=FNIx_getdolspaus((LONG Taskid), (Channel$), LONG Data)`

Input:

Taskid                    The task to which to add the channels that this function creates  
Channel\$                 A subset of virtual channels in the task that you want to work with

Return value:    0 = Function was successful

**Data:** Specifies the state of the lines in a digital output task when the task pauses.

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the state of the lines in a digital output task when the task pauses

**HTBasic Definition:**

`LONG Result=FNIx_setdolspaus((LONG Taskid), (Channel$), (LONG Data))`

Input:

Taskid                    The task to which to add the channels that this function creates  
Channel\$                 A subset of virtual channels in the task that you want to work with  
Data                       Specifies the state of the lines in a digital output task when the task pauses.

Return value:    0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the state of the lines in a digital output task when the task pauses to default value

**HTBasic Definition:**

`LONG Result=FNIx_rstdolspaus((LONG Taskid), (Channel$))`

Input:

Taskid                    The task to which to add the channels that this function creates  
Channel\$                 A subset of virtual channels in the task that you want to work with

Return value:    0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the state of the lines in a digital output task when the task completes execution

**HTBasic Definition:**

`LONG Result=FNIHx_getdolsdone((LONG Taskid), (Channel$), LONG Data)`

**Input:**

Taskid                   The task to which to add the channels that this function creates  
Channel\$                A subset of virtual channels in the task that you want to work with

**Return value:** 0 = Function was successful

**Data:** Specifies the state of the lines in a digital output task when the task completes execution.

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the state of the lines in a digital output task when the task completes execution

**HTBasic Definition:**

`LONG Result=FNIHx_setdolsdone((LONG Taskid), (Channel$), (LONG Data))`

**Input:**

Taskid                   The task to which to add the channels that this function creates  
Channel\$                A subset of virtual channels in the task that you want to work with  
Data                    Specifies the state of the lines in a digital output task when the task completes execution.

**Return value:** 0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the state of the lines in a digital output task when the task completes execution to default value

**HTBasic Definition:**

`LONG Result=FNIHx_rstdolsdone((LONG Taskid), (Channel$))`

**Input:**

Taskid                   The task to which to add the channels that this function creates  
Channel\$                A subset of virtual channels in the task that you want to work with

**Return value:** 0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Sets the digital logic family to use when the device powers up

**HTBasic Definition:**

`LONG Result=FNIHx_setdlfamput((Devicename$), (LONG Logicfamily))`

**Input:**

Devicename\$            The name of the device to which this operation applies  
Logicfamily            Specifies the logic family to set the device to when it powers up

**Return value:** 0 = Function was successful

<>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the logic family to use for generation

**HTBasic Definition:**

`LONG Result=FNHx_getdologifam((LONG Taskid), (Channel$), LONG Data)`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful  
                   Data: Specifies the logic family to use for generation  
                   <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the logic family to use for generation

**HTBasic Definition:**

`LONG Result=FNHx_setdologifam((LONG Taskid), (Channel$), (LONG Data))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with
Data	Specifies the logic family to use for generation

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the logic family to use for generation to default value

**HTBasic Definition:**

`LONG Result=FNHx_rstdologifam((LONG Taskid), (Channel$))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines if NI-DAQmx maps hardware registers to the memory space of the customer process, if possible

**HTBasic Definition:**

`LONG Result=FNHx_getdommapena((LONG Taskid), (Channel$), LONG Data)`

Input:

Taskid                    The task to which to add the channels that this function creates  
Channel\$                 A subset of virtual channels in the task that you want to work with

Return value:    0 = Function was successful

**Data:** Specifies if NI-DAQmx maps hardware registers to the memory space of the customer process, if possible

                              <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies if NI-DAQmx maps hardware registers to the memory space of the customer process, if possible

**HTBasic Definition:**

`LONG Result=FNHx_setdommapena((LONG Taskid), (Channel$), (LONG Data))`

Input:

Taskid                    The task to which to add the channels that this function creates  
Channel\$                 A subset of virtual channels in the task that you want to work with  
Data                       Specifies if NI-DAQmx maps hardware registers to the memory space of the customer process, if possible

Return value:    0 = Function was successful

                              <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets NI-DAQmx hardware registers mappings to default values

**HTBasic Definition:**

`LONG Result=FNHx_rstdommapena((LONG Taskid), (Channel$))`

Input:

Taskid                    The task to which to add the channels that this function creates  
Channel\$                 A subset of virtual channels in the task that you want to work with

Return value:    0 = Function was successful

                              <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Determines on which edge of the sample clock to generate samples

**HTBasic Definition:**

`LONG Result=FNHx_getdogeneron((LONG Taskid), (Channel$), LONG Data)`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful  
                   Data: Edge of the sample clock on which samples are generated  
                   <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies on which edge of the sample clock to generate samples

**HTBasic Definition:**

`LONG Result=FNHx_setdogeneron((LONG Taskid), (Channel$), (LONG Data))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with
Data	Specifies the edge of the sample clock on which samples are generated

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets behaviour on which edge of the sample clock to generate samples to default value

**HTBasic Definition:**

`LONG Result=FNHx_rstdogeneron((LONG Taskid), (Channel$))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines constraints to apply when the counter generates pulses

**HTBasic Definition:**

`LONG Result=FNHx_getcocongmod((LONG Taskid), (Channel$), LONG Data)`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful  
                   Data: Constraints to apply when the counter generates pulses  
                   <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Specifies constraints to apply when the counter generates pulses

**HTBasic Definition:**

`LONG Result=FNHx_setcocongmod((LONG Taskid), (Channel$), (LONG Data))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with
Data	Specifies constraints to apply when the counter generates pulses

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Resets the constraints to apply when the counter generates pulses

**HTBasic Definition:**

`LONG Result=FNHx_rstcocongmod((LONG Taskid), (Channel$))`

Input:

Taskid	The task to which to add the channels that this function creates
Channel\$	A subset of virtual channels in the task that you want to work with

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.5 and higher.



**Purpose:** Determines in seconds the amount of time to offset the exported Sample clock

**HTBasic Definition:**

LONG Result=FNIHx\_getexscdloff((LONG Taskid), REAL Data)

Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful

**Data:** Returns in seconds the amount of time to offset the exported sample clock

                              <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Specifies in seconds the amount of time to offset the exported Sample clock

**HTBasic Definition:**

LONG Result=FNIHx\_setexscdloff((LONG Taskid), (REAL Data))

Input:

Taskid                    The task to which to add the channels that this function creates

Data                        Specifies the amount of time to offset the exported Sample clock

Return value:    0 = Function was successful

                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Resets the amount of time to offset the exported Sample clock to default value

**HTBasic Definition:**

LONG Result=FNIHx\_rstexscdloff((LONG Taskid))

Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful

                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Specifies the name of the physical channel upon which this virtual channel is based

**HTBasic Definition:**

LONG Result=FNIHx\_setphyschnnam((LONG Taskid), (Channel\$), (Data\$))

Input:

Taskid                    The task to which to add the channels that this function creates

Channel\$                    A subset of virtual channels in the task that you want to work with

Data\$                        The name of the physical channel upon which this virtual channel is based

Return value:    0 = Function was successful

                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the terminal to which to route the Pause Trigger

**HTBasic Definition:**

`LONG Result=FNHx_getexptouttm((LONG Taskid), Data$, (LONG Buffersize))`

Input:

Taskid                    The task to which to add the channels that this function creates  
Buffersize                Size of Data\$ string buffer

Return value:    0 = Function was successful  
                              **Data\$:** Terminal to which to route the Pause Trigger  
                              <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.35 and higher.

**Purpose:** Specifies the terminal to which to route the Pause Trigger

**HTBasic Definition:**

`LONG Result=FNHx_setexptouttm((LONG Taskid), (Data$))`

Input:

Taskid                    The task to which to add the channels that this function creates  
Data\$                      Terminal to which the Pause Trigger is routed

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets to the default terminal to which to route the Pause Trigger

**HTBasic Definition:**

`LONG Result=FNHx_rstexptouttm((LONG Taskid))`

Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the active level of the exported Pause Trigger

**HTBasic Definition:**

LONG Result=FNHx\_getexptlaclv((LONG Taskid), LONG Data)

Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful

**Data:** Active level of the exported Pause Trigger

                              <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the active level of the exported Pause Trigger

**HTBasic Definition:**

LONG Result=FNHx\_setexptlaclv((LONG Taskid), (LONG Data))

Input:

Taskid                    The task to which to add the channels that this function creates

Data                      Specifies active level of the exported Pause Trigger

Return value:    0 = Function was successful

                              <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the active level of the exported Pause Trigger to it's default state

**HTBasic Definition:**

LONG Result=FNHx\_rstexptlaclv((LONG Taskid))

Input:

Taskid                    The task to which to add the channels that this function creates

Channel\$                 A subset of virtual channels in the task that you want to work with

Return value:    0 = Function was successful

                              <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the polarity of the exported Reference Trigger

**HTBasic Definition:**

LONG Result=FNHx\_getexrftrppo((LONG Taskid), LONG Data)

Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful

**Data:** The polarity of the exported Reference Trigger

                    <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies polarity of the exported Reference Trigger

**HTBasic Definition:**

LONG Result=FNHx\_setexrftrppo((LONG Taskid), (LONG Data))

Input:

Taskid                    The task to which to add the channels that this function creates

Data                      Specifies the polarity of the exported Reference Trigger

Return value:    0 = Function was successful

                    <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the polarity of the exported Reference Trigger to it's default state

**HTBasic Definition:**

LONG Result=FNHx\_rstexrftrppo((LONG Taskid))

Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful

                    <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the polarity of the exported Start Trigger

**HTBasic Definition:**

LONG Result=FNHx\_getexsttrppo((LONG Taskid), LONG Data)

Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful

**Data:** The polarity of the exported Start Trigger

                    <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies polarity of the exported Start Trigger

**HTBasic Definition:**

LONG Result=FNHx\_setexsttrppo((LONG Taskid), (LONG Data))

Input:

Taskid                    The task to which to add the channels that this function creates

Data                      Specifies the polarity of the exported Start Trigger

Return value:    0 = Function was successful

                    <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the polarity of the exported Start Trigger to it's default state

**HTBasic Definition:**

LONG Result=FNHx\_rstexsttrppo((LONG Taskid))

Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful

                    <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**FNHx\_getexrdfxedc** NI-DAQmx 8.3: **DAQmxGetExportedRdyForXferEventDeassertCond**

**Purpose:** Determines when the ready for transfer event deasserts

**HTBasic Definition:**

`LONG Result=FNHx_getexrdfxedc((LONG Taskid), LONG Data)`

Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful

**Data:** Returns when the ready for transfer event deasserts

                    <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**FNHx\_setexrdfxedc** NI-DAQmx 8.3: **DAQmxSetExportedRdyForXferEventDeassertCond**

**Purpose:** Specifies when the ready for transfer event deasserts

**HTBasic Definition:**

`LONG Result=FNHx_setexrdfxedc((LONG Taskid), (LONG Data))`

Input:

Taskid                    The task to which to add the channels that this function creates

Data                      Specifies when the ready for transfer event deasserts

Return value:    0 = Function was successful

                    <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**FNHx\_rstexrdfxedc** NI-DAQmx 8.3: **DAQmxResetExportedRdyForXferEventDeassertCond**

**Purpose:** Resets the ready for transfer event deassertion to it's default state

**HTBasic Definition:**

`LONG Result=FNHx_rstexrdfxedc((LONG Taskid))`

Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful

                    <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**FNHx\_getexrdyfxed**

NI-DAQmx 8.3: DAQmxGetExportedRdyForXferEventDeassertCondCustomThreshold

**Purpose:** Determines in samples the threshold below which the Ready for Transfer Event deasserts

**HTBasic Definition:**`LONG Result=FNHx_getexrdyfxed((LONG Taskid), LONG Data)`Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful

**Data:** The treshhold below which the Ready Transition deasserts

                    <>0 = Error occured

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**FNHx\_setexrdyfxed**

NI-DAQmx 8.3: DAQmxSetExportedRdyForXferEventDeassertCondCustomThreshold

**Purpose:** Specifies in samples the threshold below which the Ready for Transfer Event deasserts

**HTBasic Definition:**`LONG Result=FNHx_setexrdyfxed((LONG Taskid), (LONG Data))`Input:

Taskid                    The task to which to add the channels that this function creates

Data                      Specifies the treshhold below which the Ready Transition deasserts

Return value:    0 = Function was successful

                    <>0 = Error occured

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**FNHx\_rstexrdyfxed**

NI-DAQmx 8.3: DAQmxResetExportedRdyForXferEventDeassertCondCustomThreshold

**Purpose:** Resets the threshold below which the Ready for Transfer Event deasserts to it's default state

**HTBasic Definition:**`LONG Result=FNHx_rstexrdyfxed((LONG Taskid))`Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful

                    <>0 = Error occured

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the terminal to which to export the Data Active Event

**HTBasic Definition:**

`LONG Result=FNHx_getexdaceout((LONG Taskid), Data$, (LONG Buffersize))`

Input:

Taskid                    The task to which to add the channels that this function creates  
Buffersize                Size of Data\$ string buffer

Return value:    0 = Function was successful  
                              **Data\$:** The terminal to which to export the Data Active Event  
                              <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the terminal to which to export the Data Active Event

**HTBasic Definition:**

`LONG Result=FNHx_setexdaceout((LONG Taskid), (Data$))`

Input:

Taskid                    The task to which to add the channels that this function creates  
Data\$                      Terminal to which to export the Data Active Event

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the terminal to which to export the Data Active Event to it's default state

**HTBasic Definition:**

`LONG Result=FNHx_rstexdaceout((LONG Taskid))`

Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.



**FNHx\_getexdaevl1** NI-DAQmx 8.3: **DAQmxGetExportedDataActiveEventLvlActiveLvl**

**Purpose:** Determines the polarity of the exported Data Active Event

**HTBasic Definition:**

`LONG Result=FNHx_getexdaevl1((LONG Taskid), LONG Data)`

Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful

**Data:** The polarity of the exported Data Active Event

                    <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**FNHx\_setexdaevl1** NI-DAQmx 8.3: **DAQmxSetExportedDataActiveEventLvlActiveLvl**

**Purpose:** Specifies the polarity of the exported Data Active Event

**HTBasic Definition:**

`LONG Result=FNHx_setexdaevl1((LONG Taskid), (LONG Data))`

Input:

Taskid                    The task to which to add the channels that this function creates

Data                      Specifies the polarity of the exported Data Active Event

Return value:    0 = Function was successful

                    <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**FNHx\_rstexdaevl1** NI-DAQmx 8.3: **DAQmxResetExportedDataActiveEventLvlActiveLvl**

**Purpose:** Resets the polarity of the exported Data Active Event to it's default state

**HTBasic Definition:**

`LONG Result=FNHx_rstexdaevl1((LONG Taskid))`

Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful

                    <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the terminal to which to route the Ready for Start Event

**HTBasic Definition:**

`LONG Result=FNHx_getexrfseotm((LONG Taskid), Data$, (LONG Buffersize))`

Input:

Taskid                    The task to which to add the channels that this function creates  
Buffersize                Size of Data\$ string buffer

Return value:    0 = Function was successful  
                              **Data\$:** The terminal to which to route the Ready for Start Event  
                              <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the terminal to which to route the Ready for Start Event

**HTBasic Definition:**

`LONG Result=FNHx_setexrfseotm((LONG Taskid), (Data$))`

Input:

Taskid                    The task to which to add the channels that this function creates  
Data\$                      The terminal to which the Ready for Start Event is routed

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the terminal to which to route the Ready for Start Event to it's default state

**HTBasic Definition:**

`LONG Result=FNHx_rstexrfseotm((LONG Taskid))`

Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**FNHx\_getexrfsselal** NI-DAQmx 8.3: **DAQmxGetExportedRdyForStartEventLvlActiveLvl**

**Purpose:** Determines the polarity of the exported Ready for Start Event

**HTBasic Definition:**

`LONG Result=FNHx_getexrfsselal((LONG Taskid), LONG Data)`

Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful

**Data:** The polarity of the exported Ready for Start event

                  <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**FNHx\_setexrfsselal** NI-DAQmx 8.3: **DAQmxSetExportedRdyForStartEventLvlActiveLvl**

**Purpose:** Specifies the polarity of the exported Ready for Start Event

**HTBasic Definition:**

`LONG Result=FNHx_setexrfsselal((LONG Taskid), (LONG Data))`

Input:

Taskid                    The task to which to add the channels that this function creates

Data                      Specifies the polarity of the exported Ready for Start event

Return value:    0 = Function was successful

                  <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**FNHx\_rstexrfsselal** NI-DAQmx 8.3: **DAQmxResetExportedRdyForStartEventLvlActiveLvl**

**Purpose:** Resets the polarity of the exported Ready for Start Event to it's default state

**HTBasic Definition:**

`LONG Result=FNHx_rstexrfsselal((LONG Taskid))`

Input:

Taskid                    The task to which to add the channels that this function creates

Return value:    0 = Function was successful

                  <>0 = Error occurred

**Note:**

This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the product category of the device

**HTBasic Definition:**

LONG Result=FNHx\_getdevprdcats((Device\$), LONG Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful  
                               **Data:** The product category of the device  
                               <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates an array containing the names of the modules in the chassis

**HTBasic Definition:**

LONG Result=FNHx\_getdevchmdnm((Device\$), Data\$, (LONG Buffersize))

Input:

Device\$                    The device to be used for this function  
 Buffersize                Size of Data\$ string buffer

Return value:    0 = Function was successful  
                               **Data\$:** Names of the modules in the chassis  
                               <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates if the device supports analog triggering

**HTBasic Definition:**

LONG Result=FNHx\_getdevatrsup((Device\$), LONG Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful  
                               **Data:** Indicates if the device supports analog triggering  
                               <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates if the device supports digital triggering

**HTBasic Definition:**

LONG Result=FNHx\_getdevdtrsup((Device\$), LONG Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful  
                               **Data:** Indicates if the device supports digital triggering  
                               <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the maximum rate for an analog input task if the task contains only a single channel from this device

**HTBasic Definition:**

LONG Result=FNHx\_getdevaimscr((Device\$), REAL Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Maximum rate for an analog input task  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the maximum rate for an analog input task if the task contains multiple channels from this device

**HTBasic Definition:**

LONG Result=FNHx\_getdevaimmcr((Device\$), REAL Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Maximum rate for an analog input task  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the minimum rate for an analog input task on this device

**HTBasic Definition:**

LONG Result=FNHx\_getdevaimrat((Device\$), REAL Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Minimum rate for an analog input task  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates if the device supports simultaneous sampling

**HTBasic Definition:**

LONG Result=FNHx\_getdevasssup((Device\$), LONG Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates if the device supports simultaneous sampling  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the triggers supported by this device for an analog input task

**HTBasic Definition:**

LONG Result=FNHx\_getdevaitrus((Device\$), LONG Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates the triggers supported by this device for an analog input task

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates pairs of input voltage ranges supported by this device

**HTBasic Definition:**

LONG Result=FNHx\_getdevavrang((Device\$), REAL Data(\*), (LONG Arraysize))

Input:

Device\$                    The device to be used for this function

Arraysize                  Size of the Data array in samples

Return value:    0 = Function was successful

**Data (\*):** Pairs of input voltage ranges (low / high value etc.)

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the set of discrete internal voltage excitation values supported by this device

**HTBasic Definition:**

LONG Result=FNHx\_getdevaviedv((Device\$), REAL Data(\*), (LONG Arraysize))

Input:

Device\$                    The device to be used for this function

Arraysize                  Size of the Data array in samples

Return value:    0 = Function was successful

**Data (\*):** Set of discrete internal voltage excitation values

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates pairs of internal voltage excitation ranges supported by this device

**HTBasic Definition:**

LONG Result=FNHx\_getdevavierv((Device\$), REAL Data(\*), (LONG Arraysize))

Input:

Device\$                    The device to be used for this function

Arraysize                  Size of the Data array in samples

Return value:    0 = Function was successful

**Data (\*):** Pairs of internal voltage excitation ranges

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the pairs of current input ranges supported by this device

**HTBasic Definition:**

`LONG Result=FNHx_getdevacurra((Device$), REAL Data(*), (LONG Arraysize))`

Input:

Device\$                    The device to be used for this function  
 Arraysize                Size of the Data array in samples

Return value:    0 = Function was successful  
                               **Data(\*)**: Pairs of current input ranges (low / high value etc.)  
                               <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the set of discrete internal current excitation values supported by this device

**HTBasic Definition:**

`LONG Result=FNHx_getdevaciedv((Device$), REAL Data(*), (LONG Arraysize))`

Input:

Device\$                    The device to be used for this function  
 Arraysize                Size of the Data array in samples

Return value:    0 = Function was successful  
                               **Data(\*)**: Set of discrete internal current excitation values  
                               <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the pairs of frequency input ranges supported by this device

**HTBasic Definition:**

`LONG Result=FNHx_getdevafrqrn((Device$), REAL Data(*), (LONG Arraysize))`

Input:

Device\$                    The device to be used for this function  
 Arraysize                Size of the Data array in samples

Return value:    0 = Function was successful  
                               **Data(\*)**: Pairs of frequency input ranges (low / high value etc.)  
                               <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the input gain settings supported by this device

**HTBasic Definition:**

`LONG Result=FNHx_getdevaigain((Device$), REAL Data(*), (LONG Arraysize))`

Input:

Device\$                    The device to be used for this function  
 Arraysize                Size of the Data array in samples

Return value:    0 = Function was successful  
                               **Data(\*)**: Input gain settings  
                               <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the coupling types supported by this device

**HTBasic Definition:**

LONG Result=FNHx\_getdevaicoup((Device\$), LONG Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates the coupling types supported by this device  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the set of discrete lowpass cutoff frequencies supported by this device

**HTBasic Definition:**

LONG Result=FNHx\_getdevalcfdv((Device\$), REAL Data(\*), (LONG Arraysize))

Input:

Device\$                    The device to be used for this function  
Arraysize                  Size of the Data array in samples

Return value:    0 = Function was successful

**Data(\*):** Set of discrete lowpass cutoff frequencies  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates pairs of lowpass cutoff frequency ranges supported by this device

**HTBasic Definition:**

LONG Result=FNHx\_getdevalcfrv((Device\$), REAL Data(\*), (LONG Arraysize))

Input:

Device\$                    The device to be used for this function  
Arraysize                  Size of the Data array in samples

Return value:    0 = Function was successful

**Data(\*):** Pairs of lowpass cutoff frequency ranges  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates an array containing the names of the analog output physical channels available on the device

**HTBasic Definition:**

LONG Result=FNHx\_getdevaophch((Device\$), Data\$, (LONG Buffersize))

Input:

Device\$                    The device to be used for this function  
Buffersize                 Size of Data\$ string buffer

Return value:    0 = Function was successful

**Data\$:** Array with names of analog output physical channels  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.



**Purpose:** Indicates if the device supports the sample clock timing type for analog output tasks

**HTBasic Definition:**

LONG Result=FNHx\_getdevaosdsp((Device\$), LONG Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates if the device supports the sample clock timing type for analog output tasks

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the maximum analog output rate of the device

**HTBasic Definition:**

LONG Result=FNHx\_getdevaomart((Device\$), REAL Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates the maximum analog output rate

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the minimum analog output rate of the device

**HTBasic Definition:**

LONG Result=FNHx\_getdevaomirt((Device\$), REAL Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates the minimum analog output rate

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the triggers supported by this device for analog output tasks

**HTBasic Definition:**

LONG Result=FNHx\_getdevaotrus((Device\$), LONG Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates the triggers supported by this device for analog output tasks

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates pairs of output voltage ranges supported by this device

**HTBasic Definition:**

LONG Result=FNHx\_getdevaovorg((Device\$), REAL Data(\*), (LONG Arraysize))

Input:

Device\$                    The device to be used for this function  
 Arraysize                 Size of the Data array in samples

Return value:    0 = Function was successful  
                               **Data(\*)**: Pairs of output voltage ranges supported by this device  
                               <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates pairs of output current ranges supported by this device

**HTBasic Definition:**

LONG Result=FNHx\_getdevaocurg((Device\$), REAL Data(\*), (LONG Arraysize))

Input:

Device\$                    The device to be used for this function  
 Arraysize                 Size of the Data array in samples

Return value:    0 = Function was successful  
                               **Data(\*)**: Pairs of output current ranges supported by this device  
                               <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the output gain settings supported by this device

**HTBasic Definition:**

LONG Result=FNHx\_getdevaogain((Device\$), REAL Data(\*), (LONG Arraysize))

Input:

Device\$                    The device to be used for this function  
 Arraysize                 Size of the Data array in samples

Return value:    0 = Function was successful  
                               **Data(\*)**: Output gain settings supported by this device  
                               <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the maximum digital input rate of the device

**HTBasic Definition:**

LONG Result=FNHx\_getdevdimart((Device\$), REAL Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful  
                               **Data**: Indicates the maximum digital input rate of the device  
                               <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the triggers supported by this device for digital input tasks

**HTBasic Definition:**

LONG Result=FNHx\_getdevditrus((Device\$), LONG Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates the triggers supported by this device for digital input tasks

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the maximum digital output rate of the device

**HTBasic Definition:**

LONG Result=FNHx\_getdevdomart((Device\$), REAL Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates the maximum digital output rate of the device

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the triggers supported by this device for digital output tasks

**HTBasic Definition:**

LONG Result=FNHx\_getdevdotrus((Device\$), LONG Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates the triggers supported by this device for digital output tasks

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the triggers supported by this device for counter input tasks

**HTBasic Definition:**

LONG Result=FNHx\_getdevcitrus((Device\$), LONG Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates the triggers supported by this device for counter input tasks

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates if the device supports the sample clock timing type for counter input tasks

**HTBasic Definition:**

LONG Result=FNHx\_getdevciscsp((Device\$), LONG Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates if the device supports the sample clock timing for counter input tasks

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates in bits the size of the counters on the device

**HTBasic Definition:**

LONG Result=FNHx\_getdevcimasz((Device\$), LONG Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates in bits the size of the counters on the device

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Indicates in Hertz the maximum counter timebase frequency

**HTBasic Definition:**

LONG Result=FNHx\_getdevcimatb((Device\$), REAL Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates in Hertz the maximum counter timebase frequency

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the triggers supported by this device for counter output tasks

**HTBasic Definition:**

LONG Result=FNHx\_getdevcotrus((Device\$), LONG Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates the triggers supported by this device for counter output tasks

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the size of the counters on the device

**HTBasic Definition:**

`LONG Result=FNHx_getdevcomasz((Device$), LONG Data)`

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates the size of the counters on the device  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Indicates in Hertz the maximum counter timebase frequency

**HTBasic Definition:**

`LONG Result=FNHx_getdevcomatb((Device$), REAL Data)`

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates in Hertz the maximum counter timebase frequency  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the number of DMA channels on the device

**HTBasic Definition:**

`LONG Result=FNHx_getdevnmdchn((Device$), LONG Data)`

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates the number of DMA channels on the device  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the name of the CompactDAQ chassis that contains this module

**HTBasic Definition:**

`LONG Result=FNHx_getdevcdcdnm((Device$), Data$, (LONG BufferSize))`

Input:

Device\$                    The device to be used for this function  
BufferSize                Size of Data\$ string buffer

Return value:    0 = Function was successful

**Data\$:** Name of the CompactDAQ chassis that contains this module  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the slot number in which this module is located in the CompactDAQ chassis

**HTBasic Definition:**

LONG Result=FNHx\_getdevcdslnm((Device\$), LONG Data)

Input:

Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates the slot number in which this module is located in the CompactDAQ chassis

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates if the device(s) detected an overcurrent condition for any virtual channel in the task

**HTBasic Definition:**

LONG Result=FNHx\_getrdocuchex((LONG Taskid), LONG Data)

Input:

Taskid                    The task to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates if the device(s) detected an overcurrent condition for any virtual channel in the task

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Indicates the names of any virtual channels in the task for which an overcurrent condition has been detected

**HTBasic Definition:**

LONG Result=FNHx\_getrdovcuchn((LONG Taskid), Data\$, (LONG BufferSize))

Input:

Taskid                    The task to be used for this function  
BufferSize                Size of Data\$ string buffer

Return value:    0 = Function was successful

**Data\$:** Returns the names of any virtual channels in the task for which an overcurrent condition has been detected

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Indicates if the device(s) detected an open current loop for any virtual channel in the task

**HTBasic Definition:**

LONG Result=FNHx\_getrdoculcex((LONG Taskid), LONG Data)

Input:

Taskid                    The task to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates if the device(s) detected an open current loop for any virtual channel in the task

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Indicates the names of any virtual channels in the task for which the devices detected an open current loop

**HTBasic Definition:**

`LONG Result=FNHx_getrdopculch((LONG Taskid), Data$, (LONG Buffersize))`

Input:

Taskid                    The task to be used for this function  
Buffersize                Size of Data\$ string buffer

Return value:    0 = Function was successful  
                              **Data\$:** Returns the names of any virtual channels in the task  
  for which the devices detected an open current loop  
                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Indicates the number of devices in the task

**HTBasic Definition:**

`LONG Result=FNHx_gettsknumdev((LONG Taskid), LONG Data)`

Input:

Taskid                    The task to be used for this function

Return value:    0 = Function was successful  
                              **Data:** Indicates the number of devices in the task  
                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the action to take when the onboard memory of the device becomes empty

**HTBasic Definition:**

`LONG Result=FNHx_getsmclufbeh((LONG Taskid), LONG Data)`

Input:

Taskid                    The task to be used for this function

Return value:    0 = Function was successful  
                              **Data:** Determines the action to take when the onboard memory  
  of the device becomes empty  
                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the action to take when the onboard memory of the device becomes empty

**HTBasic Definition:**

`LONG Result=FNHx_setsmclufbeh((LONG Taskid), (LONG Data))`

Input:

Taskid                    The task to be used for this function

Return value:    0 = Function was successful  
                              **Data:** Specifies the action to take when the onboard memory  
  of the device becomes empty  
                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the action to take when the onboard memory of the device becomes empty to it's default state

**HTBasic Definition:**

`LONG Result=FNHx_rstsmclufbeh((LONG Taskid))`

Input:

Taskid                    The task to be used for this function

Return value:    0 = Function was successful  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the rate at which to clock the analog-to-digital converter

**HTBasic Definition:**

`LONG Result=FNHx_getaicvratex((LONG Taskid), (Device$), REAL Data)`

Input:

Taskid                    The task to be used for this function  
 Device\$                   The device to be used for this function

Return value:    0 = Function was successful  
                   **Data:** Determines the rate at which to clock the analog-to-digital  
    converter  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the rate at which to clock the analog-to-digital converter

**HTBasic Definition:**

`LONG Result=FNHx_setaicvratex((LONG Taskid), (Device$), (REAL Data))`

Input:

Taskid                    The task to be used for this function  
 Device\$                   The device to be used for this function  
 Data                       Specifies the rate at which to clock the analog-to-digital converter

Return value:    0 = Function was successful  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the rate at which to clock the analog-to-digital converter to it's default state

**HTBasic Definition:**

`LONG Result=FNHx_rstaicvratex((LONG Taskid), (Device$))`

Input:

Taskid                    The task to be used for this function  
 Device\$                   The device to be used for this function

Return value:    0 = Function was successful  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.



**Purpose:** Indicates the maximum convert rate supported by the task, given the current devices and channel count

**HTBasic Definition:**

`LONG Result=FNIHx_getaicvmtx((LONG Taskid), (Device$), REAL Data)`

Input:

Taskid                    The task to be used for this function  
Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates the maximum convert rate supported by the task, given the current devices and channel count

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the terminal of the signal to use as the AI Convert Clock

**HTBasic Definition:**

`LONG Result=FNIHx_getaicvsrctx((LONG Taskid), (Device$), Data$, (LONG Buffersize))`

Input:

Taskid                    The task to be used for this function  
Device\$                    The device to be used for this function  
Buffersize                Size of Data\$ string buffer

Return value:    0 = Function was successful

**Data\$:** The terminal of the signal used as the AI Convert Clock

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the terminal of the signal to use as the AI Convert Clock

**HTBasic Definition:**

`LONG Result=FNIHx_setaicvsrctx((LONG Taskid), (Device$), (Data$))`

Input:

Taskid                    The task to be used for this function  
Device\$                    The device to be used for this function  
Data\$                      The terminal of the signal used as the AI Convert Clock

Return value:    0 = Function was successful

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the terminal of the signal to use as the AI Convert Clock to it's default state

**HTBasic Definition:**

`LONG Result=FNIHx_rstaicvsrctx((LONG Taskid), (Device$))`

Input:

Taskid                    The task to be used for this function  
Device\$                    The device to be used for this function

Return value:    0 = Function was successful

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines on which edge of the clock pulse an analog-to-digital conversion takes place

**HTBasic Definition:**

`LONG Result=FNHx_getaicvaceex((LONG Taskid), (Device$), LONG Data)`

Input:

Taskid                    The task to be used for this function  
Device\$                    The device to be used for this function

Return value:    0 = Function was successful  
                              **Data:** The terminal of the signal used as the AI Convert Clock  
                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies on which edge of the clock pulse an analog-to-digital conversion takes place

**HTBasic Definition:**

`LONG Result=FNHx_setaicvaceex((LONG Taskid), (Device$), (LONG Data))`

Input:

Taskid                    The task to be used for this function  
Device\$                    The device to be used for this function  
Data                        Specifies on which edge of the clock pulse an analog-to-digital conversion takes place

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets on which edge of the clock pulse an analog-to-digital conversion takes place to it's default state

**HTBasic Definition:**

`LONG Result=FNHx_rstaicvaceex((LONG Taskid), (Device$))`

Input:

Taskid                    The task to be used for this function  
Device\$                    The device to be used for this function

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the number of AI Convert Clock Timebase pulses needed to produce a single AI Convert Clock pulse

**HTBasic Definition:**

`LONG Result=FNHx_getaicvtbdex((LONG Taskid), (Device$), LONG Data)`

Input:

Taskid                    The task to be used for this function  
Device\$                    The device to be used for this function

Return value:    0 = Function was successful  
                              **Data:** Number of AI Convert Clock Timebase pulses needed to produce a single AI Convert Clock pulse  
                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the number of AI Convert Clock Timebase pulses needed to produce a single AI Convert Clock pulse

**HTBasic Definition:**

`LONG Result=FNHx_setaicvtbdex((LONG Taskid), (Device$), (LONG Data))`

Input:

Taskid	The task to be used for this function
Device\$	The device to be used for this function
Data	Specifies the number of AI Convert Clock Timebase pulses needed to produce a single AI Convert Clock pulse

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets on which edge of the clock pulse an analog-to-digital conversion takes place to it's default state

**HTBasic Definition:**

`LONG Result=FNHx_rstaicvtbdex((LONG Taskid), (Device$))`

Input:

Taskid	The task to be used for this function
Device\$	The device to be used for this function

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the terminal of the signal to use as the AI Convert Clock Timebase

**HTBasic Definition:**

`LONG Result=FNHx_getaicvtbsex((LONG Taskid), (Device$), LONG Data)`

Input:

Taskid	The task to be used for this function
Device\$	The device to be used for this function

Return value: 0 = Function was successful  
                   **Data:** Determines the terminal of the signal to use as the AI Convert Clock Timebase  
 <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the terminal of the signal to use as the AI Convert Clock Timebase

**HTBasic Definition:**

`LONG Result=FNHx_setaicvtbsex((LONG Taskid), (Device$), (LONG Data))`

Input:

Taskid	The task to be used for this function
Device\$	The device to be used for this function
Data	Specifies the terminal of the signal to use as the AI Convert Clock Timebase

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the terminal of the signal to use as the AI Convert Clock Timebase to it's default state

**HTBasic Definition:**

`LONG Result=FNHx_rstaicvtbsex((LONG Taskid), (Device$))`

Input:

Taskid	The task to be used for this function
Device\$	The device to be used for this function

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the units of delay

**HTBasic Definition:**

`LONG Result=FNHx_getdlfscdunt((LONG Taskid), LONG Data)`

Input:

Taskid	The task to be used for this function
--------	---------------------------------------

Return value: 0 = Function was successful  
                  **Data:** Determines the units of delay  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the units of delay

**HTBasic Definition:**

`LONG Result=FNHx_setdlfscdunt((LONG Taskid), (LONG Data))`

Input:

Taskid	The task to be used for this function
Data	Specifies the units of delay

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the units of delay to it's default state

**HTBasic Definition:**

`LONG Result=FNIHx_rstdlfsdunt((LONG Taskid))`

**Input:**

Taskid                    The task to be used for this function

**Return value:**    0 = Function was successful  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the units of delay

**HTBasic Definition:**

`LONG Result=FNIHx_getdlfsdoux((LONG Taskid), (Device$), LONG Data)`

**Input:**

Taskid                    The task to be used for this function  
 Device\$                   The device to be used for this function

**Return value:**    0 = Function was successful  
                                      **Data:** Determines the units of delay  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the units of delay

**HTBasic Definition:**

`LONG Result=FNIHx_setdlfsdoux((LONG Taskid), (Device$), (LONG Data))`

**Input:**

Taskid                    The task to be used for this function  
 Device\$                   The device to be used for this function  
 Data                       Specifies the units of delay

**Return value:**    0 = Function was successful  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the units of delay to it's default state

**HTBasic Definition:**

`LONG Result=FNIHx_rstdlfsdoux((LONG Taskid), (Device$))`

**Input:**

Taskid                    The task to be used for this function  
 Device\$                   The device to be used for this function

**Return value:**    0 = Function was successful  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Get the amount of time to wait after receiving a Sample Clock edge before beginning to acquire the sample

**HTBasic Definition:**

`LONG Result=FNHx_getdlfscdlay((LONG Taskid), REAL Data)`

Input:

Taskid                    The task to be used for this function

Return value:    0 = Function was successful

**Data:** Amount of time to wait after receiving a Sample Clock edge before beginning to acquire the sample  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Set the amount of time to wait after receiving a Sample Clock edge before beginning to acquire the sample

**HTBasic Definition:**

`LONG Result=FNHx_setdlfscdlay((LONG Taskid), (REAL Data))`

Input:

Taskid                    The task to be used for this function  
Data                        Amount of time to wait after receiving a Sample Clock edge before beginning to acquire the sample

Return value:    0 = Function was successful

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the units of delay to it's default state

**HTBasic Definition:**

`LONG Result=FNHx_rstdlfscdlay((LONG Taskid))`

Input:

Taskid                    The task to be used for this function

Return value:    0 = Function was successful

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Gets the amount of time to wait after receiving a Sample Clock edge before beginning to acquire the sample

**HTBasic Definition:**

`LONG Result=FNHx_getdlfscdlex((LONG Taskid), (Device$), REAL Data)`

Input:

Taskid                    The task to be used for this function  
Device\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Determines the amount of time to wait after receiving a Sample Clock edge before beginning to aquire the sample  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Sets the amount of time to wait after receiving a Sample Clock edge before beginning to acquire the sample

**HTBasic Definition:**

`LONG Result=FNHx_setdlfscdlex((LONG Taskid), (Device$), (REAL Data))`

Input:

Taskid	The task to be used for this function
Device\$	The device to be used for this function
Data	Specifies the amount of time to wait after receiving a Sample Clock edge before beginning to acquire the sample

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the amount of time to wait after receiving a Sample Clock edge before beginning to acquire the sample to it's default state

**HTBasic Definition:**

`LONG Result=FNHx_rstdlfscdlex((LONG Taskid), (Device$))`

Input:

Taskid	The task to be used for this function
Device\$	The device to be used for this function

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines the physical channels to use for pattern matching

**HTBasic Definition:**

`LONG Result=FNHx_getdigpptsr((LONG Taskid), Data$, (LONG Buffersize))`

Input:

Taskid	The task to be used for this function
Buffersize	Size of Data\$ string buffer

Return value: 0 = Function was successful  
                  Data\$: The physical channels to use for pattern matching  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the physical channels to use for pattern matching

**HTBasic Definition:**

`LONG Result=FNHx_setdigpptsr((LONG Taskid), (Data$))`

Input:

Taskid	The task to be used for this function
Data\$	The physical channels to use for pattern matching

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.





**Purpose:** Determines if the Pause Trigger occurs when the physical channels specified with Source match or differ from the digital pattern specified with Pattern

**HTBasic Definition:**

LONG Result=FNHx\_getdigpptwhn((LONG Taskid), LONG Data)

Input:

Taskid                    The task to be used for this function

Return value:    0 = Function was successful

**Data:** Determines if the Pause Trigger occurs when the physical channels specified with Source match or differ from the digital pattern specified with Pattern to occur

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies if the Pause Trigger occurs when the physical channels specified with Source match or differ from the digital pattern specified with Pattern

**HTBasic Definition:**

LONG Result=FNHx\_setdigpptwhn((LONG Taskid), (LONG Data))

Input:

Taskid                    The task to be used for this function  
Data                        Specifies if the Pause Trigger occurs when the physical channels specified with Source match or differ from the digital pattern specified with Pattern

Return value:    0 = Function was successful

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets if the Pause Trigger occurs when the physical channels specified with Source match or differ from the digital pattern specified with Pattern to occur to it's default state

**HTBasic Definition:**

LONG Result=FNHx\_rstdigpptwhn((LONG Taskid))

Input:

Taskid                    The task to be used for this function

Return value:    0 = Function was successful

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Determines if the device(s) detected an overcurrent condition for any channel in the task

**HTBasic Definition:**

LONG Result=FNHx\_getwrovcchex((LONG Taskid), LONG Data)

Input:

Taskid                      The task to be used for this function

Return value:    0 = Function was successful

**Data:** Determines if the device(s) detected an overcurrent condition for any channel in the task

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Determines the names of any virtual channels in the task for which an overcurrent condition has been detected

**HTBasic Definition:**

LONG Result=FNHx\_getwrovccchns((LONG Taskid), Data\$, (LONG Buffersize))

Input:

Taskid                      The task to be used for this function  
Buffersize                  Size of Data\$ string buffer

Return value:    0 = Function was successful

**Data\$:** The names of any virtual channels in the task for which an overcurrent condition has been detected

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Determines if the device(s) detected an open current loop for any channel in the task

**HTBasic Definition:**

LONG Result=FNHx\_getwropcchex((LONG Taskid), LONG Data)

Input:

Taskid                      The task to be used for this function

Return value:    0 = Function was successful

**Data:** Determines if the device(s) detected an open current loop for any channel in the task

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Determines the names of any virtual channels in the task for which an open current loop has been detected

**HTBasic Definition:**

LONG Result=FNHx\_getwropccchns((LONG Taskid), Data\$, (LONG Buffersize))

Input:

Taskid                      The task to be used for this function  
Buffersize                  Size of Data\$ string buffer

Return value:    0 = Function was successful

**Data\$:** The names of any virtual channels in the task for which an open current loop has been detected

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Determines if the device(s) detected a power supply fault for any channel in the task

**HTBasic Definition:**

LONG Result=FNHx\_getwrpsfchex((LONG Taskid), LONG Data)

Input:

Taskid                      The task to be used for this function

Return value:    0 = Function was successful

**Data:** Determines if the device(s) detected a power supply fault for any channel in the task

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Determines the names of any virtual channels in the task that have a power supply fault

**HTBasic Definition:**

LONG Result=FNHx\_getwrpsfchns((LONG Taskid), Data\$, (LONG BufferSize))

Input:

Taskid                      The task to be used for this function  
BufferSize                  Size of Data\$ string buffer

Return value:    0 = Function was successful

**Data\$:** The names of any virtual channels in the task that have a power supply fault

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Determines that the next samples written are the last samples you want to generate

**HTBasic Definition:**

LONG Result=FNHx\_getwrnxwri1a((LONG Taskid), LONG Data)

Input:

Taskid                      The task to be used for this function

Return value:    0 = Function was successful

**Data:** Determines that the next samples written are the last samples you want to generate

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies that the next samples written are the last samples you want to generate

**HTBasic Definition:**

LONG Result=FNHx\_setwrnxwri1a((LONG Taskid), (LONG Data))

Input:

Taskid                      The task to be used for this function

Return value:    0 = Function was successful

**Data:** Specifies that the next samples written are the last samples you want to generate

<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the behaviour of this function to it's default state

**HTBasic Definition:**

LONG Result=FNIx\_setwrnxwrla((LONG Taskid))

Input:

Taskid                    The task to be used for this function

Return value:    0 = Function was successful  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the list of terminal configurations supported by the channel

**HTBasic Definition:**

LONG Result=FNIx\_getphchaitcf((Physchan\$), LONG Data)

Input:

Physchan\$                The device to be used for this function

Return value:    0 = Function was successful  
                                  **Data:** List of terminal configurations by the channel  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the list of terminal configurations supported by the channel

**HTBasic Definition:**

LONG Result=FNIx\_getphchaotcf((Physchan\$), LONG Data)

Input:

Physchan\$                The device to be used for this function

Return value:    0 = Function was successful  
                                  **Data:** List of terminal configurations by the channel  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates in bits the width of digital input port

**HTBasic Definition:**

LONG Result=FNIx\_getphchdipwd((Physchan\$), LONG Data)

Input:

Physchan\$                The device to be used for this function

Return value:    0 = Function was successful  
                                  **Data:** The bits width of digital input port  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

<b>FNHx_getphchdiscs</b>	<b>NI-DAQmx 8.3: DAQmxGetPhysicalChanDISampClkSupported</b>
--------------------------	---

**Purpose:** Indicates if the sample clock timing type is supported for the digital input physical channel

**HTBasic Definition:**

LONG Result=FNHx\_getphchdiscs((Physchan\$), LONG Data)

Input:

Physchan\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates if the sample clock timing type is supported for the digital input physical channel

                  <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

<b>FNHx_getphchdidsp</b>	<b>NI-DAQmx 8.3: DAQmxGetPhysicalChanDIChangeDetectSupported</b>
--------------------------	--

**Purpose:** Indicates if the change detection timing type is supported for the digital input physical channel

**HTBasic Definition:**

LONG Result=FNHx\_getphchdidsp((Physchan\$), LONG Data)

Input:

Physchan\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates if the change detection timing type is supported for the digital input physical channel

                  <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

<b>FNHx_getphchdopwd</b>	<b>NI-DAQmx 8.3: DAQmxGetPhysicalChanDOPortWidth</b>
--------------------------	--

**Purpose:** Indicates in bits the width of digital output port

**HTBasic Definition:**

LONG Result=FNHx\_getphchdopwd((Physchan\$), LONG Data)

Input:

Physchan\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** The bits width of digital output port

                  <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

<b>FNHx_getphchdoscs</b>	<b>NI-DAQmx 8.3: DAQmxGetPhysicalChanDOSampClkSupported</b>
--------------------------	---

**Purpose:** Indicates if the sample clock timing type is supported for the digital output physical channel

**HTBasic Definition:**

LONG Result=FNHx\_getphchdoscs((Physchan\$), LONG Data)

Input:

Physchan\$                    The device to be used for this function

Return value:    0 = Function was successful

**Data:** Indicates if the sample clock timing type is supported for the digital output physical channel

                  <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Indicates the terminal of the signal to use as the Sample Clock Timebase

**HTBasic Definition:**

LONG Result=FNHx\_getsctresmod((LONG Taskid), LONG Data)

Input:

Taskid                    The task to be used for this function

Return value:    0 = Function was successful

**Data:** Terminal of the signal to use as the Sample Clock Timebase  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the terminal of the signal to use as the Sample Clock Timebase

**HTBasic Definition:**

LONG Result=FNHx\_getsctresmod((LONG Taskid), (LONG Data))

Input:

Taskid                    The task to be used for this function  
Data                        Terminal of the signal to use as the Sample Clock Timebase

Return value:    0 = Function was successful

                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Resets the terminal of the signal to use as the Sample Clock Timebase to it's default state

**HTBasic Definition:**

LONG Result=FNHx\_rstsctresmod((LONG Taskid))

Input:

Taskid                    The task to be used for this function

Return value:    0 = Function was successful

                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.3 and higher.

**Purpose:** Specifies the terminal of the signal to use as the Sample Clock Timebase

**HTBasic Definition:**

LONG Result=FNHx\_getsctresmod((LONG Taskid), (LONG Data))

Input:

Taskid                    The task to be used for this function  
Data                        Terminal of the signal to use as the Sample Clock Timebase

Return value:    0 = Function was successful

                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

**Purpose:** Resets the terminal of the signal to use as the Sample Clock Timebase to it's default state

**HTBasic Definition:**

```
LONG Result=FNHx_rstscrtresmod((LONG Taskid))
```

Input:

Taskid                    The task to be used for this function

Return value:    0 = Function was successful  
                  <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.5 and higher.

<b>FNHx_getaivoacrmu</b>	<b>NI-DAQmx 8.6: DAQmxGetAIVoltageACRMSUnits</b>
--------------------------	--

**Purpose:** Determines the units to use to return voltage RMS measurements from the channel

**HTBasic Definition:**

```
LONG Result=FNHx_getaivoacrmu((LONG Taskid), (Channel$), LONG Data)
```

Input:

Taskid                    The task used in this function  
Channel\$                  The name of the channel you want to work with

Return value:    0 = Function was successful  
                  **Data**        Determines the units for the channel  
                  <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

<b>FNHx_setaivoacrmu</b>	<b>NI-DAQmx 8.6: DAQmxSetAIVoltageACRMSUnits</b>
--------------------------	--

**Purpose:** Specifies the units to use to return voltage RMS measurements from the channel

**HTBasic Definition:**

```
LONG Result=FNHx_setaivoacrmu((LONG Taskid), (Channel$), (LONG Data))
```

Input:

Taskid                    The task used in this function  
Channel\$                  The name of the channel you want to set  
Data                      Specifies the units for the channel

Return value:    0 = Function was successful  
                  <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

<b>FNHx_rstaivoacrmu</b>	<b>NI-DAQmx 8.6: DAQmxResetAIVoltageACRMSUnits</b>
--------------------------	--

**Purpose:** Resets the units to use to return voltage RMS measurements from the channel to default values

**HTBasic Definition:**

```
LONG Result=FNHx_rstaivoacrmu((LONG Taskid), (Channel$))
```

Input:

Taskid                    The task used in this function  
Channel\$                  The name of the channel you want to reset

Return value:    0 = Function was successful  
                  <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Determines the units to use to return current RMS measurements from the channel

**HTBasic Definition:**

`LONG Result=FNHx_getaicuacrmu((LONG Taskid), (Channel$), LONG Data)`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to work with

Return value: 0 = Function was successful  
                   Data Determines the units for the channel  
 <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Specifies the units to use to return current RMS measurements from the channel

**HTBasic Definition:**

`LONG Result=FNHx_setaicuacrmu((LONG Taskid), (Channel$), (LONG Data))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to set
Data	Specifies the units for the channel

Return value: 0 = Function was successful  
 <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Resets the units to use to return current RMS measurements from the channel to default values

**HTBasic Definition:**

`LONG Result=FNHx_rstaicuacrmu((LONG Taskid), (Channel$))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to reset

Return value: 0 = Function was successful  
 <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.



**Purpose:** Determines the current limit, in amperes, for the voltage channel

**HTBasic Definition:**

`LONG Result=FNHx_getaoculimit((LONG Taskid), (Channel$), REAL Data)`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to work with

Return value: 0 = Function was successful  
                   **Data** Determines the current limit in Ampere for the channel  
 <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Specifies the current limit, in amperes, for the voltage channel

**HTBasic Definition:**

`LONG Result=FNHx_setaoculimit((LONG Taskid), (Channel$), (REAL Data))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to set
Data	Specifies the current limit in Ampere for the channel

Return value: 0 = Function was successful  
 <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Resets the current limit, in amperes, for the voltage channel, to it's default value

**HTBasic Definition:**

`LONG Result=FNHx_rstaoculimit((LONG Taskid), (Channel$))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to reset

Return value: 0 = Function was successful  
 <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Determines the kind of the waveform to generate

**HTBasic Definition:**

`LONG Result=FNHx_getaofungtyp((LONG Taskid), (Channel$), LONG Data)`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to work with

Return value: 0 = Function was successful  
                   **Data**     Determines the kind of waveform to generate  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Specifies the kind of the waveform to generate

**HTBasic Definition:**

`LONG Result=FNHx_setaofungtyp((LONG Taskid), (Channel$), (LONG Data))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to set
Data	Specifies the kind of waveform to generate

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Resets the kind of the waveform to generate to it's default state

**HTBasic Definition:**

`LONG Result=FNHx_rstaofungtyp((LONG Taskid), (Channel$))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to reset

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Determines the frequency of the waveform to generate in hertz

**HTBasic Definition:**

`LONG Result=FNHx_getaofungfrq((LONG Taskid), (Channel$), REAL Data)`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to work with

Return value: 0 = Function was successful  
Data Determines the frequency of the waveform to generate  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Specifies the frequency of the waveform to generate in hertz

**HTBasic Definition:**

`LONG Result=FNHx_setaofungfrq((LONG Taskid), (Channel$), (REAL Data))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to set
Data	Specifies the frequency of the waveform to generate

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Resets the frequency of the waveform to generate to it's default value

**HTBasic Definition:**

`LONG Result=FNHx_rstaofungfrq((LONG Taskid), (Channel$))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to reset

Return value: 0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Determines the zero-to-peak amplitude of the waveform to generate in volts

**HTBasic Definition:**

`LONG Result=FNHx_getaofungamp((LONG Taskid), (Channel$), REAL Data)`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to work with

Return value: 0 = Function was successful  
                   **Data**     Determines the zero-to-peak amplitude of the waveform to generate in Volt  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Specifies the zero-to-peak amplitude of the waveform to generate in volts

**HTBasic Definition:**

`LONG Result=FNHx_setaofungamp((LONG Taskid), (Channel$), (REAL Data))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to set
Data	Specifies the zero-to-peak amplitude of the waveform to generate in Volt

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Resets the zero-to-peak amplitude of the waveform to generate to it's default value

**HTBasic Definition:**

`LONG Result=FNHx_rstaofungamp((LONG Taskid), (Channel$))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to reset

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Determines the voltage offset of the waveform to generate

**HTBasic Definition:**

`LONG Result=FNHx_getaofungoff((LONG Taskid), (Channel$), REAL Data)`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to work with

Return value: 0 = Function was successful  
                   Data Determines the voltage offset of the waveform to generate  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Specifies the voltage offset of the waveform to generate

**HTBasic Definition:**

`LONG Result=FNHx_setaofungoff((LONG Taskid), (Channel$), (REAL Data))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to set
Data	Specifies the voltage offset of the waveform to generate

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Resets the voltage offset of the waveform to generate to it's default value

**HTBasic Definition:**

`LONG Result=FNHx_rstaofungoff((LONG Taskid), (Channel$))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to reset

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Determines the square wave duty cycle of the waveform to generate

**HTBasic Definition:**

`LONG Result=FNHx_getaofungsd((LONG Taskid), (Channel$), REAL Data)`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to work with

Return value: 0 = Function was successful  
                   **Data**     Determines the square wave duty cycle of the waveform to generate  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Specifies the square wave duty cycle of the waveform to generate

**HTBasic Definition:**

`LONG Result=FNHx_setaofungsd((LONG Taskid), (Channel$), (REAL Data))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to set
Data	Specifies the square wave duty cycle of the waveform to generate

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Resets the square wave duty cycle of the waveform to generate to it's default value

**HTBasic Definition:**

`LONG Result=FNHx_rstaofungsd((LONG Taskid), (Channel$))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to reset

Return value: 0 = Function was successful  
                   <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Determines if the device generates a modulated version of the waveform using the original waveform as a carrier and input from an external terminal as the signal

**HTBasic Definition:**

`LONG Result=FNHx_getaofungmtp((LONG Taskid), (Channel$), LONG Data)`

Input:

Taskid                    The task used in this function  
Channel\$                  The name of the channel you want to work with

Return value:    0 = Function was successful  
                              **Data** Determines the modulated version of the waveform (FM/AM/no)  
                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Specifies if the device generates a modulated version of the waveform using the original waveform as a carrier and input from an external terminal as the signal

**HTBasic Definition:**

`LONG Result=FNHx_setaofungmtp((LONG Taskid), (Channel$), (LONG Data))`

Input:

Taskid                    The task used in this function  
Channel\$                  The name of the channel you want to set  
Data                        Specifies the modulated version of the waveform (FM/AM/no)

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Resets if the device generates a modulated version of the waveform using the original waveform as a carrier and input from an external terminal as the signal to it's default state

**HTBasic Definition:**

`LONG Result=FNHx_rstaofungmtp((LONG Taskid), (Channel$))`

Input:

Taskid                    The task used in this function  
Channel\$                  The name of the channel you want to reset

Return value:    0 = Function was successful  
                              <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Determines the FM deviation in hertz per volt when Type is DAQmx\_Val\_FM

**HTBasic Definition:**

`LONG Result=FNHx_getaofungfmd((LONG Taskid), (Channel$), REAL Data)`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to work with

Return value: 0 = Function was successful  
                   Data Determines the FM deviation in Hertz per Volt  
 <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Specifies the FM deviation in hertz per volt when Type is DAQmx\_Val\_FM

**HTBasic Definition:**

`LONG Result=FNHx_setaofungfmd((LONG Taskid), (Channel$), (REAL Data))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to set
Data	Specifies the FM deviation in Hertz per Volt

Return value: 0 = Function was successful  
 <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Resets the FM deviation in hertz per volt when Type is DAQmx\_Val\_FM to it's default value

**HTBasic Definition:**

`LONG Result=FNHx_rstaofungfmd((LONG Taskid), (Channel$))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to reset

Return value: 0 = Function was successful  
 <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.



**Purpose:** Determines which timing engine to use for the specified timing type. Refer to device documentation for information on supported timing engines

**HTBasic Definition:**

`LONG Result=FNHx_getsamtimeng((LONG Taskid), (Channel$), LONG Data)`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to work with

Return value: 0 = Function was successful  
                   Data Determines which timing engine to use  
 <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Specifies which timing engine to use for the specified timing type. Refer to device documentation for information on supported timing engines

**HTBasic Definition:**

`LONG Result=FNHx_setsamtimeng((LONG Taskid), (Channel$), (LONG Data))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to set
Data	Specifies which timing engine to use

Return value: 0 = Function was successful  
 <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Resets which timing engine to use for the specified timing type to it's default state

**HTBasic Definition:**

`LONG Result=FNHx_rstsamtimeng((LONG Taskid), (Channel$))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to reset

Return value: 0 = Function was successful  
 <>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Determines if you can control the physical channel externally via a manual control located on the device

**HTBasic Definition:**

`LONG Result=FNHx_getphchaomce((Channel$), LONG Data)`

Input:

Channel\$                    The name of the channel you want to work with

Return value:

0 = Function was successful  
    **Data** Determines if you can control the physical channel externally via a manual control located on the device  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Specifies if you can control the physical channel externally via a manual control located on the device

**HTBasic Definition:**

`LONG Result=FNHx_setphchaomce((Channel$), (LONG Data))`

Input:

Channel\$                    The name of the channel you want to set  
Data                         Determines if you can control the physical channel externally via a manual control located on the device

Return value:

0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Resets if you can control the physical channel externally via a manual control located on the device to it's default state

**HTBasic Definition:**

`LONG Result=FNHx_rstphchaomce((Channel$))`

Input:

Channel\$                    The name of the channel you want to reset

Return value:

0 = Function was successful  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Indicates the current value of the front panel amplitude control for the physical channel in volts

**HTBasic Definition:**

**LONG** Result=FNHx\_getpcaomcamp((Channel\$), **REAL** Data)

Input:

Channel\$                    The name of the channel you want to work with

Return value:    0 = Function was successful

**Data**    Indicates the current value of the front panel amplitude control for the physical channel  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Indicates the current value of the front panel frequency control for the physical channel in hertz

**HTBasic Definition:**

**LONG** Result=FNHx\_getpcaomcfreq((Channel\$), **REAL** Data)

Input:

Channel\$                    The name of the channel you want to work with

Return value:    0 = Function was successful

**Data**    Indicates the current value of the front panel frequency control for the physical channel in Hertz  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.6 and higher.

**Purpose:** Indicates a list of all terminals on the device

**HTBasic Definition:**

**LONG** Result=FNHx\_getdevtermin((Device\$), Data\$, (**LONG** Bufsize))

Input:

Device\$                    The name of the device you want to work with  
Bufsize                    The dimensioned length of the string variable Data\$

Return value:    0 = Function was successful

**Data\$**    A list of all terminals on the device  
<>0 = Error occurred

**Note:** This function is only available if you are using NI-DAQmx 8.7 and higher.

**Purpose:** Set buffer attributes.

**HTBasic Definition:**

`LONG Result= FNHx_setbufattrib((LONG Taskid), (LONG Attrib1), (LONG Attrib2))`

Input:

Taskid	The task used in this function
Attrib1	First attribute
Attrib2	Second attribute

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Set trigger attributes.

**HTBasic Definition:**

`LONG Result= FNHx_settriattrib((LONG Taskid), (LONG Attrib1), (LONG Attrib2))`

Input:

Taskid	The task used in this function
Attrib1	First attribute
Attrib2	Second attribute

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Set counter output pulse frequency.

**HTBasic Definition:**

`LONG Result= FNHx_setcopulfreq((LONG Taskid), (Channel$), (REAL Data))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to work with
Data	Pulse frequency

Return value: 0 = Function was successful  
<>0 = Error occurred

**Purpose:** Set counter output duty cycle.

**HTBasic Definition:**

`LONG Result= FNHx_setcopdutcyc((LONG Taskid), (Channel$), (REAL Data))`

Input:

Taskid	The task used in this function
Channel\$	The name of the channel you want to work with
Data	Duty cycle

Return value: 0 = Function was successful  
<>0 = Error occurred

## Callback services

### Introduction: Software Events

Software events provide an asynchronous notification mechanism for a set of DAQ events. Unlike hardware events, software events do not require you to use a thread to wait until data is available. Using event-based programming, you can write an application that continues to perform work while waiting for data without resorting to developing a multi-threaded application.

NI-DAQmx supports the following software events:

**Every N Samples Acquired Into Buffer Event**—Occurs when the user-defined number of samples is written from the device to the PC buffer. This event works only with devices that support buffered tasks. *Note:* The value you set for this event must be evenly divisible into the buffer size if you are using DMA as your data transfer mechanism. For instance, if the buffer size is 1,000 samples, specifying 102 for this software event generates an error. Specifying 100, however, would not generate an error. If you are using IRQ as the data transfer method, the value does not need to be evenly divisible. With IRQ, however, the Data Transfer Request Condition attribute/property can affect when this software event occurs.

**Every N Samples Transferred From Buffer Event**—Occurs when the user-defined number of samples is written from the PC buffer to the device. This event works only with devices that support buffered tasks. *Note:* The value you set for this event must be evenly divisible into the buffer size if you are using DMA as your data transfer mechanism. For instance, if the buffer size is 1,000 samples, specifying 102 for this software event generates an error. Specifying 100, however, would not generate an error. If you are using IRQ as the data transfer method, the value does not need to be evenly divisible. With IRQ, however, the Data Transfer Request Condition attribute/property can affect when this software event occurs.

**Done Event**—Occurs when the task completes execution or when an error causes the task to finish. Recoverable errors that do not cause the task to finish do not cause this event to fire.

**Signal Event**—Occurs when the specified hardware signal occurs. Supported signals include the counter output event, change detection event, sample complete event, and the sample clock. The DAQmx library supports callback services.

### HTBasic: Software Events creates ON SIGNAL events

All events are supported by the HTBasic library too. For every event handler you create using the functions described later in this chapter you must register an ON SIGNAL event handler in your HTBasic program for every context where the event may occur. In case such an event is generated in the driver HTBasic gets a SIGNAL and calls the defined ON SIGNAL handler where you can e.g. retrieve your measurement data from the buffer or transfer new data into the output buffer or whatever else might be useful.

Several example programs which demonstrates the events can be found in the subdirectory "HTB Examples\Events" inside the installation directory of the HTBasic library.

**FNHx\_regevnsmevnt**

NI-DAQmx: **DAQmxRegisterEveryNSamplesEvent**

**Purpose:** Registers a callback function to receive an event when the specified number of samples is written from the device to the buffer or from the buffer to the device. This function only works with devices that support buffered tasks.

#### HTBasic Definition:

```
LONG Result=FNHx_regevnsmevnt((LONG Taskid), (LONG Eventtype), (LONG Nsamples),  
                              (LONG Options), (LONG Htb_signal_id))
```

#### Input:

Taskid	The task used in this function
Eventtype	The type of event you want to receive
Nsamples	The number of samples after which each event should occur
Options	Use this parameter to set certain options (see DAQmx documentation)
Htb_signal_id	The ID of the HTBasic signal which should be generated (0 .. 15)

**Return value:** 0 = Function was successful  
<>0 = Error occurred

#### Example:

```
INTEGER Iresult, Debug  
LONG Taskid, Nsamples, Options, Eventtype, Htb_signal_id  
Nsamples=1000 ! callback every ... samples  
Options=0 ! no special options  
Htb_signal_id=15 ! signal 15 for callback  
Eventtype=FNHx_get_constant("DAQmx_Val_Acquired_Into_Buffer", Iresult, Debug)  
!  
Result=FNHx_regevnsmevnt((Taskid), (Eventtype), (Nsamples), (Options), (Htb_signal_id))  
!  
ON SIGNAL Htb_signal_id,15 CALL Cb_everynsample ! define ON SIGNAL routine
```

**Purpose:** Registers a callback function to receive an event when the specified hardware event occurs.

**HTBasic Definition:**

```
LONG Result=FNHx_regsignevent((LONG Taskid), (LONG Signal_id), (LONG Options),
                              (LONG Htb_signal_id))
```

Input:

Taskid	The task used in this function
Signal_id	The signal for which you want to receive results. Possible values: DAQmx_Val_SampleClock, DAQmx_Val_SampleCompleteEvent DAQmx_Val_ChangeDetectionEvent, DAQmx_Val_CounterOutputEvent
Options	Use this parameter to set certain options (see DAQmx documentation)
Htb_signal_id	The ID of the HTBasic signal which should be generated (0 .. 15)

Return value: 0 = Function was successful  
<>0 = Error occurred

**Example:**

```
INTEGER Iresult, Debug
LONG Taskid, Options, Signal_id, Htb_signal_id
Signal_id=FNHx_get_constant("DAQmx_Val_SampleCompleteEvent", Iresult, Debug)
Options=0 ! no special options
Htb_signal_id=14 ! signal 14 for callback
!
Result=FNHx_regsignevent((Taskid), (Signal_id), (Options), (Htb_signal_id))
!
ON SIGNAL Htb_signal_id,15 CALL Cb_signalevent ! define ON SIGNAL routine
```

**Purpose:** Registers a callback function to receive an event when a task stops due to an error or when a finite acquisition task or finite generation task completes execution.

**HTBasic Definition:**

```
LONG Result=FNHx_regdoneevent((LONG Taskid), (LONG Options), (LONG Htb_signal_id))
```

Input:

Taskid	The task used in this function
Options	Use this parameter to set certain options (see DAQmx documentation)
Htb_signal_id	The ID of the HTBasic signal which should be generated (0 .. 15)

Return value: 0 = Function was successful  
<>0 = Error occurred

**Example:**

```
INTEGER Iresult, Debug
LONG Taskid, Options, Htb_signal_id
Options=0 ! no special options
Htb_signal_id=13 ! signal 14 for callback
!
Result=FNHx_regdoneevent((Taskid), (Options), (Htb_signal_id))
!
ON SIGNAL Htb_signal_id,15 CALL Cb_Doneevent ! define ON SIGNAL routine
```

## Watchdog functions

**FNHx\_crwdogtimtsk**

NI-DAQmx: **DAQmxCreateWatchdogTimerTask**

**Purpose:** Creates and configures a task that controls the watchdog timer of a device

### HTBasic Definition:

```
LONG Result=FNHx_crwdogtimtsk((Devicename$), (Taskname$), (LONG Taskid), (REAL T_out),  
                               (Lines$), (LONG Expstate))
```

#### Input:

Devicename\$	The name of the device, as configured in MAX, to which this operation applies
Taskname\$	The name of the task
T_out	The amount of time, in seconds, until the watchdog timer expires
Lines\$	The digital line or port to modify. You cannot modify dedicated digital input lines. You can specify a list or range of physical channels
Expstate	The state to which to set the digital physical channel when the watchdog timer expires

Return value: 0 = Function was successful  
                  **Taskid** now contains Task ID of newly created task  
                  <>0 = Error occurred

**Please note:** This function differs from the original DAQmx function since you can only create a watchdog task for one digital line or line range with one function call. In case you need to change more than one line or line range you must create several separate watchdog timers.

**FNHx\_ctrlwdogtsk**

NI-DAQmx: **DAQmxControlWatchdogTask**

**Purpose:** Controls the watchdog task according to the action you specify

### HTBasic Definition:

```
LONG Result=FNHx_ctrlwdogtsk((LONG Taskid), (LONG Action))
```

#### Input:

Taskid	The task used in this function
Action	Specifies how to control the watchdog task

Return value: 0 = Function was successful  
                  <>0 = Error occurred

## Storage

**FNHx\_savetask**

NI-DAQmx: **DAQmxSaveTask**

**Purpose:** Saves the specified task and any local channels it contains to MAX

### HTBasic Definition:

```
LONG Result=FNHx_savetask((LONG Taskid), (Saveas$), (Author$), (LONG Options))
```

#### Input:

Taskid	The task to save
Saveas\$	Name to save the task as. If you pass an empty string ("") or NULL, the name currently assigned to the task will be used
Author\$	Name to store with the task

**Return value:** 0 = Function was successful  
<>0 = Error occurred

### Example:

```
LONG Result, Taskid, Options, Options_1, Options_2
DIM Saveas$[255], Author$[255]
INTEGER Iresult, Debug
Saveas$="MySavedTaskName"
Author$="MyName"
Options_1=FNHx_get_constant("DAQmx_Val_Save_Overwrite",Iresult,Debug)
Options_2=FNHx_get_constant("DAQmx_Val_Save_AllowInteractiveEditing",Iresult,Debug)
Options=Options_1+Options_2
Result=FNHx_savetask((Taskid), (Saveas$), (Author$), (Options))
!
IF Result=0 THEN
    !... subsequent commands...
END IF
```

**FNHx\_delsavedtask**

NI-DAQmx: **DAQmxDeleteSavedTask**

**Purpose:** Deletes the specified task from MAX

### HTBasic Definition:

```
LONG Result=FNHx_delsavedtask((Taskname$))
```

#### Input:

Taskname\$	Name of the saved task to be deleted
------------	--------------------------------------

**Return value:** 0 = Function was successful  
<>0 = Error occurred

### Example:

```
LONG Result
DIM Saveas$[255]
Saveas$="MySavedTaskName"
Result=FNHx_delsavedtask((Saveas$))
!
IF Result=0 THEN
    !... subsequent commands...
END IF
```



## HTBasic specific functions

### **FNHx\_manage\_subs**

**Purpose:** Load / unload and init / de-init HTBasic library functions

**HTBasic Definition:**

```
LONG Result=FNHx_manage_subs(Func$, OPTIONAL INTEGER Verbose, Debug)
```

Input:

```
Func$           Function code  
                "LOAD"    = Load and initialize DAQmx_HTB library  
                "UNLOAD"  = Unload and de-initialize DAQmx_HTB library  
Verbose         Verbose messages (0=none, 1=all, 2=Errors only)  
Debug
```

Return value: 0 = Function was successful  
<>0 = Error occurred

**Example:**

```
INTEGER Verbose  
IF NOT FNHx_manage_subs("LOAD", Verbose) THEN STOP ! load and init DAQmx library
```

### **FNHx\_get\_constant**

**Purpose:** Returns the value of a numeric constant from NIDAQmx.h from the NI-DAQmx Library

**HTBasic Definition:**

```
DEF FNHx_get_constant(Constant$, INTEGER Result, OPTIONAL INTEGER Debug)
```

Input:

```
Constant$       Name of DAQmx constant as defined in file "NIDAQmx.h"  
Debug           Shows debug infos (constant found / not found etc.)  
                Debug = 3 : Show all debug messages
```

Return value: 0 = Function was successful  
<>0 = Error occurred

If the function was successful it returns the constant as `INTEGER`, `LONG`, or `REAL` (according to the definition of the constant)

**Example:**

```
INTEGER Iresult, Debug  
LONG Units  
Debug=3 ! Shows comprehensive debug messages (test mode)  
Units=FNHx_get_constant("DAQmx_Val_Volts", Iresult, Debug) ! get LONG constant  
IF Iresult=0 THEN  
    !... weitere Befehle ...  
END IF
```

### **FNHx\_get\_constant\$**

**Purpose:** Returns the value of a constant from NIDAQmx.h from the NI-DAQmx Library as a string

**HTBasic Definition:**

```
Const$=FNHx_get_constant$(Constant$, INTEGER Result, Rtype$, OPTIONAL INTEGER Debug)
```

Input:

```
Constant$       Name of DAQmx constant as defined in file "NIDAQmx.h"  
Debug           Shows debug infos (constant found / not found etc.)  
                Debug = 3 : Show all debug messages
```

Return value: 0 = Function was successful  
                  **Result** = 1: Function was successful 0: Function failed  
<>0 = Error occurred

If the function was successful it returns the constant as a string and the variable type in `Rtype$` ("REAL", "SHORTHEX", "LONGHEX", "STRING", "LONG")

Use this function in order to retrieve string constants, e.g. for constants for switch topologies.

## SUB Hx\_show\_errstat

**Purpose:** Shows comprehensive information about current error status

### HTBasic Definition:

SUB Hx\_show\_errstat (LONG Result, OPTIONAL Fname\$, INTEGER Verbose)

#### Input:

Result                      Result of last DAQmx function HTBasic function  
Fname\$                      Name of function (optional)  
Verbose                      Verbose messages (0 = none, 1 = all, 2 = errors only)

Return value:    0 = Function was successful  
                    <>0 = Error occurred

## FNHx\_getdaqmxvers\$

**Purpose:** Returns version number of NI-DAQmx driver DLL "nicaiu.dll"

### HTBasic Definition:

Daqvers\$=FNHx\_getdaqmxvers\$

Return value:    If function was successful it returns version number, e.g. "8.0"  
                    If function failed the returned string is ""

## FNHx\_getdllvers\$

**Purpose:** Returns version number of DAQmx\_htb wrapper DLL "DAQmx\_htb.dll"

### HTBasic Definition:

Dllvers\$=FNHx\_getdllvers\$

Return value:    If function was successful it returns version number, e.g. "1.01.085"  
                    If function failed the returned string is ""

## FNHx\_getdevname\$

**Purpose:** Extracts device name from device string or device path (strips leading "/" if necessary)

### HTBasic Definition:

Dname\$=FNHx\_getdevname\$ (Devname\$)

#### Input:

Devname\$                    Name of device or path, e.g. "Dev1" or "/Dev1/ao0"

Return value:    If function was successful it returns device name, e.g. "Dev1"  
                    If function failed the returned string is ""

## FNHx\_lib\_version\$

**Purpose:** Returns version number of HTBasic DAQmx library "DAQmx\_htb.lib"

### HTBasic Definition:

Libvers\$=FNHx\_lib\_version\$

Return value:    If function was successful it returns version number, e.g. "1.03.116"  
                    If function failed the returned string is ""

## FNHx\_init

**Purpose:** Initializes HTBasic function library

**Note:** This function will be automatically called from FNHx\_manage\_subs if you load the function library, so normally it is *not* necessary to call this function.

### HTBasic Definition:

LONG Result=FNHx\_init

Return value:    1 = Function was successful  
                    0 = Error occurred

## **FNHx\_initlib**

**Purpose:** Initializes the DAQmx wrapper DLL "DAQmx\_htb.dll" (clears memory etc.)

**Note:** This function will be automatically called from FNHx\_manage\_subs / FNHx\_init if you load the function library, so normally it is *not* necessary to call this function.

### **HTBasic Definition:**

`LONG Result=FNHx_initlib`

**Return value:** 0 = Function was successful  
1 = Error occurred

## **SUB Hx\_enter\_kbd**

**Purpose:** Enhanced ENTER KBD function. Uses windows dialog box if Basic Plus is loaded

### **HTBasic Definition:**

`SUB Hx_enter_kbd(Text$)`

**Input:**  
Text\$                      Text to display

**Return value:** none

## **FNHx\_get\_dev\_info**

**Purpose:** Retrieves a number of interesting device parameters and returns it in a string array and (optional) in a file

### **HTBasic Definition:**

`LONG Result=FNHx_get_dev_info(Devname$, Dev_info$, OPTIONAL INTEGER Verbose, Write2File)`

**Input:**  
Devname\$                      Name of device or path, e.g. "Dev1" or "/Dev1/ao0"  
Verbose                      Verbose messages (0 = none, 1 = all, 2 = errors only)  
Write2File                      1: Write device info to file "daqmx\_devinfo.txt" in program directory  
                                    0: (or omit variable) Do not write device info to file

**Return value:** 1 = Function was successful  
                                    If function was successful it returns comprehensive device info in variable **Dev\_info\$**. Make sure that this variable is dimensioned big enough to hold info (at least 400 characters).  
0 = Error occurred  
                                    If function failed the returned string is "" or truncated string if **Dev\_info\$** is too small to hold complete device info.

## **FNHx\_gettnamwpref**

NI-DAQmx: **GetTerminalNameWithDevPrefix**

**Purpose:** Retrieves terminal name with prefix information for the given task

### **HTBasic Definition:**

`LONG Result=FNHx_gettnamwpref((LONG Taskid), Terminalname$, Triggername$)`

**Input:**  
Taskid                      The task used in this function

**Return value:** 0 = Function was successful  
                                    **Triggername\$** contains the trigger name  
                                    **Terminalname\$** contains Terminal name  
<>0 = Error occurred

## **FNHx\_numofdevices**

**Purpose:** Retrieves the number of installed devices in NI MAX (Measurement & Automation Explorer)

### **HTBasic Definition:**

`LONG Result=FNHx_numofdevices(OPTIONAL Devnames$(*))`

**Return value:** >0 = (Function was successful): Number of devices  
                                    **Devnames\$** contains the list of device names available in MAX  
0 = Error occurred

**FNHx\_getsfpath\$****Purpose:** Returns the path to special windows directories.**HTBasic Definition:**Sfpath\$=FNHx\_getsfpath\$(Foldercsid\$,Subdir\$,**INTEGER** Dflag,**OPTIONAL INTEGER** Verbose)Input: Foldercsid\$ = Path CSIDLs

<b>Foldercsid\$ value</b>	<b>English Windows</b>	<b>(German translation)</b>
CSIDL_DESKTOP	Desktop	(Arbeitsplatz)
CSIDL_INTERNET	Internet	
CSIDL_PROGRAMS	Start menu: Programs	(Startmenü: Programme)
CSIDL_CONTROLS	Control Panel	(Systemsteuerung)
CSIDL_PRINTERS	Printer	(Drucker)
CSIDL_PERSONAL	My Documents	(Eigene Dateien)
CSIDL_FAVORITES	IE: Favorites	(Favoriten)
CSIDL_STARTUP	Autostart	
CSIDL_RECENT	Recent documents	(Zuletzt benutzte Dokumente)
CSIDL_SENDTO	SendTo	(Senden an)
CSIDL_BITBUCKET	Recycle Bin	(Papierkorb)
CSIDL_STARTMENU	Startmenu	(Startmenü)
CSIDL_MYMUSIC	My Music	(Eigene Musik)
CSIDL_MYVIDEO	My Video	(Eigene Videos)
CSIDL_DESKTOPDIRECTORY	Desktop Directory	(Desktopverzeichnis)
CSIDL_DRIVES	My Computer	(Mein Computer)
CSIDL_NETWORK	Network	(Netzwerk)
CSIDL_NETHOOD	Network Neighbourhood	(Netzwerkumgebung)
CSIDL FONTS	Windows\Fonts	
CSIDL_TEMPLATES	Templates	(Vorlagen)
CSIDL_COMMON_STARTMENU	"All Users" - Startmenu	(Startmenü)
CSIDL_COMMON_PROGRAMS	"All Users" - Programs	(Programme)
CSIDL_COMMON_STARTUP	"All Users" - Autostart	
CSIDL_COMMON_DESKTOPDIRECTORY	"All Users" - Desktop	
CSIDL_APPDATA	Application Data	(Anwendungsdaten)
CSIDL_PRINTHOOD	Printhood	(Druckumgebung)
CSIDL_LOCAL_APPDATA	Local Settings\Application Data	(Lokale Einstellungen\Anwendungsdaten)
CSIDL_COMMON_FAVORITES	"All Users" - Favorites	(Favoriten)
CSIDL_INTERNET_CACHE	IE: Temorary Internet files	(Temporäre Internetdateien)
CSIDL_COOKIES	IE: Cookies	
CSIDL_HISTORY	IE: History	(Verlauf)
CSIDL_COMMON_APPDATA	"All Users" - Application Data	(Anwendungsdaten)
CSIDL WINDOWS	Windows	
CSIDL_SYSTEM	Windows\System32	
CSIDL_PROGRAM_FILES	C:\Program Files	(C:\Programme)
CSIDL_MYPICTURES	My Pictures	(Eigene Bilder)
CSIDL_PROFILE	Profile (Anwenderprofil	(Benutzername))
CSIDL_SYSTEMX86	Windows\System32	
CSIDL_PROGRAM_FILES_COMMON	Common Files	(Gemeinsame Dateien)
CSIDL_COMMON_TEMPLATES	"All Users" - Templates	(Vorlagen)
CSIDL_COMMON_DOCUMENTS	"All Users" - Documents	(Dokumente)
CSIDL_COMMON_ADMINTOOLS	"All Users" - AdminTools	(Verwaltung)
CSIDL_ADMINTOOLS	Startmenu\Program Files\AdminTools	(Startmenü\Programme\Verwaltung)

Subdir\$ = Subdirectory to be added to directory name

Dflag = Flags.

Currently defined: 1 = Create subdir if it does not exist

Return value: **Sfpath\$:** Returns the path or "" if the path does not exist

# Index

## B

Buffer Functions • 85

## C

Calibration Functions • 55

Callback Services • 11, 165

Channel configuration / Creation • 16

Concept • 2

## D

DAQmxAdjust1102Cal • 71  
DAQmxAdjust1104Cal • 63  
DAQmxAdjust1112Cal • 63  
DAQmxAdjust1122Cal • 74  
DAQmxAdjust1124Cal • 64  
DAQmxAdjust1125Cal • 64  
DAQmxAdjust1126Cal • 66  
DAQmxAdjust1141Cal • 66  
DAQmxAdjust1142Cal • 68  
DAQmxAdjust1143Cal • 68  
DAQmxAdjust1502Cal • 69  
DAQmxAdjust1503Cal • 69  
DAQmxAdjust1503CurrentCal • 70  
DAQmxAdjust1520Cal • 70  
DAQmxAdjust1521Cal • 74  
DAQmxAdjust153xCal • 71  
DAQmxAdjust1540Cal • 73  
DAQmxAdjust4204Cal • 59  
DAQmxAdjust4220Cal • 59  
DAQmxAdjust4224Cal • 60  
DAQmxAdjust4225Cal • 60  
DAQmxAdjustDSAAICal • 58  
DAQmxAdjustDSAAOCal • 58  
DAQmxAdjustDSATimebaseCal • 59  
DAQmxAOSeriesCalAdjust • 62  
DAQmxCfgAnlgEdgeRefTrig • 52  
DAQmxCfgAnlgEdgeStartTrig • 51  
DAQmxCfgAnlgWindowRefTrig • 53  
DAQmxCfgAnlgWindowStartTrig • 51  
DAQmxCfgBurstHandshakingTimingExportClock • 38  
DAQmxCfgBurstHandshakingTimingImportClock • 38  
DAQmxCfgChangeDetectionTiming • 39  
DAQmxCfgDigEdgeAdvTrig • 54  
DAQmxCfgDigEdgeRefTrig • 52  
DAQmxCfgDigEdgeStartTrig • 51  
DAQmxCfgDigPatternRefTrig • 53  
DAQmxCfgDigPatternStartTrig • 52  
DAQmxCfgHandshakingTiming • 37  
DAQmxCfgImplicitTiming • 39  
DAQmxCfgSampClkTiming • 36  
DAQmxChangeExtCalPassword • 57  
DAQmxClearTask • 14  
DAQmxClearTEDS • 83  
DAQmxCloseExtCal • 56  
DAQmxConfigureTEDS • 83  
DAQmxConnectTerms • 82  
DAQmxControlWatchdogTask • 167  
DAQmxCreateAIAccelChan • 20  
DAQmxCreateAICurrentChan • 16  
DAQmxCreateAICurrentRMSChan • 34  
DAQmxCreateAIFreqVoltageChan • 18

DAQmxCreateAIMicrophoneChan • 20  
DAQmxCreateAIPosLVDTChan • 21  
DAQmxCreateAIPosRVDTChan • 21  
DAQmxCreateAIResistanceChan • 19  
DAQmxCreateAIRTDCChan • 17  
DAQmxCreateAISTrainGageChan • 19  
DAQmxCreateAITempBuiltInSensorChan • 20  
DAQmxCreateAIThrmcplChan • 17  
DAQmxCreateAIThrmstrChanIex • 18  
DAQmxCreateAIThrmstrChanVex • 18  
DAQmxCreateAIVoltageChan • 16  
DAQmxCreateAIVoltageChanWithExcit • 19  
DAQmxCreateAIVoltageRMSChan • 34  
DAQmxCreateAOCurrentChan • 26  
DAQmxCreateAOFuncGenChan • 35  
DAQmxCreateAOVoltageChan • 26  
DAQmxCreateCIAngEncoderChan • 30  
DAQmxCreateCICountEdgesChan • 29  
DAQmxCreateCIFreqChan • 28  
DAQmxCreateCIGPSTimestampChan • 31  
DAQmxCreateCILinEncoderChan • 30  
DAQmxCreateCIPeriodChan • 28  
DAQmxCreateCIPulseWidthChan • 29  
DAQmxCreateCISemiPeriodChan • 29  
DAQmxCreateCITwoEdgeSepChan • 30  
DAQmxCreateCOPulseChanFreq • 31  
DAQmxCreateCOPulseChanTicks • 33  
DAQmxCreateCOPulseChanTime • 31  
DAQmxCreateDIChan • 27  
DAQmxCreateDOChan • 27  
DAQmxCreateTEDSAIAccelChan • 24  
DAQmxCreateTEDSAICurrentChan • 22  
DAQmxCreateTEDSAIMicrophoneChan • 25  
DAQmxCreateTEDSAIPosLVDTChan • 25  
DAQmxCreateTEDSAIPosRVDTChan • 25  
DAQmxCreateTEDSAIResistanceChan • 23  
DAQmxCreateTEDSAIRTDCChan • 22  
DAQmxCreateTEDSAISTrainGageChan • 24  
DAQmxCreateTEDSAIThrmcplChan • 22  
DAQmxCreateTEDSAIThrmstrChanIex • 23  
DAQmxCreateTEDSAIThrmstrChanVex • 23  
DAQmxCreateTEDSAIVoltageChan • 21  
DAQmxCreateTEDSAIVoltageChanWithExcit • 24  
DAQmxCreateWatchdogTimerTask • 167  
DAQmxDeleteSavedTask • 168  
DAQmxDeviceSupportsCal • 62  
DAQmxDisableAdvTrig • 53  
DAQmxDisableRefTrig • 52  
DAQmxDisableStartTrig • 51  
DAQmxDisconnectTerms • 82  
DAQmxESeriesCalAdjust • 60  
DAQmxExportSignal • 84  
DAQmxGetAdvTrigType • 94  
DAQmxGetAIAccelDBRef • 101  
DAQmxGetAIACExcitSyncEnable • 89  
DAQmxGetAIADCTimingMode • 102  
DAQmxGetAIChanCalPolyForwardCoeff • 75  
DAQmxGetAIConvActiveEdgeEx • 138  
DAQmxGetAIConvMaxRateEx • 137  
DAQmxGetAIConvRateEx • 136  
DAQmxGetAIConvSrcEx • 137  
DAQmxGetAIConvTimebaseDivEx • 138  
DAQmxGetAIConvTimebaseSrcEx • 139  
DAQmxGetAICurrentACRMSUnits • 152  
DAQmxGetAIDevScalingCoeff • 76  
DAQmxGetAIIsTEDS • 101  
DAQmxGetAILowpassCutoffFreq • 87

DAQmxGetAIIowpassEnable • 88  
 DAQmxGetAIMax • 87  
 DAQmxGetAIMin • 87  
 DAQmxGetAISoundPressuredBRef • 100  
 DAQmxGetAIThrmcplScaleType • 99  
 DAQmxGetAIVoltageACRMSUnits • 151  
 DAQmxGetAIVoltagedBRef • 98  
 DAQmxGetAnlgEdgeStartTrigHyst • 88  
 DAQmxGetAOFuncGenAmplitude • 156  
 DAQmxGetAOFuncGenFMDeviation • 160  
 DAQmxGetAOFuncGenFreq • 155  
 DAQmxGetAOFuncGenModulationType • 159  
 DAQmxGetAOFuncGenOffset • 157  
 DAQmxGetAOFuncGenSquareDutyCycle • 158  
 DAQmxGetAOFuncGenType • 154  
 DAQmxGetAOVoltageCurrentLimit • 153  
 DAQmxGetBufInputBufSize • 85  
 DAQmxGetBufOutputBufSize • 85  
 DAQmxGetBufOutputOnbrdBufSize • 86  
 DAQmxGetCalDevTemp • 55  
 DAQmxGetCIGPSSyncSrc • 91  
 DAQmxGetCOConstrainedGenMode • 112  
 DAQmxGetDelayFromSampClkDelay • 142  
 DAQmxGetDelayFromSampClkDelayEx • 142  
 DAQmxGetDelayFromSampClkDelayUnits • 140  
 DAQmxGetDelayFromSampClkDelayUnitsEx • 141  
 DAQmxGetDevAICouplings • 128  
 DAQmxGetDevAICurrentIntExcitDiscreteVals • 127  
 DAQmxGetDevAICurrentRngs • 127  
 DAQmxGetDevAIFreqRngs • 127  
 DAQmxGetDevAIGains • 127  
 DAQmxGetDevAIIowpassCutoffFreqDiscreteVals • 128  
 DAQmxGetDevAIIowpassCutoffFreqRangeVals • 128  
 DAQmxGetDevAIMaxMultiChanRate • 125  
 DAQmxGetDevAIMaxSingleChanRate • 125  
 DAQmxGetDevAIMinRate • 125  
 DAQmxGetDevAISimultaneousSamplingSupported • 125  
 DAQmxGetDevAITrigUsage • 126  
 DAQmxGetDevAIVoltageIntExcitDiscreteVals • 126  
 DAQmxGetDevAIVoltageIntExcitRangeVals • 126  
 DAQmxGetDevAIVoltageRngs • 126  
 DAQmxGetDevAnlgTrigSupported • 124  
 DAQmxGetDevAOCurrentRngs • 130  
 DAQmxGetDevAOGains • 130  
 DAQmxGetDevAOMaxRate • 129  
 DAQmxGetDevAOMinRate • 129  
 DAQmxGetDevAOPhysicalChans • 128  
 DAQmxGetDevAOSampClkSupported • 129  
 DAQmxGetDevAOTrigUsage • 129  
 DAQmxGetDevAOVoltageRngs • 130  
 DAQmxGetDevChassisModuleDevNames • 124  
 DAQmxGetDevCIMaxSize • 132  
 DAQmxGetDevCIMaxTimebase • 132  
 DAQmxGetDevCISampClkSupported • 132  
 DAQmxGetDevCITrigUsage • 131  
 DAQmxGetDevCOMaxSize • 133  
 DAQmxGetDevCOMaxTimebase • 133  
 DAQmxGetDevCompactDAQChassisDevName • 133  
 DAQmxGetDevCompactDAQSlotNum • 134  
 DAQmxGetDevCOTrigUsage • 132  
 DAQmxGetDevDigTrigSupported • 124  
 DAQmxGetDevDIMaxRate • 130  
 DAQmxGetDevDITrigUsage • 131  
 DAQmxGetDevDOMaxRate • 131  
 DAQmxGetDevDOTrigUsage • 131  
 DAQmxGetDevNumDMAChans • 133  
 DAQmxGetDevProductCategory • 124  
 DAQmxGetDevTerminals • 163  
 DAQmxGetDIAcquireOn • 105  
 DAQmxGetDIDigFltrEnable • 93  
 DAQmxGetDIDigFltrMinPulseWidth • 93  
 DAQmxGetDigLvlPauseTrigSrc • 90  
 DAQmxGetDigLvlPauseTrigWhen • 91  
 DAQmxGetDigPatternPauseTrigPattern • 144  
 DAQmxGetDigPatternPauseTrigSrc • 143  
 DAQmxGetDigPatternPauseTrigWhen • 145  
 DAQmxGetDILogicFamily • 103  
 DAQmxGetDIMemMapEnable • 104  
 DAQmxGetDOGenerateOn • 111  
 DAQmxGetDOLineStatesDoneState • 108  
 DAQmxGetDOLineStatesPausedState • 107  
 DAQmxGetDOLineStatesStartState • 106  
 DAQmxGetDOLogicFamily • 109  
 DAQmxGetDOMemMapEnable • 110  
 DAQmxGetExportedDataActiveEventLvlActiveLvl • 121  
 DAQmxGetExportedDataActiveEventOutputTerm • 120  
 DAQmxGetExportedHshkEventInterlockedDeassertDelay • 92  
 DAQmxGetExportedPauseTrigLvlActiveLvl • 115  
 DAQmxGetExportedPauseTrigOutputTerm • 114  
 DAQmxGetExportedRdyForStartEventLvlActiveLvl • 123  
 DAQmxGetExportedRdyForStartEventOutputTerm • 122  
 DAQmxGetExportedRdyForXferEventDeassertCond • 118  
 DAQmxGetExportedRdyForXferEventDeassertCondCustomThreshold • 119  
 DAQmxGetExportedRefTrigPulsePolarity • 116  
 DAQmxGetExportedSampClkDelayOffset • 113  
 DAQmxGetExportedStartTrigPulsePolarity • 117  
 DAQmxGetExtCallLastDateAndTime • 57  
 DAQmxGetExtCallLastTemp • 56  
 DAQmxGetExtCalRecommendedInterval • 58  
 DAQmxGetHshkDelayAfterXfer • 92  
 DAQmxGetMasterTimebaseRate • 95  
 DAQmxGetMasterTimebaseSrc • 95  
 DAQmxGetNthTaskChannel • 15  
 DAQmxGetPauseTrigType • 90  
 DAQmxGetPhysicalChanAITermCfgs • 148  
 DAQmxGetPhysicalChanAOManualControlAmplitude • 163  
 DAQmxGetPhysicalChanAOManualControlEnable • 162  
 DAQmxGetPhysicalChanAOManualControlFreq • 163  
 DAQmxGetPhysicalChanAOTermCfgs • 148  
 DAQmxGetPhysicalChanDIChangeDetectSupported • 149  
 DAQmxGetPhysicalChanDIPortWidth • 148  
 DAQmxGetPhysicalChanDISampClkSupported • 149  
 DAQmxGetPhysicalChanDOPortWidth • 149  
 DAQmxGetPhysicalChanDOSampClkSupported • 149  
 DAQmxGetPhysicalChanName • 77  
 DAQmxGetReadOpenCurrentLoopChans • 135  
 DAQmxGetReadOpenCurrentLoopChansExist • 134

DAQmxGetReadOvercurrentChans • 134  
 DAQmxGetReadOvercurrentChansExist • 134  
 DAQmxGetReadOverloadedChans • 89  
 DAQmxGetReadOverloadedChansExist • 89  
 DAQmxGetRefClkRate • 96  
 DAQmxGetRefClkSrc • 96  
 DAQmxGetSampClkTimebaseSrc • 97  
 DAQmxGetSampClkTimingResponseMode • 150  
 DAQmxGetSampClkUnderflowBehavior • 135  
 DAQmxGetSampClkTimingEngine • 161  
 DAQmxGetSelfCalLastDateAndTime • 57  
 DAQmxGetSelfCalLastTemp • 56  
 DAQmxGetSelfCalSupported • 62  
 DAQmxGetSwitchScanRepeatMode • 94  
 DAQmxGetSyncPulseSrc • 97  
 DAQmxGetTaskNumDevices • 135  
 DAQmxGetWriteNextWriteIsLast • 147  
 DAQmxGetWriteOpenCurrentLoopChans • 146  
 DAQmxGetWriteOpenCurrentLoopChansExist • 146  
 DAQmxGetWriteOvercurrentChans • 146  
 DAQmxGetWriteOvercurrentChansExist • 146  
 DAQmxGetWritePowerSupplyFaultChans • 147  
 DAQmxGetWritePowerSupplyFaultChansExist • 147  
 DAQmxInitExtCal • 56  
 DAQmxIsTaskDone • 15  
 DAQmxLoadTask • 13  
 DAQmxMSeriesCalAdjust • 61  
 DAQmxPerformBridgeOffsetNullingCal • 55  
 DAQmxPerformBridgeOffsetNullingCalEx • 74  
 DAQmxPerformBridgeShuntCal • 75  
 DAQmxPerformStrainShuntCal • 75  
 DAQmxReadAnalogF64 • 40  
 DAQmxReadAnalogScalarF64 • 40  
 DAQmxReadBinaryI16 • 41  
 DAQmxReadBinaryI32 • 41  
 DAQmxReadCounterF64 • 43  
 DAQmxReadCounterScalarF64 • 44  
 DAQmxReadCounterScalarU32 • 44  
 DAQmxReadCounterU32 • 44  
 DAQmxReadDigitalLines • 43  
 DAQmxReadDigitalScalarU32 • 42  
 DAQmxReadDigitalU16 • 41  
 DAQmxReadDigitalU32 • 42  
 DAQmxRegisterDoneEvent • 166  
 DAQmxRegisterEveryNSamplesEvent • 165  
 DAQmxRegisterSignalEvent • 166  
 DAQmxResetAIAccelDBRef • 101  
 DAQmxResetAIADCTimingMode • 102  
 DAQmxResetAIChanCalPolyForwardCoeff • 76  
 DAQmxResetAIConvActiveEdgeEx • 138  
 DAQmxResetAIConvRateEx • 136  
 DAQmxResetAIConvSrcEx • 137  
 DAQmxResetAIConvTimebaseDivEx • 139  
 DAQmxResetAIConvTimebaseSrcEx • 140  
 DAQmxResetAICurrentACRMSUnits • 152  
 DAQmxResetAISoundPressuredBRef • 100  
 DAQmxResetAIThrmcplScaleType • 99  
 DAQmxResetAIVoltageACRMSUnits • 151  
 DAQmxResetAIVoltagedBRef • 98  
 DAQmxResetAOFuncGenAmplitude • 156  
 DAQmxResetAOFuncGenFMDDeviation • 160  
 DAQmxResetAOFuncGenFreq • 155  
 DAQmxResetAOFuncGenModulationType • 159  
 DAQmxResetAOFuncGenOffset • 157  
 DAQmxResetAOFuncGenSquareDutyCycle • 158  
 DAQmxResetAOFuncGenType • 154  
 DAQmxResetAOVoltageCurrentLimit • 153  
 DAQmxResetCOConstrainedGenMode • 112  
 DAQmxResetDelayFromSampClkDelay • 142  
 DAQmxResetDelayFromSampClkDelayEx • 143  
 DAQmxResetDelayFromSampClkDelayUnits • 141  
 DAQmxResetDelayFromSampClkDelayUnitsEx • 141  
 DAQmxResetDevice • 78  
 DAQmxResetDIAcquireOn • 105  
 DAQmxResetDigPatternPauseTrigPattern • 144  
 DAQmxResetDigPatternPauseTrigSrc • 144  
 DAQmxResetDigPatternPauseTrigWhen • 145  
 DAQmxResetDILogicFamily • 103  
 DAQmxResetDIMemMapEnable • 104  
 DAQmxResetDOGenerateOn • 111  
 DAQmxResetDOLineStatesDoneState • 108  
 DAQmxResetDOLineStatesPausedState • 107  
 DAQmxResetDOLineStatesStartState • 106  
 DAQmxResetDOLogicFamily • 109  
 DAQmxResetDOMemMapEnable • 110  
 DAQmxResetExportedDataActiveEventLvlActiveLvl • 121  
 DAQmxResetExportedDataActiveEventOutputTerm • 120  
 DAQmxResetExportedPauseTrigLvlActiveLvl • 115  
 DAQmxResetExportedPauseTrigOutputTerm • 114  
 DAQmxResetExportedRdyForStartEventLvlActiveLvl • 123  
 DAQmxResetExportedRdyForStartEventOutputTerm • 122  
 DAQmxResetExportedRdyForXferEventDeassertCond • 118  
 DAQmxResetExportedRdyForXferEventDeassertCondCustomThreshold • 119  
 DAQmxResetExportedRefTrigPulsePolarity • 116  
 DAQmxResetExportedSampClkDelayOffset • 113  
 DAQmxResetExportedStartTrigPulsePolarity • 117  
 DAQmxResetPhysicalChanAOManualControlEnable • 162  
 DAQmxResetSampClkTimingResponseMode • 150  
 DAQmxResetSampClkUnderflowBehavior • 136  
 DAQmxResetSampTimingEngine • 161  
 DAQmxResetWriteNextWriteIsLast • 148  
 DAQmxRestoreLastExtCalConst • 57  
 DAQmxSaveTask • 168  
 DAQmxSCBaseboardCalAdjust • 62  
 DAQmxSelfCal • 55  
 DAQmxSendSoftwareTrigger • 54  
 DAQmxSetAdvTrigType • 94  
 DAQmxSetAIAccelDBRef • 101  
 DAQmxSetAIACExcitSyncEnable • 89  
 DAQmxSetAIADCTimingMode • 102  
 DAQmxSetAIChanCalPolyForwardCoeff • 76  
 DAQmxSetAIConvActiveEdgeEx • 138  
 DAQmxSetAIConvRateEx • 136  
 DAQmxSetAIConvSrcEx • 137  
 DAQmxSetAIConvTimebaseDivEx • 139  
 DAQmxSetAIConvTimebaseSrcEx • 140  
 DAQmxSetAICurrentACRMSUnits • 152  
 DAQmxSetAILowpassCutoffFreq • 87  
 DAQmxSetAILowpassEnable • 88  
 DAQmxSetAISoundPressuredBRef • 100  
 DAQmxSetAIThrmcplScaleType • 99  
 DAQmxSetAIVoltageACRMSUnits • 151  
 DAQmxSetAIVoltagedBRef • 98  
 DAQmxSetAnlgEdgeStartTrigHyst • 88

DAQmxSetAOFuncGenAmplitude • 156  
 DAQmxSetAOFuncGenFMDeviation • 160  
 DAQmxSetAOFuncGenFreq • 155  
 DAQmxSetAOFuncGenModulationType • 159  
 DAQmxSetAOFuncGenOffset • 157  
 DAQmxSetAOFuncGenSquareDutyCycle • 158  
 DAQmxSetAOFuncGenType • 154  
 DAQmxSetAOVoltageCurrentLimit • 153  
 DAQmxSetBufferAttribute • 164  
 DAQmxSetBufInputBufSize • 85  
 DAQmxSetBufOutputBufSize • 85  
 DAQmxSetBufOutputOnbrdBufSize • 86  
 DAQmxSetCIGPSSyncSrc • 91  
 DAQmxSetCIPrescaler • 77  
 DAQmxSetCOConstrainedGenMode • 112  
 DAQmxSetCOPulseDutyCyc • 164  
 DAQmxSetCOPulseFreq • 164  
 DAQmxSetDelayFromSampClkDelay • 142  
 DAQmxSetDelayFromSampClkDelayEx • 143  
 DAQmxSetDelayFromSampClkDelayUnits • 140  
 DAQmxSetDelayFromSampClkDelayUnitsEx • 141  
 DAQmxSetDIAcquireOn • 105  
 DAQmxSetDIDigFltrEnable • 93  
 DAQmxSetDIDigFltrMinPulseWidth • 93  
 DAQmxSetDigitalLogicFamilyPowerUpState • 108  
 DAQmxSetDigLvlPauseTrigSrc • 90  
 DAQmxSetDigLvlPauseTrigWhen • 91  
 DAQmxSetDigPatternPauseTrigPattern • 144  
 DAQmxSetDigPatternPauseTrigSrc • 143  
 DAQmxSetDigPatternPauseTrigWhen • 145  
 DAQmxSetDILogicFamily • 103  
 DAQmxSetDIMemMapEnable • 104  
 DAQmxSetDOGenerateOn • 111  
 DAQmxSetDOLineStatesDoneState • 108  
 DAQmxSetDOLineStatesPausedState • 107  
 DAQmxSetDOLineStatesStartState • 106  
 DAQmxSetDOLogicFamily • 109  
 DAQmxSetDOMemMapEnable • 110  
 DAQmxSetExportedDataActiveEventLvlActiveLvl • 121  
 DAQmxSetExportedDataActiveEventOutputTerm • 120  
 DAQmxSetExportedHshkEventInterlockedDeassertDelay • 92  
 DAQmxSetExportedPauseTrigLvlActiveLvl • 115  
 DAQmxSetExportedPauseTrigOutputTerm • 114  
 DAQmxSetExportedRdyForStartEventLvlActiveLvl • 123  
 DAQmxSetExportedRdyForStartEventOutputTerm • 122  
 DAQmxSetExportedRdyForXferEventDeassertCond • 118  
 DAQmxSetExportedRdyForXferEventDeassertCondCustomThreshold • 119  
 DAQmxSetExportedRefTrigPulsePolarity • 116  
 DAQmxSetExportedSampClkDelayOffset • 113  
 DAQmxSetExportedStartTrigPulsePolarity • 117  
 DAQmxSetHshkDelayAfterXfer • 92  
 DAQmxSetMasterTimebaseRate • 95  
 DAQmxSetMasterTimebaseSrc • 95  
 DAQmxSetPauseTrigType • 90  
 DAQmxSetPhysicalChanAOManualControlEnable • 162  
 DAQmxSetPhysicalChanName • 113  
 DAQmxSetRefClkRate • 96  
 DAQmxSetRefClkSrc • 96  
 DAQmxSetSampClkTimebaseSrc • 97  
 DAQmxSetSampClkTimingResponseMode • 150  
 DAQmxSetSampClkUnderflowBehavior • 135  
 DAQmxSetSampTimingEngine • 161  
 DAQmxSetSwitchScanRepeatMode • 94  
 DAQmxSetSyncPulseSrc • 97  
 DAQmxSetTrigAttribute • 164  
 DAQmxSetup1102Cal • 71  
 DAQmxSetup1104Cal • 63  
 DAQmxSetup1112Cal • 63  
 DAQmxSetup1122Cal • 73  
 DAQmxSetup1124Cal • 64  
 DAQmxSetup1125Cal • 64  
 DAQmxSetup1126Cal • 66  
 DAQmxSetup1141Cal • 66  
 DAQmxSetup1142Cal • 68  
 DAQmxSetup1143Cal • 68  
 DAQmxSetup1502Cal • 69  
 DAQmxSetup1503Cal • 69  
 DAQmxSetup1520Cal • 70  
 DAQmxSetup1521Cal • 74  
 DAQmxSetup153xCal • 71  
 DAQmxSetup1540Cal • 73  
 DAQmxSetWriteNextWriteIsLast • 147  
 DAQmxSSeriesCalAdjust • 61  
 DAQmxStartTask • 13  
 DAQmxStopTask • 14  
 DAQmxSwitchCloseRelays • 81  
 DAQmxSwitchConnect • 79  
 DAQmxSwitchConnectMulti • 79  
 DAQmxSwitchCreateScanList • 79  
 DAQmxSwitchDisconnect • 79  
 DAQmxSwitchDisconnectAll • 80  
 DAQmxSwitchFindPath • 80  
 DAQmxSwitchOpenRelays • 80  
 DAQmxSwitchSetTopologyAndReset • 80  
 DAQmxSwitchWaitForSettling • 81  
 DAQmxTaskControl • 15  
 DAQmxTristateOutputTerm • 82  
 DAQmxWaitUntilTaskIsDone • 14  
 DAQmxWriteAnalogF64 • 45  
 DAQmxWriteAnalogScalarF64 • 45  
 DAQmxWriteBinaryI16 • 46  
 DAQmxWriteBinaryI32 • 46  
 DAQmxWriteCtrFreq • 48  
 DAQmxWriteCtrFreqScalar • 49  
 DAQmxWriteCtrTicks • 50  
 DAQmxWriteCtrTicksScalar • 50  
 DAQmxWriteCtrTime • 49  
 DAQmxWriteCtrTimeScalar • 49  
 DAQmxWriteDigitalLines • 47  
 DAQmxWriteDigitalScalarU32 • 47  
 DAQmxWriteDigitalU16 • 46  
 DAQmxWriteDigitalU32 • 47  
 Device Control Functions • 78  
 Device info • 3  
 Done Event • 165  
  
**E**  
 Events & Signals • 84  
 Every N Samples Event • 165  
 Example programs • 2  
  
**F**  
 FN Hx\_chnextcalpwd • 57  
 FNHx\_adj1102cal • 71



FNHx\_adj1104cal • 63  
 FNHx\_adj1112cal • 63  
 FNHx\_adj1122cal • 74  
 FNHx\_adj1124cal • 64  
 FNHx\_adj1125cal • 64  
 FNHx\_adj1126cal • 66  
 FNHx\_adj1141cal • 66  
 FNHx\_adj1142cal • 68  
 FNHx\_adj1143cal • 68  
 FNHx\_adj1502cal • 69  
 FNHx\_adj1503cal • 69  
 FNHx\_adj1503ccl • 70  
 FNHx\_adj1520cal • 70  
 FNHx\_adj1521cal • 74  
 FNHx\_adj153xcal • 71  
 FNHx\_adj1540cal • 73  
 FNHx\_adj4204cal • 59  
 FNHx\_adj4220cal • 59  
 FNHx\_adj4224cal • 60  
 FNHx\_adj4225cal • 60  
 FNHx\_adjdsaaical • 58  
 FNHx\_adjdsaaocal • 58  
 FNHx\_adjdsatbcal • 59  
 FNHx\_aosecaladjst • 62  
 FNHx\_cfdgedsttrig • 51  
 FNHx\_cfganedrtrig • 52  
 FNHx\_cfganedstrig • 51  
 FNHx\_cfganwirtrig • 53  
 FNHx\_cfganwistrig • 51  
 FNHx\_cfgbhtimeclk • 38  
 FNHx\_cfgbhtimiclk • 38  
 FNHx\_cfgchgdettim • 39  
 FNHx\_cfgdiedatrig • 54  
 FNHx\_cfgdiedrtrig • 52  
 FNHx\_cfgdipartrig • 53  
 FNHx\_cfgdiptstrig • 52  
 FNHx\_cfghshtiming • 37  
 FNHx\_cfgimptiming • 39  
 FNHx\_cfgsamclktim • 36  
 FNHx\_clearatask • 14  
 FNHx\_clearateds • 83  
 FNHx\_closeextcal • 56  
 FNHx\_configteds • 83  
 FNHx\_connterms • 82  
 FNHx\_craiaccelchn • 20  
 FNHx\_craicurmschn • 34  
 FNHx\_craicurrchan • 16  
 FNHx\_craifrvochan • 18  
 FNHx\_craimicrchan • 20  
 FNHx\_craipslvdtdch • 21  
 FNHx\_craipsrvdtdch • 21  
 FNHx\_crairesichan • 19  
 FNHx\_crairttdchan • 17  
 FNHx\_craistrigachn • 19  
 FNHx\_craitacclchn • 24  
 FNHx\_craithcpchan • 17  
 FNHx\_craithmchiex • 18  
 FNHx\_craithmchvex • 18  
 FNHx\_craitmpbschn • 20  
 FNHx\_craivlrmschn • 34  
 FNHx\_craivochwexc • 19  
 FNHx\_craivoltchan • 16  
 FNHx\_craocurrchan • 26  
 FNHx\_craofungchn • 35  
 FNHx\_craovoltchan • 26  
 FNHx\_crcianencchn • 30  
 FNHx\_crcicntedgch • 29  
 FNHx\_crcifreqchan • 28  
 FNHx\_crcigpstschn • 31  
 FNHx\_crcilinencch • 30  
 FNHx\_crciperichan • 28  
 FNHx\_crciplswchan • 29  
 FNHx\_crcisemperch • 29  
 FNHx\_crcitwedsech • 30  
 FNHx\_crcopulchfrq • 31  
 FNHx\_crcopulchtck • 33  
 FNHx\_crcopulchtim • 31  
 FNHx\_createdichan • 27  
 FNHx\_createdochan • 27  
 FNHx\_createtask • 13  
 FNHx\_crtaicurrchn • 22  
 FNHx\_crtaimicchan • 25  
 FNHx\_crtaiplvdtdch • 25  
 FNHx\_crtaiprvdtch • 25  
 FNHx\_crtairesichn • 23  
 FNHx\_crtairtdchn • 22  
 FNHx\_crtairtrgchn • 24  
 FNHx\_crtaitchiex • 23  
 FNHx\_crtaitchvex • 23  
 FNHx\_crtaithcpchn • 22  
 FNHx\_crtaivochwex • 24  
 FNHx\_crtaivoltchn • 21  
 FNHx\_crwdogtimtsk • 167  
 FNHx\_ctrlwdogtsk • 167  
 FNHx\_delsavedtask • 168  
 FNHx\_devsuppcal • 62  
 FNHx\_disadvtrig • 53  
 FNHx\_disconnterms • 82  
 FNHx\_disstarttrig • 51  
 FNHx\_esercaladjst • 60  
 FNHx\_exportsignal • 84  
 FNHx\_get\_constant • 169  
 FNHx\_get\_constant\$ • 169  
 FNHx\_get\_dev\_info • 171  
 FNHx\_getadvtrgtyp • 94  
 FNHx\_getaiacbyref • 101  
 FNHx\_getaiacexsen • 89  
 FNHx\_getaiadctimd • 102  
 FNHx\_getaichcpfco • 75  
 FNHx\_getaicuacrmu • 152  
 FNHx\_getaicvaceex • 138  
 FNHx\_getaicvmrtex • 137  
 FNHx\_getaicvratex • 136  
 FNHx\_getaicvsrceex • 137  
 FNHx\_getaicvtbdex • 138  
 FNHx\_getaicvtbsex • 139  
 FNHx\_getaidevscco • 76  
 FNHx\_getaiisteds • 101  
 FNHx\_getailpcofrq • 87  
 FNHx\_getailwpenab • 88  
 FNHx\_getaimax • 87  
 FNHx\_getaimin • 87  
 FNHx\_getaispbyref • 100  
 FNHx\_getaitcsctyp • 99  
 FNHx\_getaivoacrmu • 151  
 FNHx\_getaivolbref • 98  
 FNHx\_getanedsthys • 88  
 FNHx\_getaoculimit • 153  
 FNHx\_getaofungamp • 156  
 FNHx\_getaofungfmd • 160  
 FNHx\_getaofungfrq • 155  
 FNHx\_getaofungmtp • 159  
 FNHx\_getaofungoff • 157  
 FNHx\_getaofungsdcc • 158  
 FNHx\_getaofungtyp • 154  
 FNHx\_getbufinpsiz • 85

FNHx\_getbufoutsiz • 85  
 FNHx\_getcadevtemp • 55  
 FNHx\_getcigpssysr • 91  
 FNHx\_getcocongmod • 112  
 FNHx\_getdaqmxvers\$ • 170  
 FNHx\_getdevaciedv • 127  
 FNHx\_getdevacurra • 127  
 FNHx\_getdevafrqrn • 127  
 FNHx\_getdevaifcou • 128  
 FNHx\_getdevaigain • 127  
 FNHx\_getdevaimmcr • 125  
 FNHx\_getdevaimrat • 125  
 FNHx\_getdevaimscr • 125  
 FNHx\_getdevaitrus • 126  
 FNHx\_getdevalcfdv • 128  
 FNHx\_getdevalcfrv • 128  
 FNHx\_getdevaocurg • 130  
 FNHx\_getdevaogain • 130  
 FNHx\_getdevaomart • 129  
 FNHx\_getdevaomirt • 129  
 FNHx\_getdevaophch • 128  
 FNHx\_getdevaoscsp • 129  
 FNHx\_getdevaotrus • 129  
 FNHx\_getdevaovorg • 130  
 FNHx\_getdevasssup • 125  
 FNHx\_getdevatrsup • 124  
 FNHx\_getdevaviedv • 126  
 FNHx\_getdevavierv • 126  
 FNHx\_getdevavrang • 126  
 FNHx\_getdevcdcdnm • 133  
 FNHx\_getdevcdslnm • 134  
 FNHx\_getdevchmdnm • 124  
 FNHx\_getdevcimasz • 132  
 FNHx\_getdevcimatb • 132  
 FNHx\_getdevciscsp • 132  
 FNHx\_getdevcitrus • 131  
 FNHx\_getdevcomasz • 133  
 FNHx\_getdevcomatb • 133  
 FNHx\_getdevcotrus • 132  
 FNHx\_getdevdimart • 130  
 FNHx\_getdevditrus • 131  
 FNHx\_getdevdomart • 131  
 FNHx\_getdevdotrus • 131  
 FNHx\_getdevdtrsup • 124  
 FNHx\_getdevname\$ • 170  
 FNHx\_getdevnmchn • 133  
 FNHx\_getdevprdcatt • 124  
 FNHx\_getdevtermin • 163  
 FNHx\_getdiacquion • 105  
 FNHx\_getdidgflmpw • 93  
 FNHx\_getdidgflten • 93  
 FNHx\_getdiglvptsc • 90  
 FNHx\_getdigpptpat • 144  
 FNHx\_getdigpptsrc • 143  
 FNHx\_getdigpptwhn • 145  
 FNHx\_getdilogifam • 103  
 FNHx\_getdimmapena • 104  
 FNHx\_getdlfscdlay • 142  
 FNHx\_getdlfscdlex • 142  
 FNHx\_getdlfscduent • 141  
 FNHx\_getdlfscdunt • 140  
 FNHx\_getdllvers\$ • 170  
 FNHx\_getdlpatrwhn • 91  
 FNHx\_getdogeneron • 111  
 FNHx\_getdologifam • 109  
 FNHx\_getdolsdone • 108  
 FNHx\_getdolspauss • 107  
 FNHx\_getdolsstrts • 106  
 FNHx\_getdommapena • 110  
 FNHx\_getecaldtati • 57  
 FNHx\_getecalrecin • 58  
 FNHx\_getexcaltemp • 56  
 FNHx\_getexdaceout • 120  
 FNHx\_getexdaevlal • 121  
 FNHx\_getexhseiddl • 92  
 FNHx\_getexptlaclv • 115  
 FNHx\_getexptouttm • 114  
 FNHx\_getexrdfxedc • 118  
 FNHx\_getexrdyfxed • 119  
 FNHx\_getexrfselal • 123  
 FNHx\_getexrfseotm • 122  
 FNHx\_getexrftrppo • 116  
 FNHx\_getexscdloff • 113  
 FNHx\_getexsttrppo • 117  
 FNHx\_gethsdelaxfr • 92  
 FNHx\_getmastbasrc • 95  
 FNHx\_getmastbrate • 95  
 FNHx\_getnhttaskch • 15  
 FNHx\_getobufoutsz • 86  
 FNHx\_getpaustrtyp • 90  
 FNHx\_getpcaomcamp • 163  
 FNHx\_getpcaomcfrq • 163  
 FNHx\_getphchaitcf • 148  
 FNHx\_getphchaomce • 162  
 FNHx\_getphchaotcf • 148  
 FNHx\_getphchdidsp • 149  
 FNHx\_getphchdipwd • 148  
 FNHx\_getphchdiscs • 149  
 FNHx\_getphchdopwd • 149  
 FNHx\_getphchdoscs • 149  
 FNHx\_getphyschnam • 77  
 FNHx\_getrdocuchex • 134  
 FNHx\_getrdoculcex • 134  
 FNHx\_getrdopculch • 135  
 FNHx\_getrdovchans • 89  
 FNHx\_getrdovchexi • 89  
 FNHx\_getrdovcuchn • 134  
 FNHx\_getrefclkrat • 96  
 FNHx\_getrefclksrc • 96  
 FNHx\_getsamtimeng • 161  
 FNHx\_getscaldtati • 57  
 FNHx\_getsclktbsrc • 97  
 FNHx\_getsctresmod • 150  
 FNHx\_getsfcaltemp • 56  
 FNHx\_getsfpath\$ • 172  
 FNHx\_getsmclufbeh • 135  
 FNHx\_getswscrpmod • 94  
 FNHx\_getsynpulsrsc • 97  
 FNHx\_gettnamwpref • 171  
 FNHx\_gettsknumdev • 135  
 FNHx\_getwrnxwrla • 147  
 FNHx\_getwropcchex • 146  
 FNHx\_getwropcchns • 146  
 FNHx\_getwrovcchex • 146  
 FNHx\_getwrovccchns • 146  
 FNHx\_getwrpsfchex • 147  
 FNHx\_getwrpsfchns • 147  
 FNHx\_init • 170  
 FNHx\_initextcal • 56  
 FNHx\_initlib • 171  
 FNHx\_istaskdone • 15  
 FNHx\_lib\_version\$ • 170  
 FNHx\_loadtask • 13  
 FNHx\_manage\_subs • 169  
 FNHx\_msercaladjst • 61  
 FNHx\_numofdevices • 171

FNHx\_perfbrishcal • 75  
 FNHx\_perfbrofncex • 74  
 FNHx\_perfbrofnlca • 55  
 FNHx\_perfstrshcal • 75  
 FNHx\_rdbinaryi16 • 41  
 FNHx\_rdbinaryi32 • 41  
 FNHx\_rdcntscal64 • 44  
 FNHx\_rdcntscal32 • 44  
 FNHx\_readanalf64 • 40  
 FNHx\_readanalsf64 • 40  
 FNHx\_readcntrf64 • 43  
 FNHx\_readcntru32 • 44  
 FNHx\_readdigits32 • 42  
 FNHx\_readdigitu16 • 41  
 FNHx\_readdigitu32 • 42  
 FNHx\_readdiglines • 43  
 FNHx\_regdoneevent • 166  
 FNHx\_regevnsmevnt • 165  
 FNHx\_regsignevent • 166  
 FNHx\_resetdevice • 78  
 FNHx\_reslexcalcon • 57  
 FNHx\_rstaiacbyref • 101  
 FNHx\_rstaiadctimd • 102  
 FNHx\_rstaichcpfco • 76  
 FNHx\_rstaicuacrmu • 152  
 FNHx\_rstaiavaceex • 138  
 FNHx\_rstaicvratex • 136  
 FNHx\_rstaicvsrcecx • 137  
 FNHx\_rstaicvtbdex • 139  
 FNHx\_rstaicvtbsex • 140  
 FNHx\_rstaispbyref • 100  
 FNHx\_rstaitcsctyp • 99  
 FNHx\_rstaivoacrmu • 151  
 FNHx\_rstaivolbref • 98  
 FNHx\_rstaoculimit • 153  
 FNHx\_rstaofungamp • 156  
 FNHx\_rstaofungfmd • 160  
 FNHx\_rstaofungfrq • 155  
 FNHx\_rstaofungmtp • 159  
 FNHx\_rstaofungoff • 157  
 FNHx\_rstaofungsdsc • 158  
 FNHx\_rstaofungtyp • 154  
 FNHx\_rstcocongmod • 112  
 FNHx\_rstdiacquion • 105  
 FNHx\_rstdigpptpat • 144  
 FNHx\_rstdigpptsr • 144  
 FNHx\_rstdigpptwhn • 145  
 FNHx\_rstdilogifam • 103  
 FNHx\_rstdimmapena • 104  
 FNHx\_rstdlfsdclay • 142  
 FNHx\_rstdlfsdclay • 143  
 FNHx\_rstdlfscduecx • 141  
 FNHx\_rstdlfscdunt • 141  
 FNHx\_rstdogeneron • 111  
 FNHx\_rstdologifam • 109  
 FNHx\_rstdolsdone • 108  
 FNHx\_rstdolspauss • 107  
 FNHx\_rstdolsstrts • 106  
 FNHx\_rstdommapena • 110  
 FNHx\_rstexdaceout • 120  
 FNHx\_rstexdaeval • 121  
 FNHx\_rstexptlaclv • 115  
 FNHx\_rstexptouttm • 114  
 FNHx\_rstexrdfxedc • 118  
 FNHx\_rstexrdyfxed • 119  
 FNHx\_rstexrfselal • 123  
 FNHx\_rstexrfseotm • 122  
 FNHx\_rstexrftrppo • 116  
 FNHx\_rstexscdloff • 113  
 FNHx\_rstexsttrppo • 117  
 FNHx\_rstphchaomce • 162  
 FNHx\_rstsamtimg • 161  
 FNHx\_rstscresmod • 150  
 FNHx\_rstsmclufbeh • 136  
 FNHx\_rstwrnxwila • 148  
 FNHx\_savetask • 168  
 FNHx\_scbbcaldadjst • 62  
 FNHx\_selfcal • 55  
 FNHx\_selfcalsupp • 62  
 FNHx\_sendsofttrig • 54  
 FNHx\_set1102cal • 71  
 FNHx\_set1104cal • 63  
 FNHx\_set1112cal • 63  
 FNHx\_set1122cal • 73  
 FNHx\_set1124cal • 64  
 FNHx\_set1125cal • 64  
 FNHx\_set1126cal • 66  
 FNHx\_set1141cal • 66  
 FNHx\_set1142cal • 68  
 FNHx\_set1143cal • 68  
 FNHx\_set1502cal • 69  
 FNHx\_set1503cal • 69  
 FNHx\_set1520cal • 70  
 FNHx\_set1521cal • 74  
 FNHx\_set153xcal • 71  
 FNHx\_set1540cal • 73  
 FNHx\_setadvtrgtyp • 94  
 FNHx\_setaiacbyref • 101  
 FNHx\_setaiacexsen • 89  
 FNHx\_setaiadctimd • 102  
 FNHx\_setaiachcpfco • 76  
 FNHx\_setaiacuacrmu • 152  
 FNHx\_setaiavaceex • 138  
 FNHx\_setaicvratex • 136  
 FNHx\_setaicvsrcecx • 137  
 FNHx\_setaicvtbdex • 139  
 FNHx\_setaicvtbsex • 140  
 FNHx\_setailpcofrq • 87  
 FNHx\_setailwpenab • 88  
 FNHx\_setaispbyref • 100  
 FNHx\_setaitcsctyp • 99  
 FNHx\_setaivoacrmu • 151  
 FNHx\_setaivolbref • 98  
 FNHx\_setanedsthys • 88  
 FNHx\_setaoaculimit • 153  
 FNHx\_setaofungamp • 156  
 FNHx\_setaofungfmd • 160  
 FNHx\_setaofungfrq • 155  
 FNHx\_setaofungmtp • 159  
 FNHx\_setaofungoff • 157  
 FNHx\_setaofungsdsc • 158  
 FNHx\_setaofungtyp • 154  
 FNHx\_setbufattrib • 164  
 FNHx\_setbufinpsiz • 85  
 FNHx\_setbufoutsiz • 85  
 FNHx\_setcigpssysr • 91  
 FNHx\_setciprescal • 77  
 FNHx\_setcocongmod • 112  
 FNHx\_setcopdutcyc • 164  
 FNHx\_setcopulfreq • 164  
 FNHx\_setdiacquion • 105  
 FNHx\_setdidgflmpw • 93  
 FNHx\_setdidgflten • 93  
 FNHx\_setdiglvptsc • 90  
 FNHx\_setdigpptpat • 144  
 FNHx\_setdigpptsr • 143

FNHx\_setdigpptwhn • 145  
 FNHx\_setdilogifam • 103  
 FNHx\_setdimmapena • 104  
 FNHx\_setdlfampust • 108  
 FNHx\_setdlfscdlay • 142  
 FNHx\_setdlfscdlex • 143  
 FNHx\_setdlfscduex • 141  
 FNHx\_setdlfscdunt • 140  
 FNHx\_setdlpatrwhn • 91  
 FNHx\_setdogeneron • 111  
 FNHx\_setdologifam • 109  
 FNHx\_setdolsdone • 108  
 FNHx\_setdolspauss • 107  
 FNHx\_setdolsstrts • 106  
 FNHx\_setdommapena • 110  
 FNHx\_setexdaceout • 120  
 FNHx\_setexdaevlal • 121  
 FNHx\_setexhseiddl • 92  
 FNHx\_setexptlaclv • 115  
 FNHx\_setexptouttm • 114  
 FNHx\_setexrdfxedc • 118  
 FNHx\_setexrdyfxed • 119  
 FNHx\_setexrfselal • 123  
 FNHx\_setexrfseotm • 122  
 FNHx\_setexrftrppo • 116  
 FNHx\_setexscdloff • 113  
 FNHx\_setexsttrppo • 117  
 FNHx\_sethsdelaxfr • 92  
 FNHx\_setmastbasrc • 95  
 FNHx\_setmastbrate • 95  
 FNHx\_setobufoutsz • 86  
 FNHx\_setpaustrtyp • 90  
 FNHx\_setphchaomce • 162  
 FNHx\_setphyschnam • 113  
 FNHx\_setrefclkrat • 96  
 FNHx\_setrefclksrc • 96  
 FNHx\_setsamtimeng • 161  
 FNHx\_setsclktbsrc • 97  
 FNHx\_setsctresmod • 150  
 FNHx\_setsmclufbeh • 135  
 FNHx\_setsynpulsrsc • 97  
 FNHx\_settriattrib • 164  
 FNHx\_setwrnxwrla • 147  
 FNHx\_ssercaladjst • 61  
 FNHx\_starttask • 13  
 FNHx\_stoptask • 14  
 FNHx\_swclosrelays • 81  
 FNHx\_swconnect • 79  
 FNHx\_swconnectmul • 79  
 FNHx\_swcrscanlist • 79  
 FNHx\_swdisconnall • 80  
 FNHx\_swdisconnect • 79  
 FNHx\_swfindpath • 80  
 FNHx\_swopenrelays • 80  
 FNHx\_swsettopares • 80  
 FNHx\_swwaitforset • 81  
 FNHx\_taskcontrol • 15  
 FNHx\_tristoutterm • 82  
 FNHx\_waituntaskdn • 14  
 FNHx\_wrctrfreq • 48  
 FNHx\_wrctrfrqscal • 49  
 FNHx\_wrctrticks • 50  
 FNHx\_wrctrtickssc • 50  
 FNHx\_wrctrtime • 49  
 FNHx\_wrctrtimscal • 49  
 FNHx\_writdigitu16 • 46  
 FNHx\_writdigitu32 • 47  
 FNHx\_writdiglines • 47

FNHx\_writdigscu32 • 47  
 FNHx\_writeanalf64 • 45  
 FNHx\_writeanasf64 • 45  
 FNHx\_writebini16 • 46  
 FNHx\_writebini32 • 46  
 Function Overview • 4

## G

Get/Set Functions • 87  
 Get/Set/Reset Functions • 77  
 Getting started • 3  
 Global variables • 11

## H

HTBasic specific functions • 169  
 HTBasic variable mapping • 12

## I

Installation • 2

## L

Library usage • 3

## M

MAX • 2

## N

NI-DAQmx C Reference Help • 2  
 NI-DAQmx Documentation • 12  
 NI-DAQmx Help • 2

## O

ON SIGNAL • 165

## R

Read Functions • 40  
 Return variables of DLL functions • 12

## S

Signal Event • 165  
 Signal Routing • 82  
 Storage • 168  
 SUB Hx\_enter\_kbd • 171  
 SUB Hx\_show\_errstat • 170  
 Switch Functions • 79

## T

Task configuration and control • 13  
 TEDS • 83  
 Timing functions • 36  
 Triggering Functions • 51

## W

Watchdog functions • 167  
 Write Functions • 45



**Tech Soft GmbH**  
Test & Measurement  
Karmeliterweg 114  
D - 13465 Berlin, Germany

ph +49 - 30 - 40608-0  
fax +49 - 30 - 40608-220  
e-mail: HTB@TechSoft.de

[www.TechSoft.de](http://www.TechSoft.de)